



BioSpec[®] / MedSpec[®] / PharmaScan[®]

MRI/MRS

***ParaVision[®]* Programming Course**

3-7 April 2006

Bruker BioSpin MRI GmbH

Rudolf-Plank-Str. 23

D-76275 Ettlingen

www.bruker-biospin.de



Table of Contents

1. General Information

- Program
- List of Participants
- City Map of Ettlingen
- City Map of Karlsruhe
- Timetable for Public Transportation
- General Questionnaire
- Customer Questionnaire

2. Courses

- Programming in *ParaVision*[®]
- Pulse Programs and Parameters
- Introduction to Pulse Programming (I+II)
- Applied Pulse Programming
- Introduction to PVM Programming
- Pipeline Filters

3. Certificate



Programming in *ParaVision*®


ParaVision Programming Course
Ettlingen, April 3-7, 2007
Franek Hennel




ParaVision Overview




Parameter editor





Spectrometer control



Scan control



Acquisition




Geometry editor



ParaVision Overview, cont.



• Workflow



Openness of ParaVision



The user can access:

- Timing control
- Parameters in memory
- Parameter files
- Raw data files
- Image files
- Data flow




Timing control



- Pulse program

```
d0 fq1:f1
d8 gatepulse 1
p1:sp1 ph0
d1 grad{ (-50) | (0) | (0) }
d2
d1 groff
d1 grad{ (50) | (0) | (0) }
50u REC_ENABLE
ADC_INIT(ph0, ph1)
d3 ADC_START
d1 groff
```



Pulse Program Commands



- Fixed delays
- Variable delays
- Incrementable delays
- Fixed RF pulses
- Variable length pulses
- Pulse shapes
- Power lists (var. flip angles)
- Gradient amplitudes
- Gradient ramps
- Gradient shapes
- Gradient functions (steps)
- Frequency lists
- Phase lists

...



Access to parameters

BRUKER

Parameter memory can be accessed from:

- Shell command line
- Shell script (marco)
- C-code
 - Automations
 - Methods

BRUKER
BIOSPIN

ParaVision Parameter Properties

BRUKER

Parameters have

- Types (int, double, user-defined, ...)
- Attributes (visibility, etc.)
- Dynamic array sizes
- Relations (reaction procedures)

Acquisition dimension	3
Dimension Description	[0-2] < 3 Spatial Spatial Spatial
acquisition size	[0-2] < 3 [0-2] [0-2] [0-2]

BRUKER
BIOSPIN

Access to parameters from shell



- Command line

```
/home/xy> pvcmd -set pvScan NR 100
```

- Shell script

```
#!/bin/bash
path=`pvcmd -a pvScan -r pvDsetPath -path PROCNO`
chmod -R u+w $path
$XWINNMRHOME/prog/bin/freco -r -n -p -d $path
pvcmd -a pvScan pvDsetSetScanStatus
chmod ugo-w $path/reco
```



Macros



- Shell scripts placed in
.../curdir/<user>/ParaVision/macros
are called MACROS.
- Macro Manager's role:
 - Executing macros
 - Recording macros
 - Copying, Deleting
 - Editing
 - Output viewing



Macro Manager



```
# Definition of test functions
function PvMacroMgrR3FuncAddCategory {
    echo -n Checking $FUNCNAME ...
    pvcmd -a JMacroManager -r JMAddCategory -category testCategory -dir /tap
    CATEGORIES=`pvcmd -a JMacroManager -r JMListCategories`
    NUM_FOUND=`echo $CATEGORIES | grep testCategory -c`
    if [ "$NUM_FOUND" -eq 1 ]; then
        printResult $OK
    else

```


- Other tools:
 - Command Observer
 - Command Sender



Macro examples



- Any sequence of PV commands can be programmed
- Graphic interfaces may be added with Tcl-Tk



Automations




- Automations are simple c-programs which:
 - Use simplified c-syntax based on c-macros
 - Get parameter values from memory
 - Set parameter values in memory
 - Send PV commands
- Simmilar funtion to macros
- Advantages:
 - Features of c-language, e.g. maths,
 - Faster access to parameters, esp. arrays



Pipeline Filters




- Pipeline Filters are special Automations which:
 - Access the data on-line
 - Modify data contents and size



Pipeline Filters: Applications

BRUKER

- Data modification:
 - Offset demodulation,
 - Resampling
 - Non-standard sorting
- Discarding data segments
 - Navigators
- Non-standard accumulation
 - K-space weighting
- Synthesising data
 - "Key-hole" strategies




BRUKER
BIOSPIN

ParaVision Methods (PVM)

BRUKER


- High-level Interface to Acquisition Parameters



Parameter
Editors



EXPERIMENT



BRUKER
BIOSPIN

PVM Features



- Methods are programmed in C
- Independent sources
- Independent binary (one file/method)
- Method-specific parameters can be defined
- Global parameters can be used
- Mechanism for parameter groups
- Pre-defined modules (e.g. DTI)



What a method does (example)



- Introduce "Echo Time"
- Include Echo Time in protocol
- React to protocol loading
 - Echo Time is legal ?
- React when Echo Time is changed
 - Restrict it,
 - Adapt other parameters (TR, etc)
 - Derive delays d1, d2
- React when other parameters change
 - New TR agrees with TE?

} "Relations"



Macros, , PVM?

BRUKER

- Pipeline filters are special. But macros, automations, methods – they all read and set parameters. What is actually the difference?
- When to write what?

BRUKER
BIOSPIN

- Batch execution of PV commands
 - “I want to run fifty scans overnight”
- Graphic interfaces (Tcl-Tk)
 - “I need a button for each experiment step”
- Starting external software from PV
 - “I want to run my own reco on a selected scan”
- Automatic post-processing
 - “When the scan is ready, I want to start ISA”
- Small modifications at base-level
 - “I am fine with FLASH, but need trim 7 to be 50%”

↳ test an idea before writing a method

BRUKER
BIOSPIN

When to write automations



- Small modifications at base-level involving maths or arrays
 - "I want to use FLASH with this powerlist:"
$$\text{Power}[n] = 20 \log (\cos^{-1}(\exp(-(d1/d2))/n)), n=1...1024$$
- Processing with access to several parameters
 - "Reading parameters with *pvcmd* is too slow for me"



When to write methods



- Any experiment with routine interface
 - "I have my macro for modifying FLASH, but each time I change geometry, my changes are lost"
 - "I want to implement my sequence so that EVERYBODY can use it"
- Projects based on PVM modules
 - "Bruker EPI readout is not bad, but I want to use it with my own excitation sequence"
- Adjustments
 - "I want the frequency to be adjusted for the central voxel for each scan"

?? ->



What we will learn



- **Monday, Tuesday:**
 - Base-level parameters and pulse programming
 - First meeting with *gre*
- **Wednesday**
 - PVM programming, part 1 (*gre* mutations)
 - Parameter definitions
 - Relations
 - Timing
 - RF pulses
- **Thursday**
 - PVM programming, part 2
 - Using modules
 - Creating modules
 - Accelerated MRI
- **Friday**
 - Pipeline filters





Pulse Program Overview

ParaVision Programming Course

April 3-7, 2006

Paul Freitag



Overview



- **Pulse Programs and Parameters**
 - General concepts
 - First Example
 - Control Parameters
- **Introduction To Pulse Programming (1)**
 - Development Tools
 - Project 1: Hello World
 - Timing Details / Transmitter Control
 - Dynamic Changes
 - Project 2: Spin Echo
- **Introduction To Pulse Programming (2)**
 - Gradient Control – Transformations
 - Project 3: Gradients and Parameters
 - Sampling Control
 - Project 4: Phase Cycling and Digitizer Control
 - Gradient Pulses
 - Project 5: Gradients and Pulse Frequency
 - Gradient Control – Ramp
 - Project 6: Gradient Cycling



Objectives


BRUKER

- understand the meaning of pulse and gradient programs
- know tools and conventions for pulse programming
- know to write simple pulse programs from scratch
- be able to understand pulse programs delivered with Bruker methods
- be able to solve problems occurring during pulse program development

BRUKER
BIOSPIN

PPG Concept

BRUKER




BRUKER
BIOSPIN


ACQP: Baselevel parameters



- Sizes/repetitions
- Timing
- Pulses
- Bandwidth
- Digitizer Control
- Gradient transformation
- Gradient trim values
- Gradient ramps
-




Pulse Program Overview



Pulse Program Example

BRUKER



```
#include "Avance.incl"
#include "DBX.include"
```

Standard Header Files

```
start, d0 fq1:f1
      (p0:sp0 ph0):f1
```

Pulse Generation

```
depa REC_ENABLE
      ADC_INIT(ph0,ph1)
deosc ADC_START
      10u ADC_END
```

Data Collection

diff cent in P13.0. → SETUP_GOTO(start)

Flow of Control

```
exit
```

Phase Definition

```
ph1=0
ph0=0
```

BRUKER
BIOSPIN

General Info

BRUKER

- ASCII Text – use arbitrary editor
- Default definition in include files:
need #include "Avance.incl" and #include "DBX.include"
- Commands are linked to durations (delays) –
* implicitly or explicitly
- Durations can be specified as parametrized delays or constants
- Constant delays are floating point numbers with time unit s(econd), m(illi), u(icro)
^ do, dt, aqg ^ by "10u"
- Semicolons introduce comments

BRUKER
BIOSPIN

RF-Pulses



- RF Pulses are generated by a pulse delay
- Parametrized (p0-p63) or constant pulse delays (e.g. 10up) can be used
- A pulse shape can be specified :sp0 - :sp31
- Frequency (SFO1), amplitude (PL[1]) and phase will be set independently for the synthesizer channel.
- Synthesizer channels (:f1) are linked to nuclei parameters



Gradients



- Gradients described in terms of subject oriented coordinate system
- Transformation matrix to change the image geometry
- Individual control of gradient ramps
- User specific gradient shapes
- Phase encoding by variable gradient functions
- Independent



Data Collection

BRUKER

- Standard sequence for data collection
depa REC_ENABLE
ADC_INIT(NOPH,ph1)
deosc ADC_START
10u ADC_END
- Acquisition delay must be slightly longer than
dwell * number of points
- Use synthesizer phase setting ph1
- Phase list must be defined at end of pulse
sequence

BRUKER
BIOSPIN

Parameter Environment

BRUKER

- Run in endless loop for setup mode (gsp):
start, d0
SETUP_GOTO(start) (goto start in Pv 3!)
exit
- In Scan Mode (gop), experiment end is linked
to acquired data (ParaVision <= 3) or loop
structure (ParaVision 4).
- Parameters control number of acquired scans
- NA (number of averaged scans)
- NI (number of data objects)

BRUKER
BIOSPIN



Introduction to Pulse Programming (1)

ParaVision Programming Course
April 3-7, 2006
Paul Freitag



Overview



- Pulse Programs and Parameters
 - General concepts
 - First Example
 - Control Parameters
- **Introduction To Pulse Programming (1)**
 - **Development Tools**
 - **Project 1: Hello World**
 - **Timing Details / Transmitter Control**
 - **Dynamic Changes**
 - **Project 2: Spin Echo**
- Introduction To Pulse Programming (2)
 - Gradient Control – Transformations
 - Project 3: Gradients and Parameters
 - Sampling Control
 - Project 4: Phase Cycling and Digitizer Control
 - Gradient Pulses
 - Project 5: Gradients and Pulse Frequency
 - Gradient Control – Ramp
 - Project 6: Gradient Cycling



Tools for Pulse Programming

BRUKER

- Editor – adapt UXNMR_editor environment (Pv3: setres tool in XWINNMR)
- Edit > Pulsprogram
- Tools > Pulsprogram Tool *only in Pv4.*
(Pv3 Display)

Interactive editor with syntax check and graphical pulse sequence display

BRUKER
BIOSPIN

Starting from Scratch

BRUKER

*try this at home.
Pv4 specific?*

- New Scan / Cancel during new protocol
- data set status new
- no method defined
- parameters initialised for onepulse
- gradient parameters have no value
- multiple slice packages cannot be defined
- The ACQP parameter PULPROG must be set to the name of the pulse program
- PULPROG parameter defines pulse program in exp/stan/nmr/lists/pp

BRUKER
BIOSPIN

/opt/PPC/PPG/basic PPG


/opt/PV4.0/exp/stm/nmr/lists/pp

Project 1: Hello World

BRUKER

By developing this project, you will learn:

- to write a minimum pulse program from scratch
- to use ParaVision tools for pulse program development
- which parameters need to be set on a new data set
- which parameters implicitly have influence on data acquisition



Processing the Text

BRUKER

- line orientation: \n separate lines
- ; comment lines with semicolon
- significant length of identifiers limited: 7 valid characters
- Standard Headers Avance.incl and DBX.include must be included for imaging applications
- include files : use C preprocessor
- special value PULPROG=*pulseprogram* uses expanded pulse program stored with dataset

↓
guarant

BRUKER
BIOSPIN

PULPROG example

BRUKER

PULPROG = simple.ppg

expl/stan/nmr/lists/pp/simple.ppg:

```
#include <deosc.mod>
start, d0
  p0
  aq
goto start
exit
```

expl/stan/nmr/lists/pp/deosc.mod:

```
define delay deosc = { $DEOSC }
"deosc = abs(deosc) * 1e-6"
```

uses pulse program stored in dataset.
PULPROG = pulseprogram

data/<usr>/nmr/<id>/1/pulseprogram:

```
# 1 "/u/exp/stan/nmr/lists/pp/simple.ppg"
# 1 "/u/exp/stan/nmr/lists/pp/simple.ppg"
# 1 "/u/exp/stan/nmr/lists/pp/deosc.mod" 1
define delay deosc = { $DEOSC }
deosc = abs(deosc) * 1e-6"
# 1 "/u/exp/stan/nmr/lists/pp/simple.ppg" 2
start, d0
  p0
  aq
goto start
exit
```

BRUKER
BIOSPIN

Module Concept

BRUKER

- *#include <modulehead>*
- *#include <modulebody>*
- *if (ParaVisionParameterCondition)*
 {
 else
 }
- *#define ADC_START* : C-Preprocessor
 macro definition possible

def'n is in dbx/include directory

BRUKER
BIOSPIN

Loop structure: Scanner control



- Number of acquired scans determined by loop parameters:
- RCU based spectrometers:
acquisition is aborted, as soon as the expected number of scans has been acquired
- Note display options for accumulation:
 - each scan
 - each accumulation
 - each phase encoding step



Loop parameters



ie. ways to accumulate the data

- NS: number of simply accumulated scans *(same the phase encoding)*
- ACQ_ns_list_size / ACQ_ns_list:
description of echo accumulation scheme
- ACQ_phase_factor: connected scans
- NI: number of image objects *number of slices or numbo. of echoes in MSlice*
- ACQ_dim / ACQ_size : description of image object
- NA: phase encoding accumulation
- NAE: image accumulation
- NR: number of repetitions




Loop structure Overview

BRUKER



- 3 NS accumulation
- 2 ACQ_phase_factor
- 3 NI Image objects / NSLICES
- 2 NA accumulation
- 3 ACQ_size[1] / ACQ_phase_factor
- ACQ_size[2]
- NR repetition



BRUKER
BIOSPIN

Timing concepts

BRUKER

- timing units
s(seconds), m(milli), u(micro)
- resolution 80MHz, 0.125u
- explicit timing <-> implicit timing
- actions start at beginning of delay
- parallel timing possible (use parentheses)
- multiple control commands possible within same delay

BRUKER
BIOSPIN

Timing parameters



- D[0-31] : general purpose delays/unit [s]
- P[0-31] : general purpose pulses/unit [u]
- VDLIST: variable delay list/unit [s]
- VPLIST: variable pulse list/unit[u]
- DE: minimum delay/unit [u] → ringing delay?
- DEOSC: minimum acquisition delay[u]
- DEPA: minimum preamp delay[u]




Timing commands



- $\langle \text{float} \rangle u, m, s$ fixed delay
- $\langle \text{float} \rangle up, mp, sp$ fixed pulse delay
- $d0-d31$ parametrized delay
- $p0-p31$ param. pulse delay
- vd variable delay
- vp variable pulse delay
- de min. ringing delay
- dw dwell time (real)
- $depa$ min. delay
- $\langle \text{delay} \rangle * \langle \text{scaling factor} \rangle$ scaled delay



Timing example



(10u) (20u) ; parallel pulse trains in parentheses

d0


p1*2

vd

deosc



Synthesizer Control Overview



Channel concept



- Eight transmit channels can be addressed
- In pulse program, channel id :f1 - :f8
- Basic frequency defined by nucleus
- Amplifier mapping part of routing
- default channel: f1



Transmitter Control Parameters



- TPQQ[0-15]/DPQQ[0-15]
 - .name: shape file in exp/stan/nmr/lists/wave
 - .power: attenuation in dB from -6dB to 120dB
 - .offset: phase shift generated offset in Hz
- SFO1-SFO8: Basic frequency in Hz
- ACQ_O1_list_size/ACQ_O1_list-
ACQ_O3_list – ACQ_O1B_list
 - offset frequency lists – B for second oscillator
- PL[0-31]: default power levels for channels
PL[1] for :f1 ...



Transmitter Control Commands

BRUKER

- *fq1:f1-fq8:f1*
set frequency from ACQ_O1_list ...,
implicit increment
- *fq1b receive* : Use second DDS for receive
- *fq8b:f1* set frequency from ACQ_O1B_list
- *pl0:f1-pl31:f1*:
set power level from PL[0]-PL[31]
- *:sp0-:sp31 [(currentpower)]*
use shape from TP(DP)QQ[0]-TP(DP)QQ[15],
use current power as option

BRUKER
BIOSPIN

Phase List Definition

BRUKER

- *ph0:f1-ph31:f1*
set phase from ph0-ph31
- *ph0-ph31=(<unit>)* list definition
- Shortcuts for repeated / shifted groups
 - *ph0= { < list function> } * <factor>*
 - *ph1= { < list function> } ^ {increment}*
- phase values as multiples of the fractional unit, default 4
- receiver (mixing) phase only fractional unit 4 available

BRUKER
BIOSPIN

Transmitter Control Details

BRUKER

*refers to
↓ the freq. channel*

- gatepulse 1: amplifier presetting
- preset off: optimize performance *(if not included, may have timing probs.)*
- shape timing: 0.1u grid
- shape files exp/stan/nmr/lists/vw
- dynamic changes limited
- setup control: linked to sliders
- TPQQ[0-7] : sp0-sp7
- DPQQ[0-7] : sp8-sp15
- TPQQ[8-15] : sp16-sp23
- DPQQ[8-15] : sp24-sp31

BRUKER
BIOSPIN

Pulse Power Setup

BRUKER

TPQQ and DPQQ connected with sliders in Spectrometer Control Tool:

The screenshot shows two overlapping windows from the Bruker software. The 'Transmitter Pulse' window on the left has tabs for 'Transmit Pulse' and 'Decoupling Pulse'. The 'Frequency Setting Control' window on the right shows parameters for 'Basic Frequency', 'Offset Freq', and 'Dc Attenuator 0' and 'Dc Attenuator 1'. The 'Dc Attenuator 0' is set to 50.0 and 'Dc Attenuator 1' is set to 0.0. There are also checkboxes for 'Lock' and 'Independent Scale'.

TPQQ: Transmit Tx
corresponds to channel 1
DPQQ: Decoupling Dc
corresponds to channel 2

BRUKER
BIOSPIN

Transmitter Control Example

BRUKER

...

10u pl1:f1

10u fq1:f1

d8 gatepulse 1

(p0:sp0(currentpower) ph1):f1

...

ph1 = {{0 1}*2}^1 ; 0 1 0 1 1 2 1 2

ph2 = (360) 0 90 180 270

BRUKER
BIOSPIN

Changing the Flow of Control

BRUKER

- `<label>`, mark position in pulse program
- `goto <label>`: jump to label
- `lo to <label> times <counter>` : repeat block times counter
- counter may be ParaVision parameter or I0-I31 corresponding to L[0]-L[31]
- `if "<condition>" goto <label>`
- `exit` : end execution

BRUKER
BIOSPIN

Increment, Decrement and Reset



- *ipp0-ipp31, rpp0-rpp31*: increment/reset phase list pointer
- *ip0-ip31, dp0-dp31, rp0-rp31*: increment/decrement/reset phase list
- *ipu0-ipu31, dpu0-dpu31, rpu0-rpu31*: increment,decrement,reset pulse delay by value defined in INP[0]-INP[31]
- *id0-id31, dd0-dd31,rd0-rd31*: increment, decrement, reset delay by value defined in IND[0]-IND[31]
- *ivd/ivp*: increment pulse/delay list pointer



Calculating in Pulse Programs



- "expression": calculation in double quotes
- C-style syntax for assignment and simple algebra and functions
- time units are preserved
- ParaVision parameters cannot be evaluated directly
- examples:
 - "d0 = d1 * 2 - p0"
 - "p0 < 100u"



User defined Commands

BRUKER

- New commands can be defined in the beginning of the pulse program (before the first delay) ✓ 7 characters only!
 - *define delay* delayname
 - *define pulse* pulsename
 - *define loopcounter* loopcountername
 - *define list*<delay> delaylistname
 - *define list*<pulse> pulselistname
 - *define list*<power> powerlistname
 - *define list*<frequency> fqlistname

BRUKER
BIOSPIN

Initialising User Commands

BRUKER

- User command lists must be initialized
 - either directly:
define list ... = { values }
 - from file:
define list ... = <file>
 - from ParaVision parameter
define = {\$ParaVisionParameterName}

BRUKER
BIOSPIN

Using and Navigating in Lists




- User defined list commands can be used like the built in commands:
 - *define list<listtype> somelist: definition*
 - *somelist*: use list at current position
 - *somelist.inc*: increment list pointer
 - *somelist.dec*: decrement list pointer
 - *somelist.res*: reset list pointer
 - *somelist[ndx]*: access list value at ndx
 - *somelist.idx*: access list pointer in relation
 - *somelist^*: use list with autoincrement



Pulse Power List Example



```
define list<power> plist = { 120 60 30 0 }  
start, 10u plist:f1 plist.inc  
  10up  
  5u  
  lo to start times 4  
exit
```




Project 2: Spin Echo



By developing this project you will learn

- to use parameters for timing
- to implement easy relations between the parameters in the pulse program
- to vary properties of rf-pulses during runtime



Introduction to Pulse Programming (2)

ParaVision Programming Course



April 3-7, 2006

Paul Freitag

Overview

- **Pulse Programs and Parameters**
 - General concepts
 - First Example
 - Control Parameters
- **Introduction To Pulse Programming (1)**
 - Development Tools
 - Project 1: Hello World
 - Timing Details / Transmitter Control
 - Dynamic Changes
 - Project 2: Spin Echo
- **Introduction To Pulse Programming (2)**
 - **Gradient Control – Transformations**
 - **Project 3: Gradients and Parameters**
 - **Sampling Control**
 - **Project 4: Phase Cycling and Digitizer Control**
 - **Gradient Pulses**
 - **Project 5: Gradients and Pulse Frequency**
 - **Gradient Control – Ramp**
 - **Project 6: Gradient Cycling**

12000
41496 Hz/cm



Gradient Control Parameters



- ACQ_patient_pos:
normalized subject coordinate system
- NSLICES: number of slices
- ACQ_grad_matrix[][][]
transformation matrix per slice
- ACQ_trim: trim values
- ACQ_read_scale, ..._phase_..., slice :
Scaling values in logical coordinates
- ACQ_phase_encoding_mode: ramp mode
- ACQ_phase_enc_start: ramp shift



Gradient Coordinate Transformation



$$\text{grad}\{(80,80,100) | (50,50,50) | (0,0,100)\}$$

	read			phase			slice		
	x	y	z	x	y	z	x	y	z
ACQ_scaling_read/..									
{ 1,1,0.5}									
ACQ_grad_matrix									
$\begin{vmatrix} 0.7 & 0.7 & 0 \\ -0.7 & 0.7 & 0 \\ 0 & 0 & -1 \end{vmatrix}$									
ACQ_patient_pos									
head_left									
PREEMP_X/Y/Z_factor									
{ 0.98, 0.98, 0.94 }									


$(80,80,100) \cdot (0.7, 0.7, 0) = (56, 56, 0)$
 $(50,50,50) \cdot (-0.7, 0.7, 0) = (-35, 35, 0)$
 $(0,0,100) \cdot (0, 0, -0.5) = (0, 0, -50)$
 $(56,56, 0) + (-35, 35, 0) + (0,0,-50) = (21,91,-50)$
 $(21,91,-50) \cdot z^*-1, xy \rightarrow yx = (91,21,50)$
 $(91,21,50) \cdot (0.98, 0.98, 0.94) = (89.18, 20.58, 47)$



Gradient Calibration

BRUKER

- PREEMP_ramp_time : technical ramp time (amplifier input ramp)
- PREEMP_ramp_mode : ramp description
- PREEMP_rise_time : physical ramp time
- PREEMP_X_factor/Y_factor/Z_factor: calibration constants



BRUKER
BIOSPIN


Gradient Commands

BRUKER

- *grad*{ <read> | <phase> | <slice> }
 - <read>/<phase>/<slice> gradient strength as percentage:
 - (0) : single trim value, range -100 to +100
 - (100,100,100) : x/y/z trim value
 - (t0)-(t99) : parameter trim value
 - <>+*trim : expression of trim values
- *groff*: shortcut for *grad*{(0)|(0)|(0)}
- delay marks **start** of gradient ramp
- effective ramp delay PREEMP_ramp_time

BRUKER
BIOSPIN

Gradient command examples



```

ramp*0.5 grad { (20) | (0) | (0) }
ramp*1.5      ; ACQ_trim[0] = {-30, -30, -30 }
ramp      grad { (t0) | (0) | (-20,-20,-20) }
ramp*3
ramp*3 groff
    
```




Project 3: First Gradients




By developing this project you will learn

- which parameters must be set to use gradients within a pulse program
- how geometric transformations influence gradients
- how static gradients can be switched in a pulse program



Data Sampling Overview

BRUKER



Acquisition Commands

BRUKER



- **REC_ENABLE:**
 - receiver enable for RX22 - at least depa
- **ADC_INIT(_B)(<mixing>, <reference>)**
 - prepare acquisition (use fq8b frequency list)
 - set mixing phase and reference phase
 - implicit delay DE - DEPA
- **ADC_START:**
 - mark acquisition of first point
 - must wait DEOSC until ADC_END
- **ADC_END:**
 - mark end of acquisition/minimum delay 10u

BRUKER BIOSPIN

Receiver Reference vs Mixing Phase

BRUKER

- ADC_INIT_B(mixing phase , reference phase)
- reference phase -> synthesizer phase
- mixing phase = channel mixing
(u,v) -> ($\pm u, \pm v$)
- attention: mixing phase prevents fast loops
- Dual DDS: TX -> RX phase continuous
RX -> TX phase coherent



BRUKER
BIOSPIN

Scan Parameters

BRUKER

- ACQ_size[0] : size of scan (real valued) alias TD
- ACQ_scan_size: combination flag
- ACQ_phase_factor: scan division
- ACQ_scan_shift: group delay compensation
- SW_h: bandwidth $dw = SW_h / TD$
- DIGTYP: digitizer selection
- DIGMOD: filter mode analog/digital
- AQ_mod: data sampling mode qdig/qsim

no group delay related to digital filter

BRUKER
BIOSPIN

Acquisition Timing


BRUKER

- Parameter Setting:
 - ACQ_size[0] = 1024
 - SW_h = 50000 => dw = 20 aq = 10240u
 - DSPFIRM = DSP_medium
 - Derive minimum delay for optimum shift
 - ACQ_scan_shift = -1 -> -26
 - DEOSC = -1 -> -10560
 - Derive optimum shift for available delay
 - DEOSC = +10400
- ACQ_scan_shift = -10

BRUKER
BIOSPIN

Acquisition Start Details

BRUKER



BRUKER
BIOSPIN

Acquisition Start Example



synthesizer b.
fq1b receive ; use ACQ_O1B_list for *channel 1*
...
10u fq8b:f1 ; set and increment ACQ_O1B_list
...
depa REC_ENABLE
ADC_INIT_B(NOPH,ph0)
deosc ADC_START
10u ADC_END
; wait at least 800u until next ADC_INIT




Project 4: Phase Cycling and Digitizer Control



By developing this project you will learn

- to implement minimum timing for data acquisition
- to select a frequency offset for acquisition
- to run phase cycling



Slice Increment



- One ACQ_grad_matrix for each slice
- Dimension cannot be changed within acqp
- Number of orientations in NSLICES
- change orientation
 - *islice* : next orientation matrix
 - *rslice* : previous orientation matrix
 - *zslice* : first orientation matrix
 - *sslice* : save slice index
 - *rslice* : restore slice index




Gradient Shapes



- Ramp functions with length definition
- `grad{ ... <function>(<trim>,<len>) | ... }`
- length must match function definition
- user defined functions possible
- delay defines step resolution


```
ufun = { 0.5, 0.8, 1.0, 0.8, 0.6, 0.4, 0.2, 0.0 } ...  
d4 grad{(ufun(100,8)) | (0) | (0) }
```



Gradient Shape Example



```
1m
1s grad{ (0) | sinp(100,30) | cos(100,30) }
1s grad{ r1d(100,10) | (0) | (0) }
1m groff
exit
```



Gradient Sequences



- Gradient sequence in parallel to pulse timing can be defined
- `gc_control {`
 - `<delay> grad{...|...|...}`
 - `<delay> grad{...|...|...}`
 - ...
- `}`




gc_control Example



```

1m gc_control {
  10m grad { (10) | (0) | (0) }
  10m grad { (0) | (10) | (0) }
  10m grad { (0) | (20) | (0) }
  10m groff
}
15m ; pulse timing independent
10mp ; from gradient timing
15m
exit
  
```

*this parallel gradient
const timing may not be displayed correctly!*




Project 5: Gradients and Pulse Frequency



By developing this project you will learn

- To use a slice selection gradient
- To vary the transmitter frequency for pulses



Phase encoding



- ramp list functions *r1d*, *r2d*, *r3d*
- linked to parameters
ACQ_phase_encoding_mode/enc_start
- further built in functions
 - *sin*
 - *cos*
 - *sinp*
 - *gauss*<truncvalue>
 - *step*
- *lgrad* <function>[<2d>|<3d>] = length
length definition/setup mode




Gradient Ramp Increment



- step through gradient ramp lists
 - *igrad* <function> : increment pointer
 - *dgrad* <function> : decrement pointer
 - *zgrad* <function> : reset pointer
 - *sgrad* <function> : save status
 - *rgrad* <function> : restore status



Gradient Encoding Example: Centered




- ACQ_phase_encoding_mode[1] = Centered
- ACQ_trim[0] = { 30, 30, 30 }
- ACQ_size[1] = 8
- ACQ_phase_enc_start[1] = -1

```
igrad r2d = ACQ_size[1]  
ramp grad{ (0) | r2d(t0) | (0) }  
igrad r2d
```



Gradient Encoding Example: RARE



```
slice, rgrad r2d  
cloop, d1 grad((0) | r2d(t0) | (0))  
groff  
igrad r2d  
lo to cloop times ACQ_phase_factor  
sgrad r2d  
lo to slice time NSLICES  
lo to start times ACQ_size[1]/ACQ_phase_factor
```




Project 6: Gradient Cycling



By developing this project you will learn

- to vary gradients during runtime
- to use parameters to control gradient variation
- to switch off gradients during setup mode







Pulse Programming Exercises

*ParaVision Programming Course
April 3-7, 2006*

Author:
Paul Freitag

Bruker BioSpin MRI GmbH



Table of Contents

1	Starting from scratch	5
1.1	Objective	5
1.2	Project description	5
1.3	Extensions	6
1.4	Hints	6
1.5	Solution	7
2	Pulses, Parameters and Relations	8
2.1	Objective	8
2.2	Project	8
2.3	Extensions	9
2.4	Hints	9
2.5	Solution	10
3	First Gradients	12
3.1	Objective	12
3.2	Project	12
3.3	Extensions	13
3.4	Hints	13
3.5	Solution	14
4	Fast Data Acquisition and Phase Cycling	16
4.1	Objectives	16
4.2	Project	16
4.3	Extensions	17
4.4	Hints	17
4.5	Solution	18
5	Gradients and Pulse Frequency	20
5.1	Objectives	20
5.2	Project	20
5.3	Extensions	21
5.4	Hints	21

5.5	Solution	22
6	Gradient Cycling	24
6.1	Objectives	24
6.2	Project	24
6.3	Extensions	25
6.4	Hints	25
6.5	Solution	26

1 Starting from scratch

1.1 Objective


By developing this project, you will learn:

- to rename an easy pulse program for your data set
- to use ParaVision tools for pulse program development
- which parameters need to be set on a new data set

1.2 Project description

In this project, you will create your first pulse program of your own.

For this purpose, you create a new data set without loading a protocol, you create a copy of the pulse program project1.ppg provided in your home directory to the pulse program library directory. You use the parameter editor to select the pulse program for your newly created data set. You use the pulse program display to inspect the pulse program and to correct a syntax error (missing comment sign) in the interactive editor. You set suitable acquisition parameters and run a first pulse-acquire experiment on your spectrometer. Make a list of parameters for which you think that they influence the acquisition. You have finished as soon as you have acquired a reasonable fid file.



1.3 Extensions

If you finish early, you may try to minimize your pulse program, i.e. try to remove as many lines in the pulse program as possible to allow you still to acquire an fid.

1.4 Hints

- Create a new data set (new scan in pvScan), using the cancel button, when asked to load a protocol.
- Pulse programs must be installed in the <PvInstDir>/exp/stan/nmr/lists/pp directory.
- Parameters to be set: Pulse Program Name (PULPROG), Relaxation delay D[0] Pulse delay P[0], Pulse Shape 0 TPQQ[0] .name, TPQQ[0] .power. A suitable pulse will be the blockpulse bp .exc.
- In order to simplify parameter input, you can open a parameter editor restricted to the parameters relevant to you:

```
pvcmd -a pvScan pvEditPars PULPROG D[0] P[0] TPQQ[0].name
TPQQ[0].power
```
- Further parameters have influence on the acquisition: The bandwidth (SW_h), the group delay compensation (ACQ_scan_shift), the amplifier gating delay D[8] (derived from CONFIG_amplifier_enable).

\uparrow
 TR

making a simple macro:

Macro Manager → record → stop (input name of macro → edit file.

1.5 Solution

```
;
; PPC 2006 Pulse Programming Course
;
; Project 1: first pulse program - simplified SINGLEPULSE
;
; Explicit parameters
;
; P0 - excitation pulse delay
; d0 - relaxation delay
; d8 - CONFIG_amplifier_enable per convention
;

;-----necessary include files-----
#include<Avance.incl>
#include <DBX.include>

; automatic pulse gating switched off for fast sequences
preset off

;label--delays--commands----gradients-----remarks-----
start, d0                                ; relaxation delay

;-----excitation-----
d8      gatepulse 1                       ; pulse gating
(p0:sp0 ph0):f1                          ; pulse TPQQ[0]

;-----data acquisition-----
depa    REC_ENABLE                       ; open receiver
ADC_INIT(ph0, ph1)                       ; prepare acquisition
Aqq     ADC_START                         ; start data sampling
5m      ; group delay comp.
10u     ADC_END                           ; end data sampling

;-----data loops-----
SETUP_GOTO(start)

exit

;-----phase definitions-----
ph0 = 0                                   ; excitation phase
ph1 = 0                                   ; reference phase
```

2 Pulses, Parameters and Relations


2.1 Objective

By developing this project you will learn

- to use parameters for timing
- to implement easy relations between the parameters in the pulse program
- to vary properties of delays during runtime

2.2 Project

In this project, we want to implement an easy spin-echo 1D method, which can be controlled by base level parameters. By acquiring a series of scans with varying echo time, the method can be used to observe T2 decay.



A spin echo method consists of an excitation pulse, a preparation delay an inversion pulse and the detection period. The spin echo will appear after the same time as the preparation delay. In order to observe the echo, the flip angle of the inversion pulse must be 180deg and the position of the acquisition delay must be adjusted such that the echo will occur in the middle of the acquisition delay. By varying the delay between the excitation and the refocusing pulse, the echo position can be varied.

To achieve this, create a new scan and create a copy of the pulse program from the first project to `project2.ppg`. Add a new pulse into the pulse program controlled by the parameters `P[1]` and `TPQQ[1]` and a phase list `ph2`. Use the delay `D[1]` to describe the echo time. Calculate filling delays in relations in the pulse program in order to position the echo correctly in the acquisition window. Use this first implementation of the pulse sequence to set up suitable parameters for the refocusing pulse and the echo time. Once you observe a reasonable echo, add a loop over the parameter `NI` to the pulse program. Add a variable

delay between the two pulses and the pulse and the acquisition window and vary the delay according to a list in say 6 steps of 50ms each.

2.3 Extensions

Use a different method, to vary the echo time. Possible methods are:

- `vd` delay and `ACQ_vd_list` parameter with `ivd`
- `define list<delay> vardel = { 0.05 0.1 0.15 0.2 0.25 0.3 }`
with `vardel.inc`
- `id1` command in combination with the `IN[1]` parameter
- Incrementing `D[1]` in a "parameter relation"

2.4 Hints

- Clone the first scan
- Copy the pulse program `project1.ppg` to the new pulses program `project2.ppg`
- Remember that the pulse should consist of the two lines
`d8 gatepulse 1`
`(p1:sp1 ph2):f1`
- Add the phase list `ph2=0` to the end of the pulse program – think which phase relative to the excitation pulse is reasonable.
- To calculate the correct relations you may define new delays
`define delay drest`
- To calculate the correct echo delay, you must take into account, that the `ADC_INIT` command has an implicit delay `de-depa`.
- When the echo delay `D[1]` is introduced between the refocusing pulse and the acquisition delay, a suitable filling delay between excitation and refocusing pulse might calculate
`"drest = d1 - de + depa + aqq/2 - p0/2 - d8"`
- You can start an adapted parameter editor
`pvcmd -a pvScan pvEditPars PULPROG NI D[0] D[1] P[0]`
`TPQQ[0].name TOQQ[0].power P[1] TPQQ[1].name TOQQ[1].power`
- Add the loop
`lo` to start times `NI`
just before the `SETUP_GOTO(start)` macro
- When adding a delay list, don't forget to increment the list pointer

2.5 Solution

```

;
; PPC 2006 Pulse Programming Course
;
; Project 2: simple spin echo
;
; Explicit parameters
;
; p0 - excitation pulse delay
; p1 - refocusing pulse delay
; d0 - relaxation delay
; d1 - echo time
; d8 - CONFIG_amplifier_enable per convention
; NI - image object loop
;

;-----necessary include files-----
#include<Avance.incl>
#include <DBX.include>

;-----internal parameter definitions and relations-----
define delay dr
"dr = aqg / 2 + dur1 - p0/2 - d8"
define list<delay> vdl = { 0.05 0.10 0.15 0.20 0.25 0.30 }

; automatic pulse gating switched off for fast sequences
preset off

;label--delays--commands---gradients-----remarks-----
start,      d0                                ; relaxation delay

;-----excitation-----
d8          gatepulse 1                       ; pulse gating
(p0:sp0 ph0):f1                                ; pulse TPQQ[0]

;-----echo time padding-----
dr          ; echo time padding
vdl        ; var. echo spacing

;-----inversion-----
d8          gatepulse 1                       ; pulse gating
(p1:sp1 ph2):f1                                ; pulse TPQQ[1]

vdl        ; var. echo spacing

;-----data acquisition-----
depa      REC_ENABLE                          ; open receiver
ADC_INIT(ph0, ph1)                            ; prepare acquisition
aqg      ADC_START                            ; start data sampling
5m       vdl.inc                              ; group delay comp.
10u      ADC_END                              ; end data sampling

;-----data loops-----
lo to start times NI                          ; NI "image" objects

SETUP_GOTO(start)
    
```

exit

```
-----phase definitions-----  
ph0 = 0 ; excitation phase  
ph1 = 0 ; reference phase  
ph2 = 0
```

3 First Gradients

*PREEMP - grad. - const
strength in Hz/cm.*

3.1 Objective


By developing this project you will learn

- which parameters must be set to use gradients within a pulse program
- how geometric transformations influence gradients
- how static gradients can be switched in a pulse program

3.2 Project

In this project, we want to implement a basic gradient echo and learn which base level parameters must be used to control the gradient system.

In order to observe a gradient echo, the spin system must be dephased after excitation and rephased during the acquisition period by applying a gradient of the opposite polarity. The echo will appear when the integral of the rephrasing gradient reaches the same integral as the dephasing gradient.



In order to implement this, create a new scan and create a copy of the pulse program from project 1 to project3.ppg, add a dephasing gradient after the end of the excitation pulse which has half of the length of the acquisition interval. Then switch immediately to the opposite polarity and start data sampling. Switch the gradient off, when data sampling is finished. Use the trim value t_0 to parameterize the gradient strength. Have in mind that the minimum delay for switching gradient states depends on the gradient system and preemphasis setting ($\text{PREEMP_ramp_time}/\text{PREEMP_rise_time}$). Before you can start data acquisition, you must make sure that parameters describing the gradient transformation (ACQ_grad_matrix) have been set to a suitable value.

3.3 Extensions

Use different logical gradient directions to create your echo

3.4 Hints

- Clone the first scan
- Copy the pulse program `project1.ppg` to the new pulse program `project3.ppg`
- A gradient with trim value e.g. in read direction can be applied with the command
`grad { (t0) | (0) | (0) }`
- Negative trim values are not allowed but expressions where a trim value is subtracted from the value (0)
`grad { (0)-(t0) | (0) | (0) }`
- You can start an adapted parameter editor
`pvcmd -a pvScan pvEditPars PULPROG D[0] D[1] P[0]
TPQQ[0].name TOQQ[0].power ACQ_trim[0] ACQ_grad_matrix`
- Import the ramp length from a ParaVision parameter: many methods set `D[4]` to the ramp time – remember that ramp parameters are specified in microseconds.
- Define `delay ramp = { PREEMP_rise_time }`
- `"ramp = ramp * 1e-6"`

3.5 Solution

```

;
; PPC 2006 Pulse Programming Course
;
; Project 3: simple gradient echo
;
; Explicit parameters
;
; p0 - excitation pulse delay
; d0 - relaxation delay
; d8 - CONFIG_amplifier_enable per convention
; t0 - trim value ACQ_trim[0]
; NI - image object loop
;

;-----necessary include files-----
#include<Avance.incl>
#include <DBX.include>

;-----internal parameter definitions and relations-----
define delay ramp = { $PREEMP_rise_time }
"ramp = ramp * 1e-6"
define delay dr
"dr = aqg / 2 "

; automatic pulse gating switched off for fast sequences
preset off

;label--delays--commands---gradients-----remarks-----
start, d0                                     ; relaxation delay

;-----excitation-----
d8 gatepulse 1                               ; pulse gating
(p0:sp0 ph0):f1                             ; pulse TPQQ[0]

;-----dephasing gradient-----
dr      grad { (0)-(t0) | (0) | (0) }
ramp    groff

;-----refocussing gradient-----
ramp    grad { (t0) | (0) | (0) }; Refocussing
;-----data acquisition-----
depa    REC_ENABLE                          ; open receiver
ADC_INIT(ph0, ph1)                          ; prepare acquisition
aqg ADC_START                               ; start data sampling
5m      groff                                ; gradients off
10u     ADC_END                              ; end data sampling


;-----data loops-----
lo to start times NI                         ; NI "image" objects

SETUP_GOTO(start)

exit
    
```

```

;-----phase definitions-----
ph0 = 0 ; excitation phase
ph1 = 0 ; reference phase
    
```



$$d \text{ rest} + \text{ramp} = \text{ramp} + a q q$$

-13000
~15 kHz

4 Fast Data Acquisition and Phase Cycling

4.1 Objectives

By developing this project you will learn


- to implement minimum timing for data acquisition
- to select a frequency offset for acquisition
- to run phase cycling

*fq1b receive
to prepare b synth channel for receiving*

4.2 Project

In this project, you will get acquainted with the details of the data collection and phase cycling process.


To apply a read offset in imaging, the LO frequency of the receiver can be tuned during data acquisition. This will lead to a shift of the acquired data in the frequency domain and can be interpreted as a spatial shift in the image domain.



To implement this project clone the scan from project 3, create a copy of the pulse program to `project4.ppg`. To achieve minimum timing in a pulse sequence, replace the fixed delay following the acquisition delay by a delay of length `DEOSC - aqc`. Relations ensure that `DEOSC` will be set to a minimum value when set to a value less than 0.

To implement the read offset, a second synthesizer channel must be used. This is implied by choosing the macro `ADC_INITB` instead of `ADC_INIT`. Usage of the synthesizer must be Prepared by the `fq1b receive` command at the beginning of the pulse program, actually frequency setting is performed by the `fq8b` command by default linked to the `ACQ_O1B_list` parameter. By setting suitable frequencies you can see the signal move in the reconstruction display.

Phase cycling is often used to remove offsets during averaging. By varying the transmitter phase in 90deg steps and combining the complex data of the scans acquired this way accordingly, offsets related to the receiver hardware will cancel out.



To implement phase cycling, add an averaging loop `lo` to start times `NA`. Use different phase lists for transmit, reference and receiver phase. Make sure that phase pointers are incremented within the averaging loop. Use a phase cycle 0 90 180 270 deg for transmitter and receiver phase and observe that the data adds up, when `NA = 4` is selected. When the receiver phase is kept constant, the signal will cancel out however. Note that the same effect seems to appear for the reference phase but in fact changing the reference phase will not cancel out offset effects, because data combination is not affected.

4.3 Extensions

- Observe the effect of the `ACQ_scan_shift` parameter on the `DEOSC` parameter

4.4 Hints

- Define the delay `deosc` in the pulse program and import it from the ParaVision parameter `DEOSC`. Rescale the value in a relation from microseconds to seconds:

```
define delay deosc = {$DEOSC}
"deosc = abs(deosc)*1e-6 - aqg"
```
- Set the parameter `DEOSC` to a negative value to have the relations calculate the minimum value for `DEOSC`.
- Initialize the parameter `ACQ_O1B_list_size` with 3 and set the values of `ACQ_O1B_list` = { -1000, 2000, 0 }.
- For phase cycling, define three phase lists `ph0`, `ph1`, `ph2` after the `exit` command.
- Use the phase command `ph0` to set the transmitter phase `p0` `ph0: f1`.
- Use the phase command `ph2` for the receiver reference and `ph1` for the receiver mixing phase.
- Don't forget to increment the phase lists in the main acquisition loop using the `ipp0`, `ipp1`, `ipp2` commands.
- Set the parameter `NA = 4` to average data
- Set the parameter `NI = 4` to observe the different phase settings simultaneously.

4.5 Solution

```

;
; PPC 2006 Pulse Programming Course
;
; Project 4: data acquisition details
;
; Explicit parameters
;
; p0 - excitation pulse delay
; d0 - relaxation delay
; d8 - CONFIG_amplifier_enable per convention
; t0 - trim value ACQ_trim[0]
; NI - image object loop
; NA - averaging loop

;-----necessary include files-----
#include<Avance.incl>
#include <DBX.include>

;-----internal parameter definitions and relations-----
define delay deosc = { $DEOSC }
"deosc = abs(deosc) * 1e-6 - aqq"
define delay ramp = { $PREEMP_rise_time }
"ramp = ramp * 1e-6"
define delay dr
"dr = aqq / 2 "

; automatic pulse gating switched off for fast sequences
preset off
; use second synthesizer for LO frequency
fq1b receive

;label--delays--commands---gradients-----remarks-----
start, d0 fq8b:f1 ; select LO freq.

;-----excitation-----
d8 gatepulse 1 ; pulse gating
(p0:sp0 ph0):f1 ; pulse TPQQ[0]

;-----dephasing gradient-----
dr grad { (0)-(t0) | (0) | (0) }
ramp groff

;-----refocussing gradient-----
ramp grad { (t0) | (0) | (0) } ; Refocussing

;-----data acquisition-----
depa REC ENABLE ; open receiver
ADC_INIT_B(ph2, ph1) ; prepare acquisition
aqg ADC_START ; start data sampling
deosc groff ; minimum eosc delay
10u ADC_END ; end data sampling
10u ipp0 ipp1 ipp2


;-----data loops-----
lo to start times NI ; NI "image" objects
    
```

```
lo to start times NA ; averaging loop
    10u      rpp0 rpp1 rpp2 ; reset phases
```

```
SETUP_GOTO(start)
```

```
exit
```

```
;-----phase definitions-----
ph0 = 0 2 1 3 ; excitation phase
ph1 = 0 ; reference phase
ph2 = 0 2 1 3 ; receiver phase
```



5 Gradients and Pulse Frequency

5.1 Objectives


By developing this project you will learn

- to use a slice selection gradient
- to vary the transmitter frequency for pulses

5.2 Project

In this project you extend the gradient echo developed in project 3 by slice selection and by varying the excitation frequency to acquire profiles for different slices.

To implement slice selection, the slice gradient must be applied during the excitation pulse and refocused afterwards. For symmetric excitation pulses the integral of the refocusing gradient must be the inverse of the integral from the middle of the pulse to the end. The profile of the excited slice depends on the gradient strength and the bandwidth of the applied excitation pulse.



To implement the project, clone the previous scan, copy the pulse program to `project5.ppg`. Switch the slice gradient on immediately before the excitation pulse and switch it to the opposite polarity immediately after the pulse. Switch the gradient off after half of the the pulse duration. Use the trim value τ_1 to control the gradient strength

In order to change the transmitter frequency, apply a frequency list at the beginning of the acquisition loop to frequency channel 1. Remember to change the frequency after each excitation. Use the `ACQ_O1_list/ACQ_O1_list_size` parameters to define different offset frequencies and set `NI` to `ACQ_O1_list_size` in order to observe the different slice profiles simultaneously.

5.3 Extensions

Use a different list syntax for frequency switching – remember that the `fq1-fq8` commands include an automatic increment of the acquisition pointer.

5.4 Hints

- You can start an adapted parameter editor

```
pvcmd -a pvScan pvEditPars PULPROG NI D[0] P[0] TPQQ[0].name  
TOQQ[0].power ACQ_n_trim ACQ_trim[0] ACQ_trim[1]  
ACQ_O1_list_size ACQ_O1_list
```
- If you don't use the built-in frequency setting commands `fq1-fq8`, do not forget to switch increment the frequency list at suitable points.

5.5 Solution

```

;
; PPC 2006 Pulse Programming Course
;
; Project 4: data acquisition details
;
; Explicit parameters
;
; p0 - excitation pulse delay
; d0 - relaxation delay
; d8 - CONFIG_amplifier_enable per convention
; t0 - trim value ACQ_trim[0]
; t1 - trim value ACQ_trim[1] : slice selection
; NI - image object loop
; NA - averaging loop

;-----necessary include files-----
#include<Avance.incl>
#include <DBX.include>

;-----internal parameter definitions and relations-----
define delay deosc = { $DEOSC }
"deosc = abs(deosc) * 1e-6 - aqq"
define delay ramp = { $PREEMP_rise_time }
"ramp = ramp * 1e-6"
define delay dr
"dr = aqq / 2 "
define delay dq
"dq = ramp - d8 - 10u"
define delay df
"df = (p0 + ramp) / 2"

; automatic pulse gating switched off for fast sequences
preset off
; use second synthesizer for LO frequency
fq1b receive

;label--delays--commands---gradients-----remarks-----
start, d0 fq8b:f1 ; select LO freq.

;-----excitation with slice selection -----
dq grad { (0) | (0) | (t1) }
3u fq1:f1 ; set transmit freq.
d8 gatepulse 1 ; pulse gating
(p0:sp0 ph0):f1 ; pulse TPQQ[0]
ramp groff

;-----slice selection refocussing-----
df grad { (0) | (0) | (0) - (t1) }
ramp groff

;-----dephasing gradient-----
dr grad { (0)-(t0) | (0) | (0) }
ramp groff

;-----refocussing gradient-----
ramp grad { (t0) | (0) | (0) } ; Refocussing
    
```

```
-----data acquisition-----
depa    REC_ENABLE                ; open receiver
ADC_INIT_B(ph2, ph1)             ; prepare acquisition
aqq ADC_START                    ; start data sampling
deosc   groff                    ; minimum eosd delay
10u    ADC_END                   ; end data sampling
10u    ipp0 ipp1 ipp2

-----data loops-----
lo to start times NI             ; NI "image" objects
lo to start times NA             ; averaging loop
      10u    rpp0 rpp1 rpp2      ; reset phases

SETUP_GOTO(start)

exit

-----phase definitions-----
ph0 = 0 2 1 3                   ; excitation phase
ph1 = 0                          ; reference phase
ph2 = 0 2 1 3                   ; receiver phase
```

6 Gradient Cycling

6.1 Objectives


By developing this project you will learn

- to vary gradients during runtime
- to use parameters to control gradient variation
- to switch off gradients during setup mode

6.2 Project


In this project you will learn how to vary the strength of a gradient linearly within a pulse sequence. This is the basic for phase encoding in many imaging sequences. You will learn to implement such a gradient “ramp” within a pulse program and understand default links to acquisition parameters.

Under certain conditions it is desirable to acquire a gradient echo not at the middle but at a different position. One way to vary the echo position is to vary the strength of the dephasing gradient. That is what we are going to do in this experiment.



We start by creating a new dataset and copying the pulse program from project 3 to project6 .ppg. Then the trim value of the dephasing gradient will be modified by adding `r2d(t0)` to the gradient value. In the pulse sequence, the length of the gradient ramp must be declared by a statement `lgrad r2d`. The ramp function must be incremented by an `igrad r2d` statement for each loop cycle. To have full access to default parameters for the `r2d` ramp, we have to increase the dimension of the experiment `ACQ_dim = 2`. E.g. `ACQ_size[1] = 4`. Then `lgrad r2d = ACQ_size[1]` and the loop over NI should be

replaced by a loop over `ACQ_size [1]`. The discrete points used for the ramp function can be influenced by the `ACQ_phase_encoding_mode [1]` parameter. The easiest way to observe this is, to select a `User_Defined_Encoding` and select discrete values in the range -1 to +1 in the parameter `ACQ_spatial_phase_1` (the size of the array is controlled by the parameter `ACQ_spatial_size_1`).



6.3 Extensions

Use the `<2d>` modifier in the ramp definition to restrict the gradient function to a single value in setup mode.

6.4 Hints

- To observe several scans simultaneously in the display, select `ACQ_phase_factor = ACQ_size [1]`
- You can start an adapted parameter editor


```
pvcmd -a pvScan pvEditPars PULPROG ACQ_dim ACQ_size[1]
ACQ_phase_factor D[0] P[0] TPQQ[0].name TOQQ[0].power
ACQ_n_trim ACQ_trim[0] ACQ_phase_encoding_mode[1]
ACQ_spatial_size_1 ACQ_spatial_phase_1
```

6.5 Solution

```

;
; PPC 2006 Pulse Programming Course
;
; Project 4: data acquisition details
;
; Explicit parameters
;
; p0 - excitation pulse delay
; d0 - relaxation delay
; d8 - CONFIG_amplifier_enable per convention
; t0 - trim value ACQ_trim[0]
; t1 - trim value ACQ_trim[1] : slice selection
; ACQ_size[1] - image object loop
; NA - averaging loop

;-----necessary include files-----
#include<Avance.incl>
#include <DBX.include>

;-----internal parameter definitions and relations-----
define delay deosc = { $DEOSC }
"deosc = abs(deosc) * 1e-6 - aqq"
define delay ramp = { $PREEMP_rise_time }
"ramp = ramp * 1e-6"
define delay dr
"dr = aqq / 2 "
define delay dq
"dq = ramp - d8 - 10u"
define delay df
"df = (p0 + ramp) / 2"

; definition of gradient ramp
lgrad r2d<2d> = ACQ_size[1]

; automatic pulse gating switched off for fast sequences
preset off
; use second synthesizer for LO frequency
fq1b receive

;label--delays--commands----gradients-----remarks-----
start, d0 fq8b:f1 ; select LO freq.

;-----excitation with slice selection -----
dq grad { (0) | (0) | (t1) }
3u fq1:f1 ; set transmit freq.
d8 gatepulse 1 ; pulse gating
(p0:sp0 ph0):f1 ; pulse TPQQ[0]
ramp groff
;-----slice selection refocusing-----
df grad { (0) | (0) | (0) - (t1) }
ramp groff

;-----dephasing gradient-----
dr grad { (0)-(t0)-r2d(t0) | (0) | (0) }
ramp groff
    
```

```
-----refocusing gradient-----
ramp grad { (t0) | (0) | (0) }      ; Refocusing
-----data acquisition-----
depa    REC_ENABLE                  ; open receiver
ADC_INIT_B(ph2, ph1)                ; prepare acquisition
aqq ADC_START                        ; start data sampling
deosc   groff                       ; minimum eosc delay
10u     ADC_END                     ; end data sampling
10u     ipp0 ipp1 ipp2 igrad r2d     ; increment ramp

-----data loops-----
lo to start times ACQ_size[1]       ; encoding steps
lo to start times NA                 ; averaging loop
10u     rpp0 rpp1 rpp2              ; reset phases

SETUP_GOTO(start)

exit

-----phase definitions-----
ph0 = 0 2 1 3                       ; excitation phase
ph1 = 0                               ; reference phase
ph2 = 0 2 1 3                       ; receiver phase
```



Applied Pulse Programming


ParaVision Programming course
April 03-07, 2006

J. VOIRON

INTRODUCTION

- Interpretation of the „gre“ pulse program (PPG)
- Project 1 : multi-frames & 3D
- Project 2.1: slice refocussing
- Project 2.2 : slice thickness



- ## PROJECT1: „gre“ ppg interpretation
- 
- initialization of 2d phase & slice loops
 - start,
 - set the emitter/receiver frequency offset


 - stabilization time
 - slice selective excitation
 - slice rephase
 - phase encoding
 - read dephase
 - frequency encoding
 - read spoiler
 - phase rewinder
 - relaxation period

 - increment of slice gradient
 - lo to start times NSLICES
 - reset of slice gradient

 - lo to start times NA

 - increment of 2d phase gradient
 - lo to start times ACQ.size[1]
 - reset of 2d phase gradient

 - lo to start times NR

 - phase cycling list
- ← Applied at the same time
(timing optimization)
- ← Applied at the same time
(timing optimization)
- 

PROJECT1: „gre“ ppg interpretation



```

; stabilization time
10u fq8b:f1 (1)
d1 fq1:f1 (2)

(1): only activated with the
ADC_init_B macro

(2): immediatly switched on
    
```

set Rx offset frequency.
set TX

```

start,
stabilization time
slice selective excitation
slice rephase
phase encoding
read dephase
frequency encoding
read spoiler
phase rewinder
relaxation period
lo to start times    NSLICES
                    NA
                    ACQ_size[1]
                    NR
    
```

necessary because the last rep. phase ended with a spoiler.



PROJECT1: „gre“ ppg interpretation



```

; slice selection
d1 fq1:f1 grad{(0)|(0)|(t0)}
d2 gatepulse 1
(p0:sp0 ph1):f1
    
```

```

start,
stabilization time
slice selective excitation
slice rephase
phase encoding
read dephase
frequency encoding
read spoiler
phase rewinder
relaxation period
lo to start times    NSLICES
                    NA
                    ACQ_size[1]
                    NR
    
```



PROJECT1: „gre“ ppg interpretation



```
; phase encoding  
  
d3 grad{(t2) | r2d(t3) | (t1)}  
d5 groff
```

```
start,  
stabilization time  
slice selective excitation  
slice rephase  
phase encoding } simultaneous  
read dephase  
frequency encoding  
read spoiler  
phase rewinder  
relaxation period  
  
lo to start times    NSLICES  
                    NA  
                    ACQ_size[1]  
                    NR
```



PROJECT1: „gre“ ppg interpretation



```
; frequency encoding  
  
denab REC_ENABLE grad{(t5) | (0) | (0)}  
    ADC_INIT_B(ph1,ph0)  
aqq  ADC_START
```

Rx phase
ref. phase.

```
start,  
stabilization time  
slice selective excitation  
slice rephase  
phase encoding  
read dephase  
frequency encoding  
read spoiler  
phase rewinder  
relaxation period  
  
lo to start times    NSLICES  
                    NA  
                    ACQ_size[1]  
                    NR
```



PROJECT1: „gre“ ppg interpretation



```
; read spoiler & phase encoding
```

```
d3 grad{ (t8) | r2d(t6) | (0) }
```

```
d6 grad{ (t8) | (0) | (0) }
```

```
d7 groff
```

```
start,
```

```
stabilization time
```

```
slice selective excitation
```

```
slice rephase
```

```
phase encoding
```

```
read dephase
```

```
frequency encoding
```

```
read spoiler
```

```
phase rewinder } simultaneous
```

```
relaxation period
```

```
lo to start times NSLICE
```

```
NA
```

```
ACQ_size[1]
```

```
NR
```



PROJECT1: gre ppg interpretation



```
; relaxation delay
```

```
d0 ADC_END
```

```
start,
```

```
stabilization time
```

```
slice selective excitation
```

```
slice rephase
```

```
phase encoding
```

```
read dephase
```

```
frequency encoding
```

```
read spoiler
```

```
phase rewinder
```

```
relaxation period
```

```
lo to start times NSLICE
```

```
NA
```

```
ACQ_size[1]
```

```
NR
```



PROJECT1: 3D generalization of „gre“

BRUKER

- Use a new *r3d* ramp function associated to the *ll* loop (to initialize and to add in the loop structure)

- Apply it in the slice selection direction:

```
d3 grad{ (t2) | r2d(t3) | (t1)+ r3d(t4) }
```

- Don't forget to add it for the phase rewinder as well:

```
d3 grad{ (t8) | r2d(t6) | r3d(t7) }
```

BRUKER
BIOSPIN

PROJECT1: 3D generalization of „gre“

BRUKER

- Test procedure:
 - Switch the spatial acquisition parameter to „3D“
 - sinc10H/1ms/30°
 - Matrix = 128x64x64
 - Slice thickness = 20mm
 - TR / TE = 30/3.5ms

BRUKER
BIOSPIN

PROJECT1: multi-frames gre




- Implement a new movie loop (I10) (*)
- Initialize L[10] and NI accordingly (*)
- Modify the object ordering parameter (ACQ_obj_order in ACQ_INFO)

(*) refer to the section 2.2 of your handout...

- Test procedure:
 - Acquire 3 slices at 3 very different positions and 4 frames per slice.
 - Check your image sorting...



PROJECT 2.1: slice refocussing




PROJECT 2.1: slice refocussing



- Follow the realization hints in section 3.4.1 of your handout (page16)...



PROJECT 2.2: slice thickness



PROJECT 2.2: slice thickness

BRUKER

- Implement the modifications according to the sequence diagram of the sequence (also, look at the [section 3.4.2](#) of your handout...)
- Change the *ACQ_grad_matrix* [i][j][k]:

$$\begin{pmatrix} 1 & 0 & 0 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{pmatrix}$$

Image



$$\begin{pmatrix} 0 & 0 & 1 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{pmatrix}$$

Image of Slice

BRUKER
BIOSPIN

PROJECT 2.2: slice thickness

BRUKER

- Export the scan into TOPSPIN
- Extract the middle FID (i.e.#64) from the 2D dataset: type „*rser*“
- Fourier transform: type „*FT*“
- rephase the spectrum (0 and 1 phase correction order) until you get the right slice profile

- Derive the slice thickness using:

$$\Delta Z = (FOV * \Delta f) / SW_h$$

BRUKER
BIOSPIN



Applied Pulse Programming

ParaVision Programming course
April 03-07, 2006

J. VOIRON

INTRODUCTION

- Interpretation of the *gre*^(*) pulse program (*gre.ppg*)
- Project 1 : multi-frames & 3D
- Project 2.1: slice refocussing
- Project 2.2 : slice thickness

(*) simplified version of a gradient echo imaging (FLASH)

PROJECT1: „gre“ ppg interpretation

Stabilization time
Slice selection
Slice rephase

Read dephase
Frequency encoding
Read spoiler

Relaxation period

PROJECT1: „gre“ ppg interpretation

- Definition of delays
- Initialization of 2d phase & slice loops
- start,
- set the emitter/receiver frequency offset
- stabilization time
- slice selective excitation
- slice rephase
- phase encoding
- read dephase
- frequency encoding
- read spoiler
- phase rewinder
- relaxation period

Applied at the same time (timing optimization)

Applied at the same time (timing optimization)

- increment of slice gradient
- lo to start times NSLICES
- reset of slice gradient
- lo to start times NA
- increment of 2d phase gradient
- lo to start times ACQ_size[1]
- reset of 2d phase gradient
- lo to start times NR
- phase cycling list

PROJECT1: „gre“ ppg interpretation

```

10u fq8b:f1 (1)
d1 fq1:f1 (2) grad( (0) | (0) | (t0) )
  
```

(1) : only activated with the ADG_initA macro

(2) : immediately switched on

start,

stabilization time
slice selective excitation

lo to start times

NSLICES
NA
ACQ_size[1]
NR

PROJECT1: „gre“ ppg interpretation

```

d2 gatepulse 1
(p0:sp0 ph1):f1
  
```

start,

stabilization time
slice selective excitation
slice rephase

lo to start times

NSLICES
NA
ACQ_size[1]
NR

PROJECT1: „gre“ ppg interpretation

```

d3 grad( (t2) |x2d(t3) | (t1) )
d5 groff
  
```

```

start,
  stabilization time
  slice selective excitation
  slice rephase
  phase encoding
  read dephase
  frequency encoding
  
```

to start times NSLICES
NA
ACQ_size[!]
NR

PROJECT1: „gre“ ppg interpretation

```

denab REC_ENABLE grad((t5) | (0) | (0) )
  ADC_INIT_B (ph1,ph0)
aqq  ADC_START
  
```

```

start,
  stabilization time
  slice selective excitation
  slice rephase
  phase encoding
  read dephase
  frequency encoding
  read spoiler
  
```

to start times NSLICES
NA
ACQ_size[!]
NR

PROJECT1: „gre“ ppg interpretation

```

d3 grad( (t8) |x2d(t6) | (0) )
d6 grad( (t8) | (0) | (0) )
d7 groff
  
```

```

start,
  stabilization time
  slice selective excitation
  slice rephase
  phase encoding
  read dephase
  frequency encoding
  read spoiler
  phase rewinder
  relaxation period
  
```

to start times NSLICE
NA
ACQ_size[!]
NR

PROJECT1: „gre“ ppg interpretation

```

d0 ADC_END
  
```

```

start,
  stabilization time
  slice selective excitation
  slice rephase
  phase encoding
  read dephase
  frequency encoding
  read spoiler
  phase rewinder
  relaxation period
  
```

to start times NSLICE
NA
ACQ_size[!]
NR

PROJECT1: „gre“ ppg interpretation

```

. #include <avaca.inc>
. #include <SBC.inc>
. project off
.
. define delay denab
.   *denab = d5 - d4 + d6pa.
.
. lgrad r2d=2d = LRF
. rgrad r2d
. lgrad r2d = NSLICES
. Zslice
.
. #include <MEDSPEC.inc>
.
. start,
.   t0u      f230(F)      grad( 0) (0) (0)
.   t1       S21(F)      *paspulse 1
.   d3       (0) (0) (0) | grad((2) (0) (0) (0) |
.   (0) (0) (0) | (0) | grad(
.   denab    REC_ENABLE   ADC_INIT_B(ph1,ph0) grad((5) (0) | (0)
.   aqq      ADC_START
.   d3       (0) (0) (0) | grad((0) (0) (0) |
.   d6       (0) (0) (0) | grad((0) (0) | (0) |
.   d7       ADC_END
.   d8
.
. to start times NSLICES  lslice
. to start times NA      rpp1  zslice
. to start times t0      rpp1  lgrad r2d
. to start times t1
. goto start
. Exit
.
.
. ph0 = 0
. ph1 = 0
  
```

PROJECT IMPLEMENTATION: Recommendations

- All implementation will be based on the **gre** method:
 - to load from the Parameter Editor for MethodClass (from a scan status „New“)
- All proposed pulse programs are stored under:
 - /opt/PPC/PPG/AppliedPPG
- In order to be launched, these PPGs must be copied under:
 - /opt/PV4.0/exp/stan/nmr/lists/pp

PROJECT1: multi-frames „gre“



- In PPG: Implement a new movie loop (l10) (*)
- In ACQP: Initialize L[10] and NI accordingly (*)
- In ACQP: Modify the object ordering parameter (*) (ACQ_obj_order in ACQ_INFO)

Note: the ACQP parameters associated to the frame loop (L[10], NI, ACQ_obj_order) are NOT implemented in the source code of the method
→ to set manually

(*) refer to the [section 2.5.1](#) of your handout...

- Test procedure:
 - Acquire 3 slices at 3 very different positions and 4 frames per slice.
 - Check your image sorting...



Note: also, change the RECO_rotate parameter: (0, 0) → (0.5, 0)



PROJECT1: multi-frames „gre“



```

Initialization of 2d phase & slice loops
- start
- set the emitter/receiver frequency offset
-
- stabilization time
- slice selective excitation
- slice rephase
- phase encoding
- read dephase
- frequency encoding
- read spoiler
- phase rewinder
- relaxation period
-
- Increment of slice gradient
- lo to start times NSLICES
- reset of slice gradient
-
- lo to start times l10
-
- lo to start times NA
-
- Increment of 2d phase gradient
- lo to start times ACQ_size[1]
- reset of 2d phase gradient
-
- lo to start times NR
-
- phase cycling list
    
```



PROJECT1: multi-frames „gre“



Image sorting:

- Image display (output):

```

slice1: fr1(0) fr2(1) fr3(2) fr4(3)
slice2: fr1(4) fr2(5) fr3(6) fr4(7)
slice3: fr1(8) fr2(9) fr3(10) fr4(11)
    
```

- Acquisition order (input):


Sequential:	Interlaced
frame1: sl1(0) sl2(4) sl3(8)	frame1: sl1(0) sl3(8) sl2(4)
frame2: sl1(1) sl2(5) sl3(9)	frame2: sl1(1) sl3(9) sl2(5)
frame3: sl1(2) sl2(6) sl3(10)	frame3: sl1(2) sl3(10) sl2(6)
frame4: sl1(3) sl2(7) sl3(11)	frame4: sl1(3) sl3(11) sl2(7)

→ ACQ_obj_order =

(0) (4) (8) (1) (5) (9) (2) (6) (10) (3) (7) (11) (0) (8) (4) (1) (9) (5) (2) (10) (6) (3) (11) (7)



PROJECT1: 3D generalization of „gre“



PROJECT1: 3D generalization of „gre“



- In PPG: Use a new r3d ramp function associated to the l1 loop (to initialize and to add in the loop structure) (*)
- In PPG: Apply it in the slice selection direction: (*)
`d3 grad((t2) | z2d (t3) | (t1) + z3d (t4))`
- In PPG: Don't forget to add it for the phase rewinder as well: (*)
`d3 grad((t8) | z2d (t6) | z3d (t7))`

Note: the ACQP parameters associated to the third dimension (L[1], ACQ_trim[4] & [7]) are already implemented in the source code of the method.

(*) refer to the [section 2.5.2](#) of your handout...




PROJECT1: 3D generalization of „gre“



Test procedure:

- Switch the spatial acquisition parameter to „3D“ in StandardInplaneGeometry class
- sinc10H/1ms/30°
- Matrix = 128x64x64
- Slice thickness = 20mm
- TR / TE = 30/3.5ms






PROJECT 2.1: slice refocussing

- Follow the realization hints in [section 3.4.1](#) of your handout (based on the sequence diagram)
- Test procedure:
 - Select 1 centered slice
 - ```

P[1]=200us
TPOQ[1].name = „bp.exc“
 .power = PVM_RefAtt - 20 log(1000us/200us) - 6[dB]
 .offset = 0Hz

```
  - GSP mode: observe the echo position when changing the slices in Z (if working in axial orientation)

BRUKER  
SIOSPIN



### PROJECT 2.2: slice thickness

- Implement the modifications according to the sequence diagram (also, look at the [section 3.4.2](#) of your handout...)
- Change the ACQ\_grad\_matrix [i][j][k]:
 
$$\begin{pmatrix} 1 & 0 & 0 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{pmatrix} \rightarrow \begin{pmatrix} 0 & 0 & 1 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{pmatrix}$$

Image Image of Slice

BRUKER  
SIOSPIN

### PROJECT 2.2: slice thickness

Signal processing procedure:

- Export the scan into TOPSPIN
- Extract the middle FID (i.e.#64) from the 2D dataset: type „rser“
- Fourier transform: type „FT“
- rephase the spectrum (0 and 1 phase correction order) until you get the right slice profile
- Derive the slice thickness using:
 
$$\Delta Z = (FOV * \Delta f) / SW\_h$$

BRUKER  
SIOSPIN

### THANK YOU FOR YOUR ATTENTION

BRUKER  
SIOSPIN



---

# Applied Pulse Programming

---

*ParaVision Programming Course  
April 03-07, 2006*

Author:

Jérôme VOIRON

Bruker BioSpin MRI GmbH



|            |                                                                     |           |
|------------|---------------------------------------------------------------------|-----------|
| <b>1</b>   | <b>Some “rehashes” about the pipeline acquisition in ParaVision</b> | <b>3</b>  |
| <b>1.1</b> | <b>Role of the pulse program</b>                                    | <b>3</b>  |
| <b>1.2</b> | <b>Loop description</b>                                             | <b>3</b>  |
| 1.2.1      | Parameter describing the dimensionality of one object               | 4         |
| 1.2.2      | Parameter describing the type of phase encoding                     | 4         |
| 1.2.3      | Parameter describing the number of reconstructed objects:           | 4         |
| <b>1.3</b> | <b>Application to a practical example</b>                           | <b>5</b>  |
| 1.3.1      | Sequence elements and loop structure for the gradient echo sequence | 5         |
| 1.3.2      | Ppg implementation strategy                                         | 6         |
| <b>2</b>   | <b>PROJECT 1: Gradient Echo sequence (“gre”)</b>                    | <b>7</b>  |
| <b>2.1</b> | <b>Objectives</b>                                                   | <b>7</b>  |
| <b>2.2</b> | <b>Project description</b>                                          | <b>7</b>  |
| <b>2.3</b> | <b>“gre” method features</b>                                        | <b>7</b>  |
| <b>2.4</b> | <b>Interpretation of the “gre” pulse program</b>                    | <b>8</b>  |
| <b>2.5</b> | <b>Realization hints</b>                                            | <b>10</b> |
| <b>2.6</b> | <b>Proposed Solution</b>                                            | <b>11</b> |
| <b>3</b>   | <b>PROJECT 2: Slice Selection Test Procedures</b>                   | <b>13</b> |
| <b>3.1</b> | <b>Objectives</b>                                                   | <b>13</b> |
| <b>3.2</b> | <b>Projects description</b>                                         | <b>13</b> |
| 3.2.1      | Project 2.1: Slice refocusing                                       | 13        |
| 3.2.2      | Project 2.2: Slice thickness                                        | 13        |
| <b>3.3</b> | <b>Method features</b>                                              | <b>14</b> |
| 3.3.1      | Project 2.1                                                         | 14        |
| 3.3.2      | Project 2.2                                                         | 15        |
| <b>3.4</b> | <b>Realization hints</b>                                            | <b>15</b> |
| 3.4.1      | Project 2.1                                                         | 15        |
| 3.4.2      | Project 2.2                                                         | 17        |
| <b>3.5</b> | <b>Proposed solutions</b>                                           | <b>18</b> |
| 3.5.1      | Project 2.1 : ppcTest1.ppg                                          | 18        |
| 3.5.2      | Project 2.2 : ppcTest2.ppg                                          | 19        |



# 1 Some “rehashes” about the pipeline acquisition in ParaVision

## 1.1 Role of the pulse program

In ParaVision, a pulse program is just one part of a complete measurement method. The role of the pulse program is to produce experimental data points which can be analyzed and displayed afterwards by the process of reconstruction and image display. Within ParaVision, these basic parts of collecting data and the reconstruction process are linked together in the acquisition & reconstruction pipeline. The pipeline expects a fixed, but parameterized flow of the data. As a consequence of this pipeline architecture, the loop structure for the pulse program is fixed: see below.

## 1.2 Loop description

```

start,
 ...
 ...
 ADC_START (of size ACQ_scan[0])
 ...
lo to start times
 ACQ_ns_list[i]
 ACQ_phase_factor
 NI
 NA
 ACQ_size [1]
 ...
 ACQ_size [n]
 NAE
NR

```

Distinguish between three types of “loop parameters”:

1.2.1 Parameter describing the dimensionality of one object

|                |                                                                                                                          |
|----------------|--------------------------------------------------------------------------------------------------------------------------|
| ACQ_size[0, 8] | Dimension of the experiment.<br><br>The standard reconstruction will apply a Fourier transformation in these dimensions. |
|----------------|--------------------------------------------------------------------------------------------------------------------------|

1.2.2 Parameter describing the type of phase encoding

|                  |                                                                 |
|------------------|-----------------------------------------------------------------|
| ACQ_phase factor | Segmentation type.<br><br>i.e. RARE factor, EPI segmentation... |
|------------------|-----------------------------------------------------------------|

1.2.3 Parameter describing the number of reconstructed objects:

|                    |                                              |
|--------------------|----------------------------------------------|
| NI =               | Number of images. Typically used for :       |
| • NSLICES          | <----> • total number of slices (NSLICES)    |
| • ACQ_ns_list_size | <----> • echo images (T2 measurements)       |
| • L[i]             | <----> • T1 weighted images (T1 measurement) |
| • L[i]             | <----> • diffusion encoding                  |
| • L[i]             | <----> • velocity encoding                   |
| • ...              | • ...                                        |

Note: "L" is a user-defined loop counter

|     |                                                                         |
|-----|-------------------------------------------------------------------------|
| NA  | Number of averages used before the phase-encoding increments.           |
| NAE | Number of averaged experiment used after the phase-encoding increments. |
| NR  | Number of repeated experiments.<br><br>i.e. time course investigations  |

|                          |                                                                                    |
|--------------------------|------------------------------------------------------------------------------------|
| <b>Particular case :</b> | <b>Multi-echo experiment (MSME)</b>                                                |
| ACQ_ns_list[0, NI-1]     | Direct data averaging of size ACQ_size[0].<br>Application: echo collection in MSME |
| NI                       | Number of echo images x number of slices<br>(ACQ_ns_list_size) (NSLICES)           |
| NA, NAE, NR              | Kept unchanged: can be still used.                                                 |

*At first glance, this pipeline structure seems to be somewhat complicated. However, we don't need to write down explicitly all of these loops in the pulse program as long as they are not required. Loops which are not required can be skipped by setting the corresponding base-level parameter to one.*

### 1.3 Application to a practical example

A 2D multi-slices gradient echo sequence will be considered in the following of this chapter.

#### 1.3.1 Sequence elements and loop structure for the gradient echo sequence

The main elements of such a gradient echo sequence are the following:

- relaxation delay within one slice
- delays between consecutive slices
- slice-selective excitation
- phase-encoding
- frequency-encoding

According to the previous loop descriptions, the required loop structure for the 2D gradient echo sequence is schematized below. Its main components are written in bold fonts and the linkage components are in italic.

Please, notice that the following loops were skipped in the following way:

- ACQ\_ns\_list[i] = 1
- ACQ\_phase\_factor = 1
- NAE = 1
- ACQ\_size[i > 1] = 0

```

Initialization of 2d phase loop
start, relaxation delay within one slice
 set the emitter/receiver frequency offset
 slice selective excitation
 phase encoding
 frequency encoding
increment of slice gradient
lo to start times NSLICES (*)
reset of slice gradient
lo to start times NA
increment of 2d phase gradient
lo to start times ACQ_size[1]
reset of 2d phase gradient
lo to start times NR

[phase cycling list]

```

(\*) implicitly related to NI (see § 1.2.3)

### 1.3.2 Ppg implementation strategy

From the previous example, one can derive the following strategy of ppg implementation:

***Combine together the described in-loop elements as well as the delays and gradient switching points.***

***Insert the loop counters and the slice frequency offsets as well as the phase gradient increment statements.***

## 2 PROJECT 1: Gradient Echo sequence (“gre”)

### 2.1 Objectives

By developing this project you will learn how to:

- interpret a complex PPG: 2D “gre” sequence
- modify it in order to obtain a PPG capable of:
  - o multi frames acquisition (typical use: cardio movies)
  - o 3D acquisition: 3D “gre” sequence


### 2.2 Project description

The first part of this project consists of analyzing a realistic 2D imaging sequence: *gre*. This is a simplified version of gradient echo imaging (FLASH) which will be used as a basis of all programming projects developed during this course.

In the next step you will modify the *gre* sequence to make it capable of acquiring movie frames. You should add an additional loop and set the corresponding ACQP parameters.

Finally, you will extend the *gre* sequence to a 3D acquisition.

### 2.3 “gre” method features



*Note: the slice selection spoiler switching points have been combined in a way allowing optimizing the sequence timing.*

- 2D gradient-echo sequence
- Multi-slices
- Multi-packages
- Multi-frames (movie mode)
- Multi-averaging
- Multi-repetitions
- No particular phase-cycling

*The acquisition loop structure (from inner to outer loops):*

- 1- slices
- 2- movie frames
- 3- accumulation
- 4- phase-encoding
- 5- repetitions

## 2.4 Interpretation of the “gre” pulse program (also look at the sequence diagram)

```
#include<Avance.incl>
#include <DBX.include>
preset off
```

```
=====
; definition of delays
;=====
```

```
define delay denab
"denab = d4 - de + depa"
```

```
=====
; declaration of 2d and 3d loop
;=====
```

```
lgrad r2d<2d> = L[0]
zgrad r2d
```

```
lgrad slice = NSLICES
```

```

zslice
#include <MEDSPEC.include>

;=====
;Delay/Pulse spec control gradients
;=====
;-----start of the main loop-----

start,10u fq8b:f1

;-----slice selection-----

 d1 fq1:f1 grad{(0) | (0) | (t0)}
 d2 gatepulse 1
 (p0:sp0 ph1):f1

;-----slice rephase, read dephase, phase encoding-----

 d3 grad{(t2) | r2d(t3) | (t1) }
 d5 groff

;-----frequency encoding-----

 denab REC_ENABLE grad{(t5) | (0) | (0)}
 ADC_INIT_B(ph1, ph0)
 aqg ADC_START

;-----read spoiler, phase encoding-----

 d3 grad{(t8) | r2d(t6) | (0)}
 d6 grad{(t8) | (0) | (0) }
 d7 groff
 d0 ADC_END

;-----slice loop-----

 1u islice
lo to start times NSLICES

;-----averaging loop-----

 1u ipp1 zslice
lo to start times NA

;-----2D loop-----

 1u rpp1 igrad r2d ;2d loop
lo to start times l0

;-----repetition loop-----

 1u
lo to start times NR
goto start
exit
;=====
;phase lists

```

```
ph0 = 0
ph1 = 0
```

---

## 2.5 Realization hints

1. **Implementation of the multi-frames loop.** First, copy gre.ppg in MovieGre.ppg without forgetting to update the PULPROG variable accordingly. ECG triggering considerations impose that all slices are acquired for each movie frame. Thus, you need:

- to implement (considering the sections 2.2 and 2.3 of this course) a new movie loop (l10) placed at adequate position in the ppg loop structure,
- to initialize L[10] and NI accordingly (for NI, refer to its definition at section 2.2 of this course),
- to modify the object ordering parameter (ACQ\_obj\_order in ACQ\_INFO) in order to get the right object sorting, namely all frames displayed for each slice in the expected order. Of course, the object ordering parameter directly depends on the object ordering mode: per default, this mode is set to Interlaced.

In order to test this new feature, acquire 3 slices at 3 very different positions (one crossing the air bubble at the top, one at the centre and the last one at the bottom of the phantom) and set the number of frames to 4. Keep the object ordering mode to Interlaced. Check if your number of acquired frames fits with the expected one, and if they are displayed in the right order (4 frames acquired for each slice in the right order).

2. **Implementation of 3D acquisition.** For that purpose, copy first gre.ppg in 3DGre.ppg without forgetting to update the PULPROG variable accordingly.

You need to implement a new r3d function and to initialize it with the L[1] parameter (associated to the l1 loop parameter). This second phase-encoding gradient (applied in the slice selection direction) will be applied during the slice rephasing lobe of gradient. Use the table below to help you changing the 2D pulse program in a 3D pulse program.

In order to test this new feature, don't forget to switch the spatial acquisition dimension parameter (StandardInplaneGeometry class) to the 3D value (by default this parameter is set to 2D), and use the following parameters to obtain a good 3D image quality:

- Sinc10H / 1ms/30°
- Matrix size = 128x64x64
- Slice thickness = 20mm
- TR = 30ms / TE = 3.5ms

| PPG's COMMENTED PART                                       | DELAY                                      | TRIMS  | COMMENTS                                                                                                    |
|------------------------------------------------------------|--------------------------------------------|--------|-------------------------------------------------------------------------------------------------------------|
| <i>slice rephase,<br/>read dephase,<br/>phase encoding</i> | d3 (slice dephasing / rephasing rise-time) | t1, t4 | During the slice rephasing time, use <i>r3d</i> to set the phase-encoding gradient in the second direction. |



|                       |                        |    |  |
|-----------------------|------------------------|----|--|
| <i>Read spoiler,</i>  | d3 (+rise-time "up")   | t7 |  |
| <i>Phase encoding</i> | d6 (+rise-time "down") |    |  |

## 2.6 Proposed Solution

```

#include<Avance.incl>
#include <DBX.include>
preset off

;=====
; definition of delays
;=====

define delay denab
"denab = d4 - de + depa"

;=====
; declaration of 2d and 3d loop
;=====

lgrad r2d<2d> = L[0]
zgrad r2d
lgrad r3d<3d> = L[1]
zgrad r3d

lgrad slice = NSLICES
zslice
#include <MEDSPEC.include>

;=====
;Delay/Pulse spec control gradients
;=====

;-----start of the main loop-----
start,10u fq8b:f1

;-----slice selection-----

 d1 fq1:f1 grad{(0) | (0) | (t0)}
 d2 gatepulse 1
 (p0:sp0 ph1):f1

;-----slice rephase, read dephase, phase encoding --
 d3 grad{(t2) | r2d(t3) | (t1)+r3d(t4)}
 d5 groff

;-----frequency encoding-----

 denab REC_ENABLE grad{(t5) | (0) | (0)}

```

```

 ADC_INIT_B(ph1, ph0)
 aqg ADC_START

;-----read spoiler, phase encoding-----

 d3 grad{(t8)|r2d(t6)|r3d(t7)}
 d6 grad{(t8)| (0) | (0) }
 d7 groff
 d0 ADC_END

;-----slice loop-----

 lu islice
lo to start times NSLICES

;-----movie loop-----

 lu zslice
lo to start times 110

;-----averaging loop-----

 lu ipp1
lo to start times NA

;-----2D loop-----

 lu rpp1 igrad r2d ;2d loop
lo to start times 10

;-----3D loop-----

 lu igrad r3d ;3d loop
lo to start times 11

;-----repetition loop-----

lo to start times NR
goto start
exit
;=====
;phase lists

ph0 = 0
ph1 = 0
;=====

```

*Note: the lines written in bold fonts correspond to "gre" ppg's modifications according to the Project1 specification.*

## 3 PROJECT 2: Slice Selection Test Procedures

Before running your pulse program in a routine mode, it can be necessary to test for instance the slice selection you have implemented. Two properties of the slice selection pulses can be easily verified simply by modifying your 2D Gradient Echo pulse program. We will thus check:

- the slice rephasing property
- the slice thickness property

### 3.1 Objectives

By developing this project you will learn how to:

- introduce a new  $\pi$  non-selective pulse in your “gre” pulse program file
- modify some ACQP parameters influencing your pulse program and in particular:
  - ACQ\_trim parameters
  - ACQ\_grad\_matrix parameters

### 3.2 Projects description

#### 3.2.1 Project 2.1: Slice refocusing

Based on the former 2D gradient echo (“gre”) pulse program, you will implement a ( $\pi/2$  selective) – ( $\pi$  non-selective) sequence. The non-selective  $\pi$  pulse must be a block pulse (bp) placed at TE/2. Because we are only interesting in checking the reliability of the slice refocusing pulse, the readout and phase-encoding gradients will be ignored.

From a theoretical point of view, if the phase dispersion produced by the slice selection gradient has been properly rephased, then the echo position must not be dependent from any local magnetic field (B0). Practically and assuming the slice selection to be along the z-axis, it means that changing the values of the z-shim must not change the echo position (only the echo shape is expected to change). We will check this property.

#### 3.2.2 Project 2.2: Slice thickness

By slightly modifying the timing of the 2D gradient echo (“gre”) pulse program and setting the readout gradient in the same direction than the slice gradient (modify the ACQ\_grad\_matrix array parameter), you will be able to generate an excitation profile.

The slice thickness will be checked by exporting the acquired 2D image, extracting a single row of the 2D dataset and processing the corresponding signal into TopSpin. The slice thickness will be derived from the slice profile according to the following formula:


$$\Delta Z = (\text{FOV} * \Delta f) / \text{SW}_h \quad (1)$$

If the slice thickness differs from the selected value, check either your method files or the bandwidth of your excitation pulse.

*Note: after implementing these changes in the pulse program, you can no longer expect having right calculated values for the parameters that are derived from the “gre” method relations (TE, TR...).*

### 3.3 Method features


#### 3.3.1 Project 2.1



Special features: modified 2D gradient echo sequence

- ( $\pi/2$  selective)  $-(\pi$  non-selective) sequence
- no gradients applied in the readout and phase-encoding directions

### 3.3.2 Project 2.2



Special features: modified 2D gradient echo sequence

- slice selection and readout gradients in the same direction
- phase-encoding gradient

## 3.4 Realization hints

### 3.4.1 Project 2.1

1. Copy gre.ppg in ppcTest1.ppg and follow the next implementation steps to modify it. Don't forget to update the PULPROG variable accordingly.
2. Remove all the gradients in the readout and phase-encoding directions by setting the corresponding trim values to 0.
3. Define a new delay del in ppcTest1.ppg and insert it after the variable delay d5 :

- `define delay del "del = p0/2+d3+d5+d2-denab-aqq/2"`

*Note: make sure that  $d5$  ( $\approx TE$ ) is big enough to allow a reliable  $d\epsilon 1$  calculation.*

4. Insert a block pulse “bp” just before the delay  $d\epsilon 1$  (use the pulse duration  $p1$ , the shape pulse  $sp1$ , and the phase  $ph2$ ). Define the  $ph2$  frequency list at the end of your ppg :
  - $ph2 = 1$

*Note: the new  $\pi$  pulse parameter value is independent on the “gre” method. Thus you will need to initialize its corresponding ACQP parameters. The Pulse Program Tool can be used to display your modified ppg only if the ACQP parameters have been initialized.*

5. It is also necessary to insert a gate pulse command directly between the  $d5$  duration and the block pulse “bp” (use the  $d2$  delay).
6. You are now able to test your new sequence by keeping the predefined shims. Follow the following recommendations :
  - select only one slice and center it
  - $TE = 50\text{ms} / SWH = 50000\text{Hz}$ , Matrix size =  $1024 \times 128$   
( $\Rightarrow d5 = 38\text{ms} / aqq = 20.5\text{ms}$ )
  - In the ACQP class:
    - Set the PULPROG parameter according to your new ppg’s name
    - $ACQ\_O1B\_list = 3000\text{Hz}$  (setting the receiver frequency offset off-resonance will make the echo monitoring easier)
    - $P[1] = 200\mu\text{s}$
    - $TPQQ [1]$ 
      - .name = “bp”
      - .power = to calculate according the formula below (2)
      - .offset = 0Hz

*Note: use the following formula to derive the attenuation value (in dB) for the block-pulse  $180^\circ / 200\mu\text{s}$ :*

$$TPQQ[1].power = PVM\_RefAtt - 20 \log(1000\mu\text{s}/200\mu\text{s}) - 6 \text{ [dB]} \quad (2)$$

*with  $PVM\_RefAtt$  : reference attenuation calculated for a  $90^\circ/1\text{ms}$  block-pulse shape*

7. Run your sequence in GSP, and observe the echo position by changing the shim values in the Z direction (use the shim tool for that purpose). Furthermore, don’t forget to reset your changes before closing the shim tool. If your slice selection gradient is well calculated, no shift of the echo position is expected.
8. If you change the slice rephasing gradient trim value by about 5% (use the  $ACQ\_trim[1]$  parameter), you can now simulate a bad trim calculation. In such a case, expect a shift of the echo position is to expect.

### 3.4.2 Project 2.2

1. Copy gre.ppg in ppcTest2.ppg and follow the next implementation steps to modify it. Don't forget to update the PULPROG ACQP variable accordingly.
2. Move the dephasing lobe of the readout gradients (t2 trim) after the delay d5. You will need to add a new d7 delay for the ramp time down of the dephasing gradient lobe.
3. After updating the PULPROG variable according to ppcTest2.ppg, visualize the pulse program with the Pulse Program Tool. Then, by considering the pulse program diagram related to the project 2.2 change the values of the ACQ\_grad\_matrix parameter accordingly.

Assuming the axial orientation as default spatial orientation, the required changes can be described as followed:

$$\text{ACQ\_grad\_matrix}_{\text{init}} = \begin{bmatrix} 1 & 0 & 0 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{bmatrix}$$

⇓

$$\text{ACQ\_grad\_matrix}_{\text{final}} = \begin{bmatrix} 0 & 0 & 1 \\ 0 & 1 & 0 \\ 0 & 0 & 1 \end{bmatrix}$$

Check your pulse program again by re-opening the Pulse Program Tool. You will need to compile your ppg in X/Y/Z in order to really see what happened along the 3 main axis.

4. You can now measure the slice thickness. Use the following parameters:
  - sinc10H / 2ms
  - TR = 1000ms
  - Position your slice outside the inhomogeneous areas

Export the scan into TopSpin and determine the half height of the phase profile after Fourier transformation. For that purpose, type into TopSpin the following commands:

- rser for extracting a fid number from the 2D dataset : select the middle fid # 64,
- FT for applying the Fourier transformation,
- use the phasing icon (TopSpin panel) for phasing the spectrum: adjust the zero and first phase correction order by pressing the 0 and 1 buttons (panel at the top of the current spectrum window) to improve the phasing.

By applying the formula (1), derive the value of the slice thickness and compare it to the expected value.

### 3.5 Proposed solutions

#### 3.5.1 Project 2.1 : ppcTest1.ppg

```

;=====
; definition of delays
;=====

define delay denab
"denab = d4 - de + depa"
define delay del
"del = p0/2 +d3 +d5 +d2 - denab -aqq/2"

;=====
; declaration of 2d and 3d loop
;=====

#include <MEDSPEC.include>

;=====
; D/P spec control gradients
;=====
start,10u fq8b:f1
 d1 fq1:f1 grad{(0) | (0) | (t0)}
 d2 gatepulse 1
 (p0:sp0 ph1):f1
 d3 grad{(0) | (0) | (t1)}
 d5 groff
 d2 gatepulse 1
 (p1:sp1 ph2):f1
 del
 denab REC_ENABLE
 ADC_INIT_B(ph1, ph0)
 aqq ADC_START
 d3
 d6
 d7
 d0 ADC_END
goto start
exit
;=====
;phase lists

ph0 = 0
ph1 = 0
ph2 = 1
;=====

```



## 3.5.2 Project 2.2 : ppcTest2.ppg

```

#include<Avance.incl>
#include <DBX.include>
preset off

;=====
; definition of delays
;=====

define delay denab
"denab = d4 - de + depa"

;=====
; declaration of 2d and 3d loop
;=====

lgrad r2d<2d> = L[0]
zgrad r2d

lgrad slice = NSLICES
zslice


#include <MEDSPEC.include>

;=====
; D/P spec control gradients
;=====
start, 10u fq8b:f1
 d1 fq1:f1 grad{(0) | (0) | (t0)}
 d2 gatepulse 1
 (p0:sp0 ph1):f1
 d3 grad{(0) | r2d(t3) | (t1)}
 d5 groff
 d3 grad{(t2) | (0) | (0)}
 d7 groff
 denab REC_ENABLE grad{(t5) | (0) | (0)}
 ADC_INIT_B(ph1, ph0)
 aqq ADC_START
 d3 grad{(t8) | r2d(t6) | (0) }

 d6 grad{(t8) | (0) | (0) }
 d7 groff
 d0 ADC_END
 lu islice
lo to start times NSLICES
 lu ipp1 zslice
lo to start times NA
 lu rpp1 igrad r2d
lo to start times 10
 lu
lo to start times NR



```

```
goto start
exit
;=====
;phase lists
ph0 = 0
ph1 = 0
;=====
```




ParaVision Methods  
(PVM)

ParaVision Programming Course  
April 3–7, 2006



Development in *ParaVision*


- Base-level programming
  - Pulse program
  - Acquisition parameters
- High-level programming (Methods)
  - Intuitive parameterization
  - Control of consistency
  - Graphic interface (geometry)



## METHODS

BRUKER

- Experiment description



BRUKER  
BIOSPIN

## Advantages of PVM vs. IMND

BRUKER

- One binary file per method
- Separated source code
  - one directory per method
- Easy and safe installation
  - one file to be exported
  - distribution with/without sources
  - methods can be private for one user
- New features
  - Parameter redirections
  - Local parameters
  - Toolbox library
  - Methods have full control of base-level parameters

BRUKER  
BIOSPIN


## Basics of PVM Programming



- Parameter relations
- Method structure
- Tools for user development
  - Toolbox library
  - Modules



## Parameter relations




```
int parameter
{
 display_name "I'm positive";
 relations checkMe;
} Me;
```

```
void checkMe(void)
{
 if(Me < 1)
 Me = 1;
}
```




## Chain of relations



```
void relA(void)
{
 ParxRelsParRelations("parB", Yes);
}
```




## Relation redirection




```
relations X LocalRelX;
```



# Group redirection



# METHOD STRUCTURE



## Method source code:



### Selected files in the directory FLASH/

need not  
modify

|                          |                            |
|--------------------------|----------------------------|
| <b>Makefile</b>          | compilation recipe         |
| <b>FLASH.c</b>           | glue file                  |
| <b>parsDefinitions.h</b> | parameter declarations     |
| <b>parsLayout.h</b>      | definition of MethodClass  |
| <b>callbackDefs.h</b>    | redirections               |
| <b>parsRelations.c</b>   | backbone & local relations |
| <b>initMeth.c</b>        | function initMeth()        |
| <b>FLASH.ppg</b>         | pulse program              |



## Location of methods



### SOURCES: Bruker:

**`$PvHome/prog/parx/src/`**

Local, typically:

**`$HOME/methods/src/`**

### BINARIES: Bruker:

**`$PvHome/prog/parx/methods`**

Local, usable for all:

**`$PvHome/prog/parx/pub`**

Local private:

**`$ParxMethodSearchPath, typically:`**

**`$HOME/methods/bin/`**





## Tools for user development



- Example methods with source code
  - FLASH, RARE, GEFC,
  - SINGLEPULSE, CSI,
  - FISP, EPI, FAIR\_EPI
  - DtiStandard, DtiEpi (licensed)
  - gre - basis for this course
- Toolbox libraries
- Modules
- copyMethod script



## PVM Toolbox



- Documented library
  - .../docu/english/pvman/pvm\_docu/index.html
- Validity guaranteed for future versions (interfaces will never change)
- Organisation
  - ATB\_ acquisition
  - STB\_ sequence
  - PTB\_ paravision
  - MRT\_ MR utilities - eg. calc gradient strength for a particular resolution
  - UT\_ general utilities
  - CFG\_ configuration always better to rely on CFG files to maintain forward compatibility



## Toolbox documentation



```
double MRT_MinSliceThickness(double sliceThickHz,
 double sliceGradRatio,
 double limitingSliceGrad,
 double limitingRephaseGrad,
 double gradCalConst)
```

Calculate the minimum slice thickness  
There are no error conditions for this function.

Toolbox Classification:  
FUNCTION

Toolbox Functions Called: MRT\_MaxSliceGrad, MRT\_SliceThickness.

Parameter Relations Called:  
NONE

Parameters:  
sliceThicknessHz - Thickness of slice (Hz).  
sliceGradRatio - The ratio of the strengths of the slice ...



## PVM Toolbox – examples



```
STB_InitRFPulse(&ExcPulse, "gauss", 1., 90.);

gconst = CFG_GradCalConst();


Acqtime = MRT_AcqTimeEffSWh(effSWh, size[0]);

ATB_SetNI(nslices*nechoes);
```




# PVM programming: first steps

## gre the template method



PARCOMP - DEBUG = )



### Relations and range checking, cont.

```
int parameter { parsDefinitions.h
 relations myParRels;
} myPar;


void myParRels () { parsRelations.c
 myParRange ();
 backbone ();
}

void myParRange () {
 if (!ParxRelsParHasValue ("myPar"))
 myPar = 5; Initial value
 else
 myPar = MAX_OF (myPar, 0); Legal range
}
```

*built in fcn.*

BRUKER BIOSPIN

# spGre: Introducing a spoiler



↑  
new addition



## RF Pulses in PVM



MyInPulse

BRUKER

|                 |                  |                          |
|-----------------|------------------|--------------------------|
| Length          | 1.000 us         | High Level Attributes    |
| Bandwidth       | 1010.0 Hz        |                          |
| FlipAngle       | 100.0 Deg        |                          |
| Attenuation     | 30.0 dB          | Customization Attributes |
| TrimBandwidth   | 100.00 %         |                          |
| TrimAttenuation | 0.00 dB          |                          |
| TrimRephase     | 100.00 %         |                          |
| Classification  | LB INVERSION     | Low Level Attributes     |
| Filename        | *shape.mv*       |                          |
| BandwidthFactor | 1010.000000 Hzms |                          |
| IntegralRatio   | 0.618172         |                          |
| RephaseFactor   | 0.000000 %       |                          |
| MinPulseLength  | 0.102100 us      |                          |
| ShapeType       | COMBINATION      |                          |

**Wave file fields:**

- SHAPE\_EXMODE
- SHAPE\_BWFAC
- SHAPE\_INTEGFAC
- SHAPE\_REPHFAC
- NPOINTS
- SHAPE\_TYPE

*not exactly correct but can experiment to determine what is the correct classification for each pulse*

*looks editable but is not within a structure that can't tell if it is editable immediately*

BRUKER BIOSPIN

## Definition of an RF pulse

BRUKER

```
PVM_RF_PULSE_TYPE parameter
{
 display_name "Excitation Pulse";
 relations ExcPulseRel;
}ExcPulse;
```

- PVM\_RF\_PULSE\_TYPE is a structure representing physical properties of an RF pulse (length, bandwidth, etc.)
- You can define several pulses
- You can define pulse arrays

BRUKER BIOSPIN

## RF-pulse parameter relations



```
void ExcPulseRel(void)
{
 UT_SetRequest("ExcPulse");
 ExcPulseRange();
 backbone();
}
```

← Needed for  
a correct update  
(correct event handling  
on parameter edit)

```
void ExcPulseRange(void)
{
 individual restrictions (legal flip-angle, classification, etc.)
 ...
 STB_CheckRFPulse(&ExcPulse);
}
```

← given as a pointer to the function  
General check

↑ convenient range check - library, few



needs to change its contents

## Update of RF Pulses



To be called in backbone()

```
void updateRFPulses(char *nucleus)
{
 YesNo refAvailable;
 double refAtt=30;
 refAviable = CFG_GetGlobRefAtt(nucleus, &refAtt);
 STB_UpdateRFPulse("ExcPulse", &ExcPulse,
 refAvailable, refAtt);
}
```

← returns "yes" if attenuation has been set yet?

← convenient range checker, contained in library



## RF Pulses: transfer to base level



Called in BaseLevelRelations()

```
void setPpgParameters(void)
{
 ...
 TPQQ[0].power = ExcPulse.Attenuation;
 TPQQ[0].offset = 0.0;
 ParxRelsParRelations("TPQQ", Yes);

 P[0] = ExcPulse.Length * 1000;
 ParxRelsParRelations("P", Yes);
 ...
}
```

Conversion  
ms →  $\mu$ s



## RF Pulse Enums



- Easy access to available pulse shapes
- Automatic update from .../lists/wave

|                         |         |
|-------------------------|---------|
| Excitation Pulse Choice | gauss   |
| Excitation Pulse        | expand  |
| Refocusing Pulse Choice | hermite |
| Refocusing Pulse        | expand  |

Enums are related to pulses





## Using RF-Pulse Enums



definition

```
PV_PULSE_LIST parameter
{
 display_name "Excitation Pulse Choice";
 relations ExcPulseEnumRel;
} ExcPulseEnum;
```

initialisation

```
STB_InitExcPulseEnum("ExcPulseEnum");
```

Exc, Rfc, Inv, depending on classification  
of the related pulse



## Using RF-Pulse Enums, cont.



Connection PULSE → ENUM (after pulse update)


```
STB_UpdateExcPulseEnum(
 "ExcPulseEnum",
 &ExcPulseEnum,
 ExcPulse.FileName,
 ExcPulse.Classification);
```

Connection ENUM → PULSE (in enum's rels)

```
STB_UpdateExcPulseName("ExcPulseEnum",
 &ExcPulseEnum,
 ExcPulse.FileName,
 &ExcPulse.Classification);
```

Exc, Rfc, Inv







## Bandwidth Scaling in Spin Echo



An echo of the 1kHz-90 – 1kHz-180 sequence is not 1kHz broad.



## Bandwidth Scaling in SE, cont.



Scaling factor depends on shape combination:


|                    |       |
|--------------------|-------|
| Gauss – Gauss:     | 75.8% |
| Hermite – Hermite: | 85.8% |
| Sinc3 – Sinc3:     | 85.5% |

Other combinations not measured.


See Advanced Users Manual, chapter 5.




# ItGre: Timing modification



# spinEcho: Introducing an RF pulse





BRUKER

Modular design

- Find elements which often repeat
- Make them building blocks for the user
- Focus:
  - Interface for the user
  - Interface to the rest of the sequence
- Examples:
  - Contrast preparation (fat sup, IR, MT)
  - Diffusion preparation
  - EPI readout

*P/4.0: slice selective inversion*

BRUKER  
BIOSPIN


## Module components



- Piece of pulse program
- Group of parameters
- Library functions:
  - Initialisation
  - Update
  - Calculation of base-level parameters
- Interface to the main code
  - Global parameters (e.g. module duration)
  - Update arguments (e.g. nucleus)



## Usage example: Inversion-Recovery



### PULSE PROGRAM

```
start, d0
#include<inversion.mod>
Excitation
Acquisition
10 to start times Ny
```

### Method code

```
Update fixedTime
Update Inversion group:
STB_UpdateInversion(fixedTime);


Update repetition time:
minTR = PVM_InvModuleTime + ...;
TR = MAX_OF(TR, minTR);

Set base-level parameters:
ATB_SetInversionBaseLevel();
```



don't turn on too many modules at once!

### Saturation Slices Module



**Parameter Editor for Sat. Slices Parameters**

|                                 |           |
|---------------------------------|-----------|
| Number of saturation slices     | 0         |
| PVM_Fov/SatSlice Drive          | 0.0 x1    |
| Saturation slice thickness      | 10.195 mm |
| Min. saturation slice thickness | 2.548 mm  |
| Saturation slice offset         | 0.000 mm  |
| Saturation slice vector         | _expand   |
| Saturation Pulse List           | ref0      |
| Saturation Pulse                | _expand   |
| PVM_Fov/Set/Ref                 | _expand   |
| Gradient spill duration         | 0.360 ms  |
| Sybil gradient strength         | 50.000    |
| Set of Slice Module Duration    | 1.000 ms  |

OK Apply Cancel Reset Cancel


**BIOSPIN**

### greTag: Creating a tagging module

New parameter group:  
Tagging


Contains

- Tag\_Spacing
- Tag\_Gradient
- Tag\_GradientTime
- Tag\_RFPulse
- Tag\_ModuleTime



**BIOSPIN**

## greTag, cont.




## Request Handling



- Problems for group updating:


- Who was changed?
- Which structure field?
- Which array element?



## Request Handling, cont.



- Using UT\_Get/SetRequest()



## Request Handling, cont.



```
void myParRels(void)
{
 UT_SetRequest("myPar");
 myParRange();
 backbone();
}
```

```
void updateMyGroup(void)
{
 const char *const parlist = "myPar,anotherPar";
 const char *const structlist = "amplitude,duration";
 int requestDetected, parNo, structNo;

 requestDetected = UT_GetRequest(
 parlist,&parNo,
 structlist,&structNo);
}
```





## Project: Accelerated Imaging (encGre)

ParaVision Programming Course  
April 3–7, 2006  
Sascha Köhler

## New parameter group

### Encoding:

- Parallel Imaging
- Linear and centric k-space sampling
- Segmentation
- Partial Fourier encoding
- Zero filling




## Parallel Imaging: Advantages




- Reduce scan time, but keep image quality

Time-of-flight angiography of rat brain:  
matrix size = 256 x 256 x 128, resolution: 117 x 148 x 250  $\mu\text{m}$

10 min




5 min, Acceleration = 2



identical image quality



## Encoding group



## Your Task

BRUKER

### Implementation of new parameter group *Encoding* into PVM method *gre*.

Steps:

- 1) Use `copyMethod` to copy *gre* to *encGre*.
- 2) Include *Encoding* in the file `parsLayout.h`.
- 3) In `callbackDefs.h` the group handler must be directed to the method's *backbone* function:

- relations PVM\_EncodingHandler  
backbone;

BRUKER  
BIOSPIN

## Your Task

BRUKER

Steps:

- 5) In the *backbone* function in `parsRelations.c` the *Encoding* group must be updated:

```
STB_UpdateEncoding(PTB_GetSpatDim(), /* dimension */
 PVM_Matrix, /* image size */
 PVM_AntiAlias, /* a-alias */
 &seg_size, /* segment size */
 SEG_SEQUENTIAL, /* segmenting mode */
 Yes, /* PI in 2nd dim
allowed */
 Yes, /* PI ref lines in 2nd dim
allowed */
 No); /* partial ft in 2nd dim
allowed */
```

BRUKER  
BIOSPIN

## Your Task



### Steps:

- 6) Use *PVM\_EncMatrix* for experiment time calculation in **parsRelations.c**.

(PVM\_EncMatrix describes the acquisition matrix and is defined as product of image size and anti-aliasing factor)

- 7) Use *PVM\_EncMatrix* for setting *ACQ\_size* in **BaseLevelRelations.c**.

- 8) Set ACQP parameters for *User Encoding* in **BaseLevelRelations.c**.



## Your Task



### Steps:

- 9) Setting base-level parameters for controlling a multiple receiver acquisition:

- toolbox function *ATB\_SetMultiRec()*

➤ Yes

- special pulse program required (only for older versions of spectrometer hardware)

➤ No (for Avance II)

```
If(Yes == ATB_SetMultiRec())
```

```
ATB_SetPulprog("encGre.4ch");
```




## Your Task



Steps:

10) Multi receiver data require a special reconstruction mode which is defined in the **RecoRelations.c**:

- copy **RecoRelations.c** from the master solution receiving separately



## Additional Exercise



Segmented acquisition

- Use *PVM\_RareFactor* for updating the Encoding group

The purpose of this exercise is to demonstrate the functionality of the new introduced Encoding group:

- only a few steps are necessary to realize segmented acquisition





---

# Method Programming in PVM

## Part 1

---

*ParaVision Programming Course*  
*April 3-7, 2006*

Authors:

Franciszek Hennel

Sascha Köhler

Markus Wick

Bruker BioSpin MRI GmbH







## Table of Contents

|                                                                        |           |
|------------------------------------------------------------------------|-----------|
| <b>Introduction</b>                                                    | <b>5</b>  |
| <b>1 First Steps - A copy of the gre method (project myGre)</b>        | <b>7</b>  |
| 1.1 Objectives                                                         | 7         |
| 1.2 Project description                                                | 7         |
| 1.3 Specific method parameters                                         | 8         |
| 1.4 Code organization in gre                                           | 8         |
| 1.5 Your modifications                                                 | 9         |
| 1.5.1 Preparing the programming environment:                           | 9         |
| 1.5.2 Using copyMethod                                                 | 9         |
| 1.5.3 Activating debug messages                                        | 10        |
| 1.5.4 Typical compilation errors                                       | 10        |
| 1.5.5 A typical runtime error                                          | 11        |
| <b>2 Method Parameters: Adding a spoiler to gre (project spGre)</b>    | <b>12</b> |
| 2.1 Objectives:                                                        | 12        |
| 2.2 Project description                                                | 12        |
| 2.3 Specific method parameters                                         | 13        |
| 2.4 Realisation hints                                                  | 13        |
| 2.4.1 Starting point                                                   | 13        |
| 2.4.2 Steps                                                            | 13        |
| 2.5 Proposed Solution                                                  | 14        |
| <b>3 Timing Control: Long echo-time version of gre (project ltGre)</b> | <b>18</b> |
| 3.1 Objectives                                                         | 18        |
| 3.2 Project description                                                | 18        |
| 3.3 Specific method parameters                                         | 20        |
| 3.4 Special features                                                   | 20        |
| 3.5 Realisation hints                                                  | 20        |
| 3.5.1 Starting point                                                   | 20        |
| 3.5.2 Steps                                                            | 20        |
| 3.6 Possible pitfalls                                                  | 21        |
| 3.7 Possible extensions                                                | 21        |
| 3.8 Proposed Solution                                                  | 21        |
| <b>4 RF Pulses in PVM: A spin-echo method (project spinEcho)</b>       | <b>24</b> |
| 4.1 Objectives                                                         | 24        |
| 4.2 Project description                                                | 24        |
| 4.3 Specific method parameters                                         | 25        |
| 4.3.1 New method parameters                                            | 25        |
| 4.3.2 Important method parameters                                      | 25        |
| 4.3.3 Important baselevel parameters                                   | 26        |
| 4.4 Important toolbox routines                                         | 26        |
| 4.5 Realisation hints                                                  | 26        |
| 4.5.1 Starting point                                                   | 26        |
| 4.5.2 Refocusing pulse in the pulse program                            | 26        |
| 4.5.3 Refocusing pulse in the method                                   | 27        |
| 4.5.4 Pulse Choice Enum (optional)                                     | 27        |
| 4.5.5 Controlling the spin-echo position                               | 28        |
| 4.6 Possible extensions                                                | 28        |

|       |                                                                   |    |
|-------|-------------------------------------------------------------------|----|
| 4.6.1 | Different amplitudes of excitation and refocusing slice gradients | 28 |
| 4.7   | Possible pitfalls                                                 | 29 |
| 4.8   | Proposed Solution                                                 | 30 |

## Introduction

As you know from the previous course, the acquisition and reconstruction are controlled by a pulse program and by numerous “base-level” parameters grouped mainly in the ACQP and RECO classes. In principle, it is possible to set up an experiment by editing the pulse program setting these parameters directly. On the other hand, one usually prefers to specify the measurement in terms of “high level” parameters which may not have direct representation in ACQP or RECO. For example, it is more convenient to specify the echo time than to set delays used in the pulse program. The mechanism linking the high-level requests with the base-level parameters is provided by the **ParaVision Methods (PVM)**. When you type in the scan editor: “Echo time = 50 ms”, it is the currently active method which checks if this value is allowed and calculates the corresponding delays. Similarly, when you move a slider controlling the field of view in the geometry editor, the active method takes care of not letting you go below the minimum and calculates the necessary gradient amplitudes. In general, the **role of a method** is to:

- provide high-level parameters for experiment description,
- give the parameters initial values in a new scan,
- decide which parameters should be visible, editable and saved in protocol files,
- react to loading a protocol,
- react to changes of parameters made by the user in the geometry editor, scan editor, method editor, or by a shell command,
- check if the new parameter values are legal, and if all parameters agree with each other,
- derive base-level parameters from high-level parameters.

During this course you will learn programming PVM methods by examples. We will start by an overview of an existing **template method gre** and by preparing your working environment for PVM programming. The first programming example will be a simple copy of *gre* just to see how methods are compiled and installed.

In the next project (*spGre*) a user-defined spoiling gradient pulse will be introduced to the *gre* sequence. This example will show you how to add **new parameters** to a method, control their legal range, and program dependencies with other parameters (relations).

Another project (*spinEcho*) will you teach you the handling of **radio frequency pulses** in PVM methods. A special structure type which describes RF pulse properties will be introduced together with toolbox functions for controlling such structures. You will understand the mechanism which derives a necessary attenuation for a given pulse flip angle and shape, pulse bandwidth for a given duration and vice versa, etc. As you will see, all this is achieved by calling a proper updating function for a user defined RF pulse parameter. Additionally, you will also learn using pulse enums, auxiliary parameters allowing an easy selection of a valid pulse shape.

In the same project (*spinEcho*) you will also see how to use **array parameters**, change their sizes and set their elements. This will be shown on the basis of a local frequency list for the refocusing RF pulse.

Each **exercise project** is described step by step in its "**Realization hints**" section. These hints are sufficient (we believe) to complete the projects but still remain open to your invention. Some projects contain suggestions for further extensions which you can follow if the time allows. Should you need more detailed guidance, each project description also contains elements of source-code listing which can simply be copied to your files. Finally, if you prefer analyzing a **ready solution** rather than programming yourself, all projects have been pre-programmed and you can find their complete code in the /opt/PPC/PVM directory.

# 1 First Steps - A copy of the gre method (project myGre)

## 1.1 Objectives

This projects explains:


- How to prepare your login settings for PVM programming
- How to use the *copyMethod* script to produce a new method based on an existing example
- How to compile a new method and make it active in ParaVision
- How the source code of the template PVM method *gre* is organized.

## 1.2 Project description

This project is our first step in PVM programming. We will prepare the method programming environment in your home directory and produce a simple copy of the method *gre*. This is a simplified gradient echo method provided in ParaVision 3.0.2 as a programming template. It will be used throughout this course as a starting point for your projects.

We will have a look at the organization of the code of the new method (*myGre*) and discuss the purpose of all files. Debugging messages will be activated and we will be able to follow the function calling chain resulting from a parameter change, from a new scan selection, or loading a protocol. This exercise will help you understand the role of methods in ParaVision.

Finally, we will purposely produce a few typical compilation errors and runtime errors in the source code so that you can recognize and avoid them easily in future.



### 1.3 Specific method parameters

The method contains several parameters defining the sequence timing, such as the duration of phase encoding or readout dephasing. They are defined in the diagram above.

### 1.4 Code organization in gre

The source code is distributed in several files:

- Makefile                      The control of the compilation
- parsDefinitions.h            Definitions of main-level parameters. Other definition files are included in this file by #include
- parsLayout.h                 Definition of the MethodClass. Other definition files are included in this file by #include
- parsRelations.c             Relations of parameters from the main level (e.g. PVM\_EchoTime) and the backbone function. The backbone is called from all parameter relations and is responsible for inter-dependence of parameters. At its end the backbone calls SetBaseLevelParm() and SetRecoParams().
- BaseLevelRelations.c        Contains SetBaseLevelParm() which sets all ACQP parameters based on PVM parameter values established by the backbone.

FreqEncPars.h, FreqEncLayout.h, FreqEncPars.c

Code responsible for the group of parameters describing the frequency encoding. The h-files contain parameter definitions and the definition of the FrequencyEncoding parclass; they are included in the parsDefinitions.h and parsLayout, respectively. The c-file contains

- range checkers of all parclass members (e.g. ReadDepthTimeRange() for ReadDepthTime)
- relations for all parclass members (e.g. ReadDepthTimeRels())
- InitFreqEncoding() – a function giving all group members default values, used in initMeth()
- FreqEncodingLimits() – a function calculating the minimum FOV in the frequency encoding direction. Used in backbone before the update of geometry.
- FreqEncodingGradients() – a function calculating all gradient amplitudes used on the readout channel. Used in the backbone after the update of geometry

SliceSelPars.h, SliceSelLayout.h, SliceSel.c

Files which manage the slice selection channel in a similar way as the ones for the Frequency Encoding.

Phase2DPars.h, Phase2DLayout.h, Phase2D.c,

Phase2DPars.h, Phase2DLayout.h, Phase2D.c

RFPulsePars.h, RFPulseLayout.h, RFPulse.c

Similar for the 2D and 3D phase encoding gradients and for RF pulses.

## 1.5 Your modifications

### 1.5.1 Preparing the programming environment:

You need to create a directory for the method source code sub-directories and another one for the binary method files. The typical locations are:

```
$HOME/methods/src
$HOME/methods/bin
```

Create these directories.

Now you have to tell ParaVision where to search for method binary files by defining the environment variable *ParxMethodSearchPath*. Edit the file *.bashrc* in your home directory and add the following line:

```
export ParxMethodSearchPath=$HOME/methods/bin
```

Log out and in again to make this definition active.

### 1.5.2 Using copyMethod

Use copyMethod to produce a copy of *gre*, called *myGre* in your source directory. Type

/opt/PV4.0/prog/bin/copyMethod

and follow the dialog. The script will automatically start the compilation of the new method. Restart ParaVision and check if the new method is available by selecting it in a new scan. Try to display its pulse program (do not forget to copy it in the pp directory). Acquire some images with *myGre* and get familiar with the meaning of its parameters.

### 1.5.3 Activating debug messages

You will notice that each function in the *myGre* source code includes a pair of **DB\_MSG** macro calls:

```
void backbone(void)
{
 DB_MSG(("----> backbone"))
 ...
 DB_MSG(("<--- backbone"))
}
```

These are the so called debugging messages used to show the function calling chain during the execution of parameter relations. The **DB\_MSG** macro behaves exactly like the `printf` function – it can be given a format string followed by a list of numeric arguments, or just a simple text. The output goes to the startup shell window of ParaVision. The only difference is the necessary double pair of brackets.

To activate the debug macros you need to

- define the macro **DEBUG** to 1 in each c file in which the messages should be active
- define the macros **DB\_MODULE** and **DB\_LINE\_NR** to one if you want the debug messages to be accompanied by file names and line numbers, respectively.

A typical definition is

```
#define DEBUG 1
#define DB_MODULE 1
#define DB_LINE_NR 0
```

*DB-MSG specification same as printf*

*from C DEBUG FLAGS remove DNDEBUG*

**Your task:** Remove the **DNDEBUG** option from the Makefile, activate the debug messages in `initMeth.c`, `parsRelations.c` and `frecEnc.c`, and de-activate them in `BaseLevelRelations.c`.

Recompile the method by typing:

```
make clean cproto depend
make install
```

Restart Paravision and produce the following situations:

- produce a new empty scan and select the *myGre* method in Method editor
- modify the value of *Read Spoiler Duration*
- select a different scan, and go back to the one containing *myGre*

Each time, note what debug messages are printed in the startup window. Reconstitute the function calling chain. Check if it corresponds to the graph of a recommended method structure shown in Advanced Manual, Method Programming.

### 1.5.4 Typical compilation errors

Syntax error

Remove a colon at the end of a code line (anywhere) and recompile by *make install*. Find the error by reading the compiler messages. If you are using the emacs editor, you can browse



through the compilation errors with the F6 key and be directed to the corresponding code positions automatically.

Missing prototype warning

Open `parsRelations.c` and at the end of the file write the following function

```
void dummy(void) {}
```

Call this function anywhere in the backbone. Compile by typing

```
make install
```

A warning message will appear that the function *dummy* was called without a prototype. In C every function must be declared (by its prototype) before it is used. However, in method programming you do not need to care about prototypes. The only thing you need each time you define a function (or change a function's interface) is to recompile by typing

```
make clean cproto depend
```

followed by

```
make install
```

Now your modified `myGre` will compile OK. Next time, when your code changes do not concern function definitions, it is enough to compile with the second command.

### 1.5.5 A typical runtime error NA STT error!

Open the `initMeth.c` file and find the line

```
ParxRelMakeNonEditable("Phase3DGradDur, ReadDepthTime");
```

Change slightly the name of the first parameter and recompile. The compiler will accept the change because the new name is still a valid string. However, when you restart ParaVision, it will be impossible to edit the `myGre` method (all parameters will disappear except for Method). The reason will be explained in the startup window: the parameter handling mechanism will not be able to change the attributes of a non-existing parameter.

**In conclusion**, whenever your method or method parameters disappear from the screen, this is due to an error in the method which has not been detected by the compiler. Look at the startup shell messages to figure out the exact reason.

**Finally, remember to remove all introduced errors!** You will need a working `myGre` for next projects. To be quite sure, remove `myGre` and generate it with `copyMethod` again.

## 2 Method Parameters: Adding a spoiler to gre (project spGre)

### 2.1 Objectives:

By developing this project you will learn how to


- define new parameters
- program a variable gradient pulse

### 2.2 Project description

We will start from the single gradient echo method (gre) and modify it to add a variable gradient spoiler in the slice gradient direction before the excitation RF pulse.

The most important modification consists of creating of two new parameters to control the amplitude and the duration of the spoiler. The repetition time of the sequence should be modified and a new gradient trim point will be created to define the spoiler amplitude.

## Method "spGre"



## 2.3 Specific method parameters

SliceSpoilGrad – spoiler amplitude in percent of the maximum gradient.

SliceSpoilGradDur – duration of the gradient spoiler in ms.

## 2.4 Realisation hints

### 2.4.1 Starting point

The method should be developed based on the *gre* method (discussed in the first project).

### 2.4.2 Steps

1. Use `copyMethod` to copy *gre* to *spGre*.
2. Define the new double parameters *SliceSpoilGrad* and *SliceSpoilGradDur* (in *SliceSelPars.h*). Declare their relations as *SliceSpoilGradRel* and *SliceSpoilGradDurRel*.
3. Add the the two parameters in the *SliceSelection* parameter class (in *SliceSelLayout.h*).

4. Write (in `SliceSel.c`) the `SliceSpoilGradRel` and the `SliceSpoilGradDurRel` relations. Each relation should consist of calling a range checking function for the respective parameter followed by a call of the backbone. Write these range-checking functions (`SliceSpoilGraRange` and the `SliceSpoilGradDurRange`). They should give the respective parameters a default value if they have none, and constrain them. The minimum value of the gradient spoiler duration should be  $2.0 * \text{CFG\_GradientRiseTime}()$ .
5. Set the visibility of the new parameters in the `SliceSelectionParsVisibility` relation (in `SliceSel.c`) to make them appear in the protocol.
6. Call the range checking functions of the new parameters in the initialisation of the `SliceSelection` parameter group (`InitSliceSelection`, in `SliceSel.c`).
7. Add the spoiler duration to the minimum repetition time calculation (in the function `UpdateRepetitionTime` located in `parsRelations.c`).
8. Add the new trim value `t9` for the `SliceSpoilGrad` parameter and the delay  $D[10] = (\text{SliceSpoilGradDur} - \text{riseTime}) / 1000$  in the function `SetBaseLevelParam` (in `BaseLevelRelations.c`).
9. Add the new line `d10 grad{(0)|(0)|(t9)}` in the `spGre.ppg` pulse program.

## 2.5 Proposed Solution

**SliceSelPars.h** define new parameters

```
double parameter
{
 display_name "Slice Spoiler";
 units "%";
 format "%f";
 relations SliceSpoilGradRel;
}SliceSpoilGrad;

double parameter
{
 display_name "Slice Spoiler Duration";
 units "ms";
 format "%f";
 relations SliceSpoilGradDurRel;
}SliceSpoilGradDur;
```

**SliceSelLayout.h** insert the new parameters in the `SliceSelection` class

```
parclass
{
 SliceGradStabTime;
 ExcSliceRephTime;
 ExcSliceGrad;
 ExcSliceGradLim;
 ExcSliceRephGrad;
 ExcSliceRephGradLim;
 SliceSpoilGrad;
 SliceSpoilGradDur;
```

```
}SliceSelection;
```

### SliceSel.c parameter relations

```
/*-----
 * Default relations and range checker for the spoiler
 -----/

void SliceSpoilGradRange()
{
 DB_MSG(("-->SliceSpoilGradRange"));

 if(!ParxRelsParHasValue("SliceSpoilGrad"))
 {
 SliceSpoilGrad = 40.0;
 }
 else
 {
 SliceSpoilGrad =
 MIN_OF(SliceSpoilGrad,100.0);
 SliceSpoilGrad =
 MAX_OF(SliceSpoilGrad,0.0);
 }
 DB_MSG(("<--SliceSpoilGradRange"));
}

void SliceSpoilGradRel(void)
{
 DB_MSG(("-->SliceSpoilGradRel"));

 SliceSpoilGradRange();
 backbone();

 DB_MSG(("<--SliceSpoilGradRel"));
}

void SliceSpoilGradDurRange()
{
 double min;

 DB_MSG(("-->SliceSpoilGradDurRange"));
 min = 2.0*CFG_GradientRiseTime();

 if(!ParxRelsParHasValue("SliceSpoilGradDur"))
 {
 SliceSpoilGradDur = MAX_OF(min,1.0);
 }
 else
 {
 SliceSpoilGradDur =
 MAX_OF(SliceSpoilGradDur,min);
 }
 DB_MSG(("<--SliceSpoilGradDurRange"));
}

void SliceSpoilGradDurRel(void)
{
 DB_MSG(("-->SliceSpoilGradDurRel"));

 SliceSpoilGradDurRange();
 backbone();
}
```

```

 DB_MSG(("<--SliceSpoilGradDurRel"));
}

```

### SliceSel.c initialization of the parameters

```

void InitSliceSelection(YesNo showAllPars)
{
 DB_MSG("-->InitSliceSelection");

 SliceGradStabTimeRange();
 ExcSliceRephTimeRange();
 ExcSliceGradRange();
 ExcSliceGradLimRange();
 ExcSliceRephGradRange();
 ExcSliceRephGradLimRange();
 SliceSpoilGradRange();
 SliceSpoilGradDurRange();

 SliceSelectionParsVisibility(showAllPars);

 DB_MSG("<--InitSliceSelection");
}

```

### parsRelations.c include SliceSpoilGradDur in the minimum repetition time

```

PVM_MinRepetitionTime =
 nslices *
 (
 0.011 +
 SliceGradStabTime +
 SliceSpoilGradDur +
 CFG_AmplifierEnable() +
 ExcPulse.Length/2 +
 PVM_EchoTime +
 PVM_AcquisitionTime *(1.0 - PVM_EchoPosition/100) +
 ReadSpoilGradDur +
 CFG_InterGradientWaitTime()
);

```

### BaseLevelRelations.c add new trim t9 and delay D[10]

```

ATB_SetAcqTrims(10,
 ExcSliceGrad, /* t0 */
 -ExcSliceRephGrad, /* t1 */
 -ReadDepthGrad, /* t2 */
 Phase2DGrad, /* t3 */
 -Phase3DGrad, /* t4 */
 ReadGrad, /* t5 */
 -Phase2DGrad, /* t6 */
 Phase3DGrad, /* t7 */
 ReadSpoilGrad, /* t8 */
 SliceSpoilGrad /* t9 */
);

```

(...)

```
D[10] = (SliceSpoilGradDur - riseT)/1000.0;
```

**SpGre.ppg**    add the line with d10

```
start, 10u fq8b:f1
 d10
 d1 fq1:f1 grad{(0) | (0) | (t9)}
 d2 gatepulse 1 grad{(0) | (0) | (t0)}
 (p0:sp0 ph1):f1
 ...
```

## 3 Timing Control: Long echo-time version of gre (project ItGre)

### 3.1 Objectives

By developing this project you will learn how to

- modify the timing of the pulseprogram
- control the TE and TR calculation of the method
- control the maximum gradient settings dependent on the pulse program timing
- test the correct method timing by dedicated subroutines


### 3.2 Project description

We will start from the known single gradient echo method (*gre*) and modify it to get a different timing of the measurement sequence. We will separate the slice gradient rephasing from the readout dephasing and phase-encoding (in *gre* these three gradient pulses are simultaneous) and generate a TE filling delay at a different position. In the new sequence the readout dephasing will be done immediately before the read gradient is switched on. This will reduce the flow artifacts at long echo times. The necessary modification of the pulse program consists of a separate gradient switching command for the slice refocusing gradient and the introduction of a new delay for a separate control the duration of the slice refocusing.

The additional delay introduced for slice refocusing must be taken into account in the calculation of the echo time and the repetition time in the method code (TR and TE are calculated by subroutines called in the backbone). The slice refocusing time must be made independent of the 2D phase-encoding time. Additionally, since only two gradient channels are switched on simultaneously during phase-encoding (and not three as in *gre*) the maximum gradient amplitudes may be increased for the 2D phase-encoding gradient. This can be made depending on the dimensionality of the experiment (for 3D acquisition the 3D phase-encoding gradient is switched simultaneously).



Method "ltGre"



### 3.3 Specific method parameters

ExcSliceRephTime – the duration of the slice refocusing pulse

Phase2DGradLim, Phase3DGradLim, ReadDephGradLim, ExcSliceRephGradLim:  
Maximum gradient amplitude for the 2D, 3D read diphase and excitation slice rephase gradient respectively.

### 3.4 Special features

None in this project.

### 3.5 Realisation hints

#### 3.5.1 Starting point

The method should be developed based on the *gre* method provided with Para Vision 3.0.2

#### 3.5.2 Steps

10. Use copyMethod to copy *gre* to *ltGre*.
11. Modify the pulse program *ltGre.ppg* (created by copy method). In a first step, copy the line “`d3 grad{ (t2) | r2d(t3) | (t1)+r2d(t4) }`” and paste it behind the “`d5 groff`” command. Remove the gradients which should not be switched at these two timepoints. Copy this pulse program to the pp directory and test it "manually".
12. Use a new delay (e.g. *d8*) for the slice rephasing gradient statement. The method *ltGre* should control this delay with parameter *ExcSliceRephTime*.
13. Check *ExcSliceRephTime* to appear editable in the method editor: The editability of parameters are controlled by special functions for the different parameter groups. In case of the slice selection group the routine *SliceSelectionParsVisibility* (c module *SliceSel.c*) is used. Inspect the way the parameter editability is handled.
14. Find the place in the backbone routine (*parsRelations.c*) that connects the *ExcSliceRephTime* with the *Phase2DgradDur* and modify this line appropriately.
15. (Optional) Modify the setting of *Phase2DGradLim*, *Phase3DGradLim*, *ReadDephGradLim*, *ExcSliceRephGradLim* (*initMeth.c* *parsRelations.c* – backbone). The maximum gradient amplitude has 2 major constraints:
  - The number of gradient pulses switched simultaneously at different gradient channels. Note, the norm of the vector switched in read, phase and slice direction must be less or equal to 1 (one) to prevent components for the x y and z gradient coil to be above 100% gradient amplitude in case of oblique slices.
  - The amplitude difference of two adjacent gradient switching events switched with the system gradient rise time must not exceed 100% for each channel (r,p,s).
  - Consider differences of the maximum Phase2D gradient dependent on the imaging dimensionality.

$$57\% = \frac{100}{\sqrt{3}} \quad \text{for oblique slices.}$$

16. Modify the UpdateEchoTime routine called in backbone to include the new timing. Check the TR calculation in UpdateRepetitionTime: Is it necessary to change it too?
17. Inspect the structure of SetBaseLevelParam routine (BaseLevelRelations.c) and find out which subroutine is used to set up the delays. Set delay D[8] to the appropriate value.
18. Inspect the routine PrintTimingInfo. Adapt this testing routine and activate debug messages for the routines in BaseLevelRelations.c.

### 3.6 Possible pitfalls

The gradient amplitudes may exceed 100% if the limits calculation is not appropriate.

### 3.7 Possible extensions

None.

### 3.8 Proposed Solution

ItGre.ppg (complete listing; modified/new elements are printed boldface)

```
#include<Avance.incl>
#include <DBX.include>
preset off

;=====
; definition of delays
;=====

define delay denab
"denab = d4 - de + depa"

"l3 = l0 + ds"

;=====
; declaration of 2d and 3d loop
;=====

lgrad r2d<2d> = L[0]
zgrad r2d
lgrad r3d<3d> = L[1]
zgrad r3d

if(DS>0)
{
 dsl, dgrad r2d
 lo to dsl times DS
}
```

```

}

lgrad slice = NSLICES
zslice

#include <MEDSPEC.include>

;=====
=
; D/P spec control gradients
;=====
=
start, 10u fq8b:f1
 d1 fq1:f1 grad{(0) | (0) | (t0)}
 d2 gatepulse 1
 (p0:sp0 ph1):f1
 d8 grad{(0) | (0) | (t1)}
 d5 groff
 d3 grad{(t2) | r2d(t3) |
r3d(t4)}
 denab REC_ENABLE grad{(t5) | (0) | (0)}
 ADC_INIT_B(ph1, ph0)
 aqq ADC_START
 d3 grad{(t8) | r2d(t6) |
r3d(t7)}
 d6 grad{(t8) | (0) | (0) }
 d7 groff
 d0 ADC_END
 1u islice
lo to start times NSLICES
 1u ipp1 zslice
lo to start times NA
 1u rpp1 igrad r2d
lo to start times l3 ;2d loop
 1u igrad r3d
 "l3=10"
lo to start times l1 ;3d loop
lo to start times NAE
goto start
exit
;=====
;phase lists

ph0 = 0
ph1 = 0 2 1 3
;=====

```

**initMeth.c** - Necessary modification in the initMeth() function:

```

ReadDepthGradLim =
 ExcSliceRephGradLim = ExcGradLim = 50.0;

```

(directly switching from positive to negative gradient, max. step is 100% during 1 risetime).

**parsRelations.c**

Necessary modification in the backbone routine:


```
(...)
ReadDepthTime = Phase3DGradDur = Phase2DGradDur;
(...)
dim = PTB_GetSpatDim();

if(dim == 3)
{
 Phase2DGradLim = Phase3DGradLim = 61.0 /* close to 0.5*sqrt(3/2) */
}
else
{
 Phase2DGradLim = 86.0; /* close to 0.5*sqrt(3/2) */
 Phase3DGradLim = 0.1; /* not used in this case but > 0 for safety
reasons */
}
```


Necessary modification in the UpdateEchoTime routine:

```
(...)
*mintel =
 ExcPulse.Length/2 +
 ExcSliceRephTime +
 igwt ;

*minte2 =
 riseTime +
 Phase2DGradDur +
 igwt ;
(...)
```



## 4 RF Pulses in PVM: A spin-echo method (project spinEcho)



### 4.1 Objectives

By developing this project you will learn the following topics:

- How to introduce new RF pulses into the sequence,
- How RF pulses are represented in PVM (struct RF\_PULSE\_PULSE\_TYPE)
- How to handle RF pulses with toolbox functions,
- How to use RF pulse enums (struct PV\_RF\_PULSE\_LIST) for an easy pulse shape selection,
- How to handle array parameters,
- How to calculate and use a new frequency list in the method code and in the pulse program.

### 4.2 Project description


We will create a new method spinEcho based on ItGre, a variant of gre in which the slice rephasing and readout de-phasing periods are separated (see previous project). This arrangement is convenient for introducing the refocusing RF pulse since it reduces signal losses due to diffusion and flow. The project is organized in a series of short steps which can be programmed and tested progressively.

First, we introduce a refocusing group consisting of a refocusing RF pulse and a slice gradient of the same amplitude as used for the excitation. The RF pulse will be surrounded by fixed spoiling delays. We will focus on the RF pulse handling in this step and neglect the modified timing for a moment. Our goal will be to familiarize with the RF pulse tools provided in PVM.

In the second step we will introduce variable spoiler delays and take care of the correct echo positioning. This part will show you how timing issues are handled in PVM methods.

Finally, as an optional extension, we propose to introduce individual slice-selection gradient amplitudes for both RF pulses. This will require two separate frequency offset lists for both pulses – an occasion to see how array parameters are used in the method code.

## Method "spinEcho"



## 4.3 Specific method parameters

### 4.3.1 New method parameters

**RfcPulse** – a struct parameter of predefined type to be defined in RFPulsePars.h (to be implemented in the first step).

**RfcPulseEnum** – a pulse list parameter of predefined type to facilitate the selection of available refocusing pulses, controls the name of the pulshape file.

**RfcSliceGrad** – a double parameter defining the refocusing slice selection gradient amplitude to be defined in SliceSelPars.h

**RefSliceSpoilTime** – a double parameter defining the time around the refocusing pulse to spoil unwanted coherences (FID of the excitation and refocusing pulse), to be defined in SliceSelPars.h

### 4.3.2 Important method parameters

**ExcPulse** – struct parameter already used to define the excitation pulse

**ExcPulseEnum** – a pulse list parameter of predefined type to facilitate the selection of available excitation pulses.

**BwScale** – a double parameter already defined to handle the transfer bandwidth of a multiple RF pulse slice selective preparation. Note that single RF pulse excitations as in case of the gradient echo variants of gre do not really need this parameter. It is predefined


Bandwidth factor of a pulse refers to FWHM of the shaped pulse

for educational purposes to facilitate the gradient calculations in case additional RF pulses are implemented.

#### 4.3.3 Important baselevel parameters

TPQQ – array of struct parameters to define RF pulse properties on the machine level

ACQ\_O1\_list ACQ\_O2\_list – array parameters holding the offset frequencies for slice selective excitation. Note that these parameters are bound to fq1 and fq2 commands in the pulseprogram.



#### 4.4 Important toolbox routines

The following toolbox routines are used to facilitate the handling of the different elements introduced in this project. The user may refer to the HTML toolbox documentation for detailed information:


STB\_UpdateRfcPulseName – used to update the RF pulse name according to the state of the refocusing pulse list.

STB\_InitRFPulse – type based initialisation routine for RF pulse struct parameters to assure valid ranges of the struct members.

#### 4.5 Realisation hints

##### 4.5.1 Starting point

The method should be developed based on the *ltGre* method.



##### 4.5.2 Refocusing pulse in the pulse program

Use copyMethod to copy *ltGre* to *spinEcho*. Modify the pulse program *spinEcho.ppg*:

- Introduce a new RF-Pulse command surrounded by two fixed delays of 1ms. Use P1 and sp1 and a new phase list ph2. Find out the correct place in the pulse program to do this.
- Write the phase list ph2. One simple possibility is to keep the refocusing pulse orthogonal to the excitation pulse. Remaining lists will need no modification in this case.
- During the first fixed delay switch on the slice gradient ( $t_0$ ). Switch it off after the second delay.
- Copy the pulse program into the standard PP directory
- Start ParaVision and create a scan with status “new”, select the new *spinEcho* method.
- Set the correct TPQQ entry to appropriate values in the ACQP editor and provide an appropriate pulse duration. Which TPQQ and P element needs to be modified?
- Start a GSP and move the TX sliders to get the signal and to verify it is a spin echo.
- Open the method editor again and have a look at the RF pulses parameter class: Open the RF pulse struct and have a look at the struct members. Modify the shapename and keep track to the reaction of the excitation pulse list parameter. You will access in the method code some of the struct elements. Modify the classification field and check what may be selected. This reaction is related to a switch statement in the range



checking routine of the ExcPulse parameter. You will need to modify this for the refocusing pulse.

#### 4.5.3 Refocusing pulse in the method

Shut down ParaVision and introduce the new RF pulse into the method code by performing the following steps:

- a. In RFPulsePars.h: Copy the ExcPulse parameter definition and call it RfcPulse. Remember the necessary renaming of display names and relations.
- b. Make the new RF pulse visible by extension of the RF\_pulses parameter class defined in RFPulseLayout.h.
- c. In RFPulses.c, write the relations of the new pulse (RfcPulseRel) by copy-and-paste of the existing ExcPulseRel function. Remember about renaming of "ExcPulse" to "RfcPulse" in the new function!
- d. Write the range checker (RfcPulseRange) by copy-and-paste of ExcPulseRange. Make sure this range checker is called in the relations of RfcPulse. Again, remember about renaming of ExcPulse to RfcPulse. Adapt the classification restriction (only LIB\_REFOCUS, PVM\_REFOCUS, and USER\_PULSE should be allowed) as well as the flip angle limit.
- e. Initialization RF Pulses Group (function InitRfPulses) – add a call of the range checker of RfcPulse. As you remember the range checkers also initialize parameters when they have no value.
- f. Update of RF Pulses Group (function UpdateRFPulse). Add another call of the toolbox function STB\_UpdateRfPulse with the new pulse as argument.
- g. Transfer the state of RfcPulse to baselevel parameters. Find in BaseLevelRelations the subroutine which sets pulse delays and TPQQ settings and add the proper entries for the refocusing pulse.

**Recompile** the method and start ParaVision. The refocusing pulse should be editable and the correct flip-angle/gains and bandwidth/duration calculations should take place. You can try the sequence with signal, however, the image will still have mixed properties of spin- and gradient-echoes because both echoes are not yet centred – this will be corrected later.

#### 4.5.4 Pulse Choice Enum (optional)

In addition to the struct parameter describing the refocusing pulse, a special enum parameter can be added to allow choosing between available pulse shapes (you can also select any shape by typing its name directly in the Fileneme field of the struct, but this is less convenient). Such an enum parameter already exists for the excitation pulse (ExcPulseEnum). We will add another one for the refocusing (RfcPulseEnum) by copying proper parts of the code:

- a. In RFPulsePars.h: Add a definition of RfcPulseEnum,
- b. In RFPulseLayout.h: Add the RfcPulseEnum to the definition of RF\_Pulses parclass.
- c. In RfPulses.c: Write relations of the new enum (RfcPulseEnumRel). In the relations, call STB\_UpdateRfcPulseName with proper arguments. This assures the change of the pulse's Filename when the enum is changed. Range-checking of the enum is not needed (because of the next step)
- d. In the function UpdateRFPulses (file RFPulses.c): call STB\_UpdateRfcPulseEnum with proper arguments. This changes the value of the enum when the Filename of the related pulse is changed.

**Recompile** the method and restart PV. Check the interaction between RfcPulse.FileName and RfcPulseEnum.

#### 4.5.5 Controlling the spin-echo position

Our next goal is to get the gradient echo in the centre of the spin echo and to make sure the echo time is correctly calculated. We will also introduce a parameter to control the duration of the spoilers on both sides of the refocusing pulse.

- a. in the pulse program use delays d9, d10, d11 in the refocusing group, as in the sequence diagram. The echo time and the centering will be controlled by d5 and d11.
- b. Introduce a new parameter "RfcSliceSpoilerDuration" as a member of the SliceSelection group. Give it a minimum of CFG\_GradientRiseTime() in its range check. It will affect delays d9 and d10. Call the range check in the InitRFPulses function.
- c. In parsRelations.c, edit the function UpdateEchoTime. It calculates two variables *minte1*, and *minte2*, which hold the sum of all fixed delays from the excitation to the refocusing, and from the refocusing to the echo, respectively. You should add the spoiler duration and one half of the refocusing pulse to both of them. Additionally, one rise time should be added to *minte2* as a minimum value of d11. Finally, change the restriction of the echo time (PVM\_EchoTime). The minimum is now given by two times *minte1* or two times *minte2* whichever one is longer. (Attention: these two parameters are passed by pointers; the asterisk is used to access their values).
- d. Verify the function UpdateRepetitionTime. Is any modification needed?
- e. Pass the values of *minte1* and *minte2* from the backbone via BaseLevelRelations() to SetPpgParameters(), where they will be needed to set delays d5 and d11. For this purpose, two double arguments must be added to the interfaces BaseLevelRelations and SetPpgParameters.
- f. In SetPpgParameters, calculate d5/d11 as half the difference between PVM\_EchoTime and *minte1/minte2*, plus rise time (remember about ms-to-sec conversion). Set d9 to RfcSliceSpoilerDuration and d10 to RfcSliceSpoilerDuration – CFG\_GradientRiseTime(). This apparent asymmetry will compensate the delay of the gradient system.

**Recompile** the method and restart PV. Check the reaction of the echo time to modifications of spoiler duration, acquisition size, pulse durations. Check in GSP mode whether the spin echo coincides with the gradient echo. For this purpose, turn the read gradient off by setting ACQ\_scaling\_read to zero. You should now see the pure spin-echo envelope and its maximum should correspond to the gradient echo centre (you may make it narrower with shims).

## 4.6 Possible extensions

### 4.6.1 Different amplitudes of excitation and refocusing slice gradients

So far, both RF pulses are applied with the same amplitude of the slice selection gradient. The effective slice thickness is determined by the pulse of a narrower bandwidth. In the next step we will adjust the slice selection gradient individually to each pulse so that they always act at the same slice thickness. As a consequence both pulses will also need individual frequency offset lists.

1. Implementation of slice gradient and slice spoiling: (stop ParaVision)
  - Define parameters RfcGrad in SliceSelPars.h. We use backbone as default relation of RfcGrad. Any input to this parameter should be replaced by a value calculated there. This is a common approach for parameters whose values are of interest but which should not be modified directly by user input.
  - Handle visibility of RfcGrad. The visibility of all slice selection parameters is handled by the locally defined routine SliceSelectionParsVisibility. Modify this routine appropriately
  - Consider the refocusing slice selection gradient for the calculation of the minimum slice thickness. You should:
    - provide an additional argument for routine SliceSelectionLimits, namely a pointer to an RF pulse (PVM\_RF\_PULSE\_TYPE \*const rfcPulse).
    - Declare a new double variable “min” to hold the minimum slice thickness with regard to the refocusing slice gradient
    - (optional) perform a trivial type based range check for the argument rfcPulse using STB\_CheckRFPulse. If called appropriately in backbone, the argument will not be changed. For safety reasons (to avoid division by zero for illegal pulse bandwidth) this call is nevertheless a good practice.
    - Calculate the minimum slice thickness so that RfcSliceGrad does not exceed 90% gradient amplitude using routine MRT\_MinSliceThickness (hint: the slice ratio is 1.0 in this example, why?). Store return value in variable “min”
    - Consider the “min” in the return value of SliceSelectionLimits
  - Calculate the refocusing slice gradient:
    - Define a new argument const double rfcPulseBW (effective bandwidth used for the slice gradient calculation) for routine UpdateSliceSelectionGradients. Check this argument to be greater than zero (0).
    - Use MRT\_SliceGrad to calculate RfcSliceGrad
  - Adapt backbone routine
    - Adapt the call of SliceSelectionLimits (think about the new argument RfcPulse!)
    - Adapt the call of UpdateSliceSelectionLimits rfcPulseBW (here you may consider the bandwidth scaling factor!)
  - Calculate an additional frequency list
    - Modify routine SetFrequencyParameters (BaseLevelRelations.c)
      - Change the dimension of ACQ\_O2\_list to the total number of slices (note the total number of slices has been calculated for the setting of the offsets for the slice excitation. Why do we have to set ACQ\_O2\_list (hint: see syntax in pulseprogram).
      - Calculate the ACQ\_O2\_list using the ACQ\_O1\_list. What is the relation between offsets for the excitation and for the refocusing slice? Use UT\_ScaleDoubleArray to do the job.

Recompile the method, restart ParaVision and check the extended version of spinEcho.

#### 4.7 Possible pitfalls

- Forgotten UT\_SetRequest in the relations of RF pulses – the backward requests (e.g. Bandwidth and Length, Attenuation and FlipAngle) will not be accepted. Bandwidth and Attenuation appear noneditable. The request mechanism will be explained in more detail in a later project (greTag).

- No signal for offcenter slice acquisition:
  - Wrong frequency list.
  - Frequency list not used (check pulseprogram)
- Wrong amplitudes for slice gradients:
  - error in handling of bandwidth scaling factor
  - wrong handling of minimum slice thickness: Correct return value of SliceSelectionLimits, Correct arguments for this routine (check call in backbone)
- Compiler errors concerning usage of unknown functions or wrong function prototypes:
  - Forgot to do a make cproto depend before the compile (especially if function arguments have been changed!)

## 4.8 Proposed Solution

**spinecho.ppg** (complete listing; modified/new elements are printed boldface)

```

#include<Avance.incl>
#include <DBX.include>
preset off

;=====
; definition of delays
;=====

define delay denab
"denab = d4 - de + depa"

"l3 = l0 + ds"

;=====
; declaration of 2d and 3d loop
;=====

lgrad r2d<2d> = L[0]
zgrad r2d
lgrad r3d<3d> = L[1]
zgrad r3d

if(DS>0)
{
 dsl, dgrad r2d
 lo to dsl times DS
}

lgrad slice = NSLICES
zslice

#include <MEDSPEC.include>

;=====
; D/P spec control gradients
;=====
start, 10u fq8b:f1
 dq1:f1 grad{(0) | (0) | (t0)}

```

```

d2 gatepulse 1
(p0:sp0 ph1):f1
d8
d5 groff fq2:f1 grad{(0) | (0) |(t1)}
d9 grad{(0) |(0) |(t9)}
d2 gatepulse (p1:sp1 ph2):f1
d10
d11 groff
d3 grad{(t2) | r2d(t3) | r3d(t4)}
denab REC_ENABLE grad{(t5) | (0) | (0)}
 ADC_INIT_B(ph1, ph0)
aqq ADC_START
d3 grad{(t8) | r2d(t6) | r3d(t7) }

d6 grad{(t8) | (0) | (0) }
d7 groff
d0 ADC_END
lu islice
lo to start times NSLICES
lu ippl zslice
lo to start times NA
lu rpp1 igrad r2d
lo to start times l3
lu igrad r3d
"l3=10"
lo to start times l1
lo to start times NAE
SETUP_GOTO(start)
exit
;=====
;phase lists

ph0 = 0
ph1 = 0 2 1 3
ph2 = 1 3 2 0 ← phase cycling to avoid imperfections in spin echo
;=====

```

RFPulsePars.h:

```

/*
 * refocusing pulse parameter
 */

PV_PULSE_LIST parameter
{
 display_name "Refocusing Pulse Choice";
 relations RfcPulseEnumRel;
}RfcPulseEnum;

PVM_RF_PULSE_TYPE parameter
{
 display_name "Refocusing Pulse";
 relations RfcPulseRel;
}RfcPulse;

```

RFPulseLayout.h:

```

parclass
{
 BwScale;

```

```

ExcPulseEnum;
ExcPulse;
RfcPulseEnum;
RfcPulse;
} RF_Pulses;

```

### RFPulses.c

```

void InitRFPulses(void)
{
 DB_MSG(("-->InitRFPulses"));

 BwScaleRange();
 ExcPulseRange();
 RfcPulseRange();
 DeriveGainsRange();
 STB_InitExcPulseEnum("ExcPulseEnum");
 STB_InitRfcPulseEnum("RfcPulseEnum");

 DB_MSG(("<--InitRFPulses"));

 return;
}

YesNo UpdateRFPulses(YesNo deriveGains, char *nucleus)
{
 YesNo referenceAviable;
 double referenceAttenuation=30;

 DB_MSG(("-->UpdateRFPulses"));

 if(deriveGains == Yes)
 referenceAviable =
 CFG_GetGlobRefAtt(nucleus, &referenceAttenuation);
 else
 referenceAviable = No;

 STB_UpdateRFPulse("ExcPulse",
 &ExcPulse,
 referenceAviable,
 referenceAttenuation);

 STB_UpdateExcPulseEnum("ExcPulseEnum",
 &ExcPulseEnum,
 ExcPulse.Filename,
 ExcPulse.Classification);

 STB_UpdateRFPulse("RfcPulse",
 &RfcPulse,
 referenceAviable,
 referenceAttenuation);

 STB_UpdateRfcPulseEnum("RfcPulseEnum",
 &RfcPulseEnum,
 RfcPulse.Filename,
 RfcPulse.Classification);

 DB_MSG(("<--UpdateRFPulses"));
}

```

```
 return referenceAviable;
}

(...)

void RfcPulseEnumRel(void)
{
 DB_MSG("-->RfcPulsesEnumRel");

 /* set the name and clasification of RfcPulse: */

 STB_UpdateRfcPulseName("RfcPulseEnum",
 &RfcPulseEnum,
 RfcPulse.Filename,
 &RfcPulse.Classification);

 /* call the method relations */
 backbone();

 DB_MSG("<--RfcPulseEnumRel");
}

void RfcPulseRange(void)
{
 DB_MSG("-->RfcPulseRange");

 if(ParxRelsParHasValue("RfcPulse") == No)
 {
 STB_InitRFPulse(&RfcPulse,
 "gauss.rfc",
 1.0,
 180.0);
 }

 /* allowed clasification */

 switch(RfcPulse.Classification)
 {
 default:

 RfcPulse.Classification = LIB_REFOCUS;
 break;
 case LIB_REFOCUS:
 case PVM_REFOCUS:
 case USER_PULSE:
 break;
 }

 /* allowed angle for this pulse */

 RfcPulse.FlipAngle = MIN_OF(180.0,RfcPulse.FlipAngle);
 RfcPulse.FlipAngle = MAX_OF(90.0,RfcPulse.FlipAngle);

 /* general verifiation of all pulse atributes */

 STB_CheckRFPulse(&RfcPulse);

 DB_MSG("<--RfcPulseRange");
}
```

```

void RfcPulseRel(void)
{
 DB_MSG("-->RfcPulseRel");

 /*
 * Tell the request handling system that the parameter
 * RfcPulse has been edited
 */

 UT_SetRequest("RfcPulse");

 /* Check the values of RfcPulse */

 RfcPulseRange();

 /*
 * call the backbone; further handling will take place there
 * (by means of STB_UpdateRFPulse)
 */

 backbone();

 DB_MSG("<--RfcPulseRel");
}

```

### SliceSelPars.h

```

double parameter
{
 display_name "Ref. Slice Gradient";
 relations backbone;
 units "%";
 format "%f";
}RfcSliceGrad;

double parameter
{
 display_name "Ref Slice Spoiling Time";
 relations RefSliceSpoilTimeRel;
 units "ms";
 format "%f";
}RefSliceSpoilTime;

```

### SliceSelLayout.h

```

parclass
{
 SliceGradStabTime;
 ExcSliceRephTime;
 RefSliceSpoilTime;
 ExcSliceGrad;
 RfcSliceGrad;
 ExcSliceGradLim;
 ExcSliceRephGrad;
 ExcSliceRephGradLim;
}SliceSelection;

```

### SliceSel.c:



```

void SliceSelectionParsVisibility(YesNo showAllPars)
{
 const char *const editable = "ExcSliceRephTime,"
 "RefSliceSpoilTime";

 const char *const nonedit = "ExcSliceGrad,"
 "SliceGradStabTime,"
 "ExcSliceGradLim,"
 "ExcSliceRephGrad,"
 "ExcSliceRephGradLim,"
 "RfcSliceGrad";

 DB_MSG(("-->SliceSelectionParsVisibility"));

 (...)

void InitSliceSelection(YesNo showAllPars)
{
 ...
 ExcSliceRephGradLimRange();
 RefSliceSpoilTimeRange(1.0);

 SliceSelectionParsVisibility(showAllPars);

 DB_MSG(("<--InitSliceSelection"));
}

double SliceSelectionLimits(PVM_RF_PULSE_TYPE *const excPulse,
 PVM_RF_PULSE_TYPE *const rfcPulse,
 const double gradStabTime,
 const double gradCalConst,
 double *const sliceRatio)

{
 (...)

 /* range check of arguments */

 STB_CheckRFPulse(excPulse);
 STB_CheckRFPulse(rfcPulse);

 (...)

 min = MRT_MinSliceThickness(rfcPulse->Bandwidth,
 1.0,
 90.0,
 90.0,
 gradCalConst);

 DB_MSG(("<--SliceSelectionLimits"));
 return MAX_OF(minSlThk,min);
}

YesNo UpdateSliceSelectionGradients(const double slthk,
 const double sliceRatio,
 const double excPulseBW,
 const double rfcPulseBW,
 double gradCalConst)

```

```

{
 DB_MSG(("-->UpdateSliceSelectionGradients"));

 if(slthk <= 0.0)
 {
 UT_ReportError("UpdateSliceSelectionGradients: "
 "Illegal value of argument 1\n");
 return No;
 }

 if(excPulseBW < 0.0 || rfcPulseBW <0.0)
 {
 UT_ReportError("UpdateSliceSelectionGradients: "
 "Illegal value of argument 2,3\n");
 return No;
 }

 (...)
 RfcSliceGrad = MRT_SliceGrad(rfcPulseBW,slthk, gradCalConst);
 (...)

 DB_MSG(("<--UpdateSliceSelectionGradients"));
 return Yes;
}

(...)

void RefSliceSpoilTimeRange(double gradstab)
{
 double min;
 DB_MSG(("-->RefSliceSpoilTimeRange"));

 min = CFG_GradientRiseTime();

 min += gradstab > 0.0 ? 0.0:gradstab;
 min = MIN_OF(min,10.0);

 if(!ParxRelsParHasValue("RefSliceSpoilTime"))
 {
 RefSliceSpoilTime = min;
 }
 else
 {
 RefSliceSpoilTime = MAX_OF(MIN_OF(RefSliceSpoilTime,10.0),min);
 }

 DB_MSG(("<--RefSliceSpoilTimeRange"));
}

void RefSliceSpoilTimeRel(void)
{
 DB_MSG(("-->RefSliceSpoilTimeRel"));

 RefSliceSpoilTimeRange(0.0);
 backbone();

 DB_MSG(("<--RefSliceSpoilTimeRel"));
}

```

**ParsRelations.c**

```
void backbone(void)
{
 (...)

 minSliceThick = SliceSelectionLimits(&ExcPulse,
 &RfcPulse,
 GradStabTime,
 PVM_GradCalConst,
 &sliceGradRatio);

 /*
 * update geometry parameters
 */
 (...)

 /*
 * calculate gradients in logical directions
 */

 FreqEncodingGradients(PVM_Fov[0], readGradRatio, PVM_GradCalConst);
 UpdatePhase2DGradients(PVM_Matrix[1], PVM_Fov[1], PVM_GradCalConst);

 if(dim == 3)
 {
 UpdatePhase3DGradients(PVM_Matrix[2], PVM_Fov[2], PVM_GradCalConst);
 }
 else
 {
 Phase3DGrad = 0.0;
 }

 UpdateSliceSelectionGradients(PVM_SliceThick,
 sliceGradRatio,
 ExcPulse.Bandwidth*BwScale/100.0,
 RfcPulse.Bandwidth*BwScale/100.0,
 PVM_GradCalConst);

 /*
 * calculate frequency offsets
 */

 LocalFrequencyOffsetRels();

 /*
 * update sequence timing
 */

 UpdateEchoTime(&mintel, &minte2);
 UpdateRepetitionTime();

 PVM_NEchoImages = 1;

 SetBaseLevelParam(mintel, minte2);
 SetRecoParam();
}
```

```

 DB_MSG(("<--backbone"));
 return;
}

void UpdateEchoTime(double *const mintel, double *const minte2)
{
 double riseTime, igwt;

 DB_MSG("-->UpdateEchoTime");

 riseTime = CFG_GradientRiseTime();
 igwt = CFG_InterGradientWaitTime();

 *mintel =
 ExcPulse.Length/2 +
 ExcSliceRephTime +
 igwt + /* min TE1/2 filling delay */
 RefSliceSpoilTime +
 RfcPulse.Length/2;

 *minte2 =
 RfcPulse.Length/2 +
 RefSliceSpoilTime +
 igwt + /* min TE2/2 filling delay */
 Phase2DGradDur +
 PVM_AcqStartWaitTime +
 PVM_AcquisitionTime * PVM_EchoPosition/100;

 PVM_MinEchoTime = 2*MAX_OF(*mintel, *minte2);
 PVM_EchoTime = MAX_OF(PVM_MinEchoTime, PVM_EchoTime);
 PVM_EchoTime1 = PVM_EchoTime2 = PVM_EchoTime;

 DB_MSG("<--UpdateEchoTime");
 return;
}

(...)

```

### BaseLevelRelations.c

```

void SetBaseLevelParam(double mintel, double minte2)
{
 DB_MSG("-->SetBaseLevelParam");

 SetBasicParameters();

 (...)

 SetPpgParameters(mintel, minte2);

 (...)
}

```

```

 DB_MSG(("<--SetBaseLevelParam"));
}

void SetFrequencyParameters(void)
{
 int nslices;

 DB_MSG(("-->SetFrequencyParameters"));

 ATB_SetNucl(PVM_Nucleus1);

 (...)

 nslices = GTB_NumberOfSlices(PVM_NSpacks, PVM_SPackArrNSlices);
 ATB_SetAcqO1List(nslices,
 PVM_ObjOrderList,
 PVM_SliceOffsetHz);

 PARX_change_dims("ACQ_O2_list",nslices);

 UT_ScaleDoubleArray(nslices,
 ACQ_O1_list,
 RfcSliceGrad/ExcSliceGrad,
 ACQ_O2_list);

 ATB_SetAcqO1BList(nslices,
 PVM_ObjOrderList,
 PVM_ReadOffsetHz);

 (...)

 DB_MSG(("<--SetFrequencyParameters"));
}

void SetGradientParameters(void)
{
 int spatDim, dim, i;

 DB_MSG(("-->SetGradientParameters"));

 (...)

 ATB_SetAcqTrims(10,
 ExcSliceGrad, /* t0 */
 -ExcSliceRephGrad, /* t1 */
 -ReadDepthGrad, /* t2 */
 Phase2DGrad, /* t3 */
 -Phase3DGrad, /* t4 */
 ReadGrad, /* t5 */
 -Phase2DGrad, /* t6 */
 Phase3DGrad, /* t7 */
 ReadSpoilGrad, /* t8 */
 RfcSliceGrad /* t9 */
);

 (...)

 DB_MSG(("<--SetGradientParameters"));
}

```

```

void SetPpgParameters(double minte1,double minte2)
{
 double riseT,igwT;
 int slices;
 DB_MSG(("-->SetPpgParameters"));

 (...)

 D[0] = ((PVM_RepetitionTime - PVM_MinRepetitionTime)/slices
 + igwT)/1000.0;
 D[1] = (SliceGradStabTime + riseT)/1000.0;
 D[2] = CFG_AmplifierEnable()/1000.0;
 D[3] = (Phase2DGradDur - riseT) / 1000.0;
 D[5] = (PVM_EchoTime/2 -minte1 + riseT + igwT)/1000.0;
 D[11] = (PVM_EchoTime/2 -minte2 + riseT + igwT)/1000.0;

 D[4] = (riseT + PVM_AcqStartWaitTime)/1000.0;
 D[6] = (ReadSpoilGradDur - Phase2DGradDur)/1000.0;
 D[7] = riseT/1000.0;
 D[8] = (ExcSliceRephTime - riseT)/1000.0;

 D[9] = RefSliceSpoilTime/1000.0;
 D[10]= (RefSliceSpoilTime-riseT)/1000.0;

 /* set shaped pulses */
 sprintf(TPQQ[0].name,ExcPulse.FileName);
 if(PVM_DeriveGains == Yes)
 {
 TPQQ[0].power = ExcPulse.Attenuation;
 }
 TPQQ[0].offset = 0.0;

 sprintf(TPQQ[1].name,RfcPulse.FileName);
 if(PVM_DeriveGains == Yes)
 {
 TPQQ[1].power = RfcPulse.Attenuation;
 }
 TPQQ[1].offset = 0.0;

 ParxRelParRelations("TPQQ",Yes);

 /* set duration of pulses */
 P[0] = ExcPulse.Length * 1000;
 P[1] = RfcPulse.Length * 1000;
 ParxRelParRelations("P",Yes);

 DB_MSG(("<--SetPpgParameters"));
}

```









---

# Method Programming in PVM

## Part 2

---

*ParaVision Programming Course*  
*April 3 – 7, 2006*

Authors:

Franciszek Hennel

Sascha Köhler

Markus Wick


Bruker BioSpin MRI GmbH





## Table of Contents

|                                                                           |           |
|---------------------------------------------------------------------------|-----------|
| <b>Introduction</b>                                                       | <b>5</b>  |
| <b>1 Using Modules: Saturation Slices (satGre)</b>                        | <b>6</b>  |
| 1.1 Objectives                                                            | 6         |
| 1.2 Project description                                                   | 6         |
| 1.3 The Saturation Slices Parameter Group                                 | 7         |
| 1.4 Toolbox functions associated to the Saturation Slices Parameter Group | 8         |
| 1.5 Pulse Program Subroutine SatSlices()                                  | 8         |
| 1.6 Realization Hints: Usage of Sat_Slices_Parameters parclass.           | 8         |
| 1.6.1 Creation of satGre                                                  | 8         |
| 1.6.2 Modification of Method Sources                                      | 9         |
| 1.7 Realization Hints: Modification of Global Module Functionality        | 9         |
| 1.7.1 Definition of a parameter controlling the orientation constrain     | 10        |
| 1.7.2 Modifications in backbone                                           | 10        |
| 1.8 Possible Pitfalls:                                                    | 11        |
| 1.9 Proposed Solution                                                     | 11        |
| <b>2 Creating modules: Tagging (greTag)</b>                               | <b>15</b> |
| 2.1 Objectives:                                                           | 15        |
| 2.2 Project description                                                   | 15        |
| 2.3 The greTag method                                                     | 16        |
| 2.4 Specific method parameters                                            | 17        |
| 2.5 Special features                                                      | 17        |
| 2.6 Your task                                                             | 18        |
| 2.7 Realisation hints                                                     | 18        |
| 2.7.1 Starting point                                                      | 18        |
| 2.7.2 Steps                                                               | 18        |
| 2.8 Proposed solution                                                     | 19        |
| <b>3 Accelerated Imaging (encGre)</b>                                     | <b>23</b> |
| 3.1 Objectives                                                            | 23        |
| 3.2 Project description                                                   | 23        |
| 3.3 Specific method parameters                                            | 24        |
| 3.3.1 New parameter group: Encoding                                       | 24        |
| 3.4 Realization hints                                                     | 25        |
| 3.4.1 Starting point                                                      | 25        |
| 3.4.2 Steps                                                               | 25        |
| 3.5 Proposed Solution                                                     | 27        |
| 3.6 Additional Exercise: Segmented acquisition                            | 32        |
| 3.6.1 Objectives                                                          | 32        |
| 3.6.2 Project description                                                 | 32        |
| 3.6.3 Realization hints                                                   | 33        |
| 3.6.4 Proposed Solution                                                   | 33        |



## Introduction

The first project for today (*satGre*) demonstrates the **usage of modules** in PVM. We will extend the gradient echo method *gre* by a series of slice-selective saturation RF pulses. Saturation slices are commonly used to avoid aliasing of signals from outside the field of view. Positions, orientation and thickness of the saturation slices will be controlled by the Geometry editor. All we will need to add this quite complex functionality to *gre* is to include a special module in the pulse program and provide its support in the PVM code. As you will see, a pulse program module comes along with a dedicated group of PVM parameters and with toolbox functions for initializing and updating of the group. A similar modularity is proposed in PVM for several other applications: the echo-planar imaging readout, diffusion weighting, water or fat suppression, and several other preparation sequences.

The goal of the next project (*greTag*) is to show how you can **develop your own PVM modules**. We will program a pulse program module for spin tagging and include it in the gradient echo sequence. The group of parameters controlling this module will be handled in a similar way as the predefined groups used for the saturation or geometry. We will see that this modular design of method elements is very useful when similar features need to be implemented in various methods.

As already mentioned, every parameter group has an updating function and we will also write one for our tagging module. The essential problem of such functions is to react properly to modifications of group members. When members A and B are to be related, should A be derived from B, or the other way round? This should be dependent on which parameter has last been modified by the user. We will make use of the **request handling mechanism** which provides exactly this type of information.

Finally, we will explore a new feature of PVM introduced in ParaVision 4.0 allowing a simple implementation of **accelerated imaging methods**. A variant of gradient echo will be programmed (*encGre*) using the *Encoding* parameter group to handle all common schemes of accelerated phase encoding, such as parallel imaging, partial-FT and zero-filling. We will also show how this parameter group can be used for a segmented k-space scanning.

# 1 Using Modules: Saturation Slices (satGre)

## 1.1 Objectives

By developing this project you will learn how to


- use a preparation module based on a globally defined parameter class
- modify the functionality of a preparation module
- handle multidimensional array parameters
- use a pulseprogram subroutine

## 1.2 Project description

The template gre will be copied to a new method satGre. First the saturation slices module will be introduced into the method. This module is implemented as a standard preparation module consisting of a dedicated set of PVM parameters and a parclass providing a user interface. In addition this module is connected to the routine geometry editor to assure a saturation slice placement based on reference images.

Parameter initialisation, -update and the setting of the baselevel parameters that are associated to the pulseprogram subroutine are based on toolbox functions. These functions are used also to control the visibility of the parameters as shown in the method editor and the geometry editor respectively.

In a second step the saturation slices are constrained to lie always parallel to the imaging slices. Here the gradient matrix (as defined in the standard slice geometry parameter group) is used to provide a gradient vector that is transferred to the related parameter of the saturation slices group.



### 1.3 The Saturation Slices Parameter Group

Parameters, that define the Saturation Slices are grouped in the Sat\_Slices\_Parameters parclass:

PVM\_FovSatNSlices: The number of saturation slices represented as int.

PVM\_FovSatSliceOrient: The major orientation of saturation slice (axial, sagittal, coronal) implemented as array of corresponding enums.

PVM\_FovSatThick: The thickness (defined by the FWHM of the excitation/inversion profile), implemented as double array with PVM\_FovSatNSlices entries.

PVM\_FovSatOffset: The distance to the gradient isocenter, implemented as double array with PVM\_FovSatNSlices entries.

PVM\_FovSatSliceVec: The slice orientation for each saturation slice, implemented as an array of 3 element vectors with PVM\_FovSatNSlices entries.

PVM\_SatSlicesPulseEnum: The pulse list parameter for the RF pulse used for slice selective saturation/inversion.

PVM\_SatSlicesPulse: The pulse struct parameter to define the specifications of the shaped RF pulse. Dependent on the Classification field it may be specified whether an excitation pulse (for saturation purposes, LIB\_EXCITATION) or a slice selective inversion pulse is used. The pulse list parameter will represent either a list of excitation pulses or inversion pulses according to the value specified in this field.

PVM\_SatSlicesDeriveGainMode: Dependent on this parameter, the derivation of a RF pulse amplitude is handled according to the state of PVM\_DeriveGains (YesNo parameter in the main

method parclass) or according to the availability of a reference attenuation to force RF amplitude derivation independent on the state of PVM\_DeriveGains.

PVM\_FovSatGrad: The amplitude of the slice gradient implemented as a double array with PVM\_FovSatNSlices entries.

PVM\_FovSatSpoilTime: Duration of the spoiling gradient used to suppress unwanted coherences invoked by the slice selective excitation/inversion.

PVM\_FovSatSpoilGrad: Gradient amplitude of the spoiling gradient.

PVM\_FovSatModuleTime: Total duration of the SatSlices module (needed to calculate the repetition time of the sequence).

The group handler of the saturation slices parameter group is PVM\_SatSlicesHandler.

Parameter PVM\_FovSatOnOff is used to switch saturation slices on or off. Only if set to On, the parameter group is visible in the method editor and saturation slices may be defined in the geometry editor.

#### 1.4 Toolbox functions associated to the Saturation Slices Parameter Group

Parameter initialisation: STB\_InitSatSlicesModule, used to provide legal starting values. This function is intended to be called in the initMeth routine of the method.

Parameter update: STB\_UpdateSatSlicesModule, used to provide a consistent state of the parameter group. This function must be called in the backbone of the method.

Baselevel parameter setting: ATB\_SetSatSliceBaseLevel, used to derive the baselevel parameters used in the pulse program part according to the state of the saturation slice parameter group.

#### 1.5 Pulse Program Subroutine SatSlices()

A new mechanism to handle common used functionality on pulse program level is the concept of subroutines. Subroutines are defined in the header of a pulse program and may be called inside the pulseprogram body. The related pulseprogram part for this project is available after an include of PrepModulesHead.mod in the pulseprogram header. This is a common include file for most of the preparation modules delivered with the software.

To apply the saturation slices, the subroutine SatSlices() must be called at a dedicated part of the major pulse program.

#### 1.6 Realization Hints: Usage of Sat\_Slices\_Parameters parclass.

##### 1.6.1 Creation of satGre

To perform the exercises create your version of method satGre and a reasonable protocol of this method by following these steps:



- Copy method gre from the /opt/PV4.0/prog/parx/src directory in your method development environment.
- Install the pulse program satGre.ppg as created by copyMethod into the pulse program directory of your ParaVision® installation (/opt/PV4.0/exp/stan/nmr/lists/pp).

### 1.6.2 Modification of Method Sources

The following approach is typical for the handling of modules defined by a global parameter group:

1. Define a parclass "Preparations" in parsLayout.h. Members of this class are PVM\_FovSatOnoff and the parclass Sat\_Slices\_Parameters. Place this parclass inside the methodClass. The parclass is now visible in the method editor.
2. Initialize the parameter group by calling STB\_InitSatSliceModule (please refer to the related online documentation) in initMeth (before the backbone call).
3. Update the parameter group with STB\_UpdateSatSlicesModule. This toolbox function requires the final value of PVM\_Nucleus1 (to derive the gyromagnetic ratio for gradient calculations) and must be called in the backbone routine (parsRelations.c) after the update of the nuclei. Since the minimum TR depends now on the duration of this module, the update must be called before the TR calculation. To prepare also the second part (modification of module functionality), the slice geometry parameters should have been updated before the call of STB\_UpdateSatSlicesModule.
4. Consider the duration of the module in the TR calculation (modification of local routine UpdateRepetitionTime).
5. Use ATB\_SetSatSlicesBaseLevel in routine SetBaseLevelParam to derive the acquisition parameters for the pulseprogram include.
6. Modify callbackDefs.h to assure a redirection of the group handler PVM\_SatSlicesHandler to the methods backbone routine.
7. Modify the pulseprogram satGre.ppg: include PrepModulesHead.mod and call the subroutine SatSlices in the pulseprogram body.
8. Test the method: Restart ParaVision, switch on saturation slices and place them inside the geometry editor. Acquire Images to assure the saturation is done as planned in the geometry editor

### 1.7 Realization Hints: Modification of Global Module Functionality

In this second part of the project the direction of the saturation slice vectors should be constrained to be along the orientation of the first slice package. This functionality should be dependent on a new (locally defined) parameter SatParallel of type YesNo. This parameter should be visible in the Preparations parclass only if the saturation slices module is switched on. Before the parameter group is updated, the last entry of the gradient matrix (parameter

PVM\_SpackArrayGradOrient[0][2]) is copied in each entry of PVM\_SatSliceVec in case SatParallel is set to Yes.

### 1.7.1 Definition of a parameter controlling the orientation constrain

- Define a YesNo parameter SatParallel (parsDefinition.h), the default relation should be SatParallelRel.
- Make it visible in the new Preparations subclass of the method (parsLayout.h)
- Write a standard range checking function SatParallelRange and the default relation SatParallelRel in parsRelations.c
- Call the range checking function in the default relation as well as in initMeth() of the method (force parameter initialization)

### 1.7.2 Modifications in the backbone

These modifications are related to the handling of multidimensional double array parameters. Parameter PVM\_SpackArrayGradOrient is implemented as double [i][j][k], where i is an index running from 0 to PVM\_NSpacks-1, j and k is running from 0 to 2. Whereas j relates to the read (j=0), phase (j=1) and slice (j=2) direction. The index k denotes the x (k=0), y (k=1) and z (k=2) component of the gradient vector in the patient coordinate system. The implementation of multidimensional arrays allows an easy access to the row vectors of the gradient matrix:

PVM\_SpackArrGradOrient[0][2] is the slice orientation vector of the first slice package. Parameter PVM\_FovSatSliceVec is implemented as double [i][j] array, where i is an index running from 0 to PVM\_FovSatNslices-1, j is the index for the components of the saturation slice vectors.

The c-function memcpy may be used to transfer the slice vector defined in the gradient matrix to the slice vectors of the saturation slices as follows:

```
memcpy(PVM_FovSatSliceVec[i],PVM_SpackArrGradOrient[0][2],3*sizeof(double));
```

The orientation transfer should be done before STB\_UpdateSatSlicesModule is called (why?):

1. Write an if-statement, testing whether SatParallel==Yes
2. Define in the body of the if-statement int variables to hold the actual dimension of the saturation slices gradient vector array (dim) as well as a running index i to implement a loop. Define a double array vec[3] to save the slice vector defined by the orientation matrix.
3. Copy the slice vector of the orientation matrix into the local variable vec (using memcpy).
4. Use PARX\_get\_dim("PVM\_FovSatSliceVec",1) to get the dimension of this array.
5. Loop over all vectors of PVM\_FovSatSliceVec and copy the orientation stored in the local variable vec into each vector.

6. Control the visibility of SatParallel according to the state of PVM\_FovSatOnOff. (use function ParxRelsHideInEditor)

### 1.8 Possible Pitfalls:

- The parameter group does not react on user input: Relation redirection in callbackDefs.h is missing?
- TR is not calculated correctly: Adaption of UpdateRepetitionTime is missing?
- Saturation is performed always perpendicular to the slice orientation: Is the correct index of PVM\_SpackArrGradOrient used to provide the saturation orientation?

### 1.9 Proposed Solution

(modifications in boldface):

#### File parsLayout.h:

```
parclass
{
 PVM_FovSatOnOff;
 SatParallel;
 Sat_Slices_Parameters;
}Preparation;

parclass
{
 Method;
 PVM_EchoTime;
 PVM_MinEchoTime;
 PVM_RepetitionTime;
 PVM_NAverages;
 PVM_ScanTimeStr;
 PVM_DeriveGains;
 RF_Pulses;
 Nuclei;
 SequenceDetails;
 Preparation;
 StandardInplaneGeometry;
 StandardSliceGeometry;
} MethodClass;
```

#### File parsDefinition.h :

```
YesNo parameter
{
 display_name "Sat. Parallel to Slice";
 relations SatParallelRel;
}SatParallel;
```

**File initMeth.c**

```
(...)
/*
 * initialize sat slices module
 */
SatParallelRange();
STB_InitSatSlicesModule();
```

**File callbackDefs.h**

```
/* saturation slices */
relations PVM_SatSlicesHandler backbone;
```

**File parsRelations.c**

```
void backbone(void)
{
 (...)

 /*
 * calculate frequency offsets
 */

 LocalFrequencyOffsetRels();

 /*
 * update saturation slices
 */

 if(PVM_FovSatOnOff==Off)
 {
 ParxRelsHideInEditor("SatParallel");
 }
 else
 {
 ParxRelsShowInEditor("SatParallel");
 }

 if(SatParallel==Yes)
 {
 int i,siz;
 double vec[3];

 memcpy(vec,PVM_SPackArrGradOrient[0][2],3*sizeof(double));
 dim=(int)PARX_get_dim("PVM_FovSatSliceVec",1);
 for(i=0;i<dim;i++)
 {
 memcpy(PVM_FovSatSliceVec[i],vec,3*sizeof(double));
 }
 }
}
```

```

 STB_UpdateSatSlicesModule(PVM_Nucleus1);

 (...)
}

void UpdateRepetitionTime(void)
{
 int nslices, dim;
 double TotalTime, mindur, riset;

 DB_MSG("-->UpdateRepetitionTime");
 (...)
 PVM_MinRepetitionTime =
 nslices *
 (
 0.011 +
 PVM_FovSatModuleTime +
 SliceGradStabTime +
 CFG_GradientRiseTime() +
 CFG_AmplifierEnable() +
 ExcPulse.Length/2 +
 PVM_EchoTime +
 PVM_AcquisitionTime * (1.0 - PVM_EchoPosition/100) +
 ReadSpoilGradDur +
 CFG_InterGradientWaitTime()
);
 (...)
}

void SatParallelRange(void)
{
 if(!ParxRelsParHasValue("SatParallel"))
 {
 SatParallel = No;
 }
 else
 {
 if(SatParallel != Yes && SatParallel != No)
 {
 SatParallel = No;
 }
 }
}

void SatParallelRel(void)
{
 SatParallelRange();
 backbone();
}

```

```
void SetBaseLevelParam()
{
 DB_MSG("-->SetBaseLevelParam");
 (...)

 ATB_SetSatSlicesBaseLevel();

 PrintTimingInfo();

 DB_MSG("<--SetBaseLevelParam");
}
}
```

## 2 Creating modules: Tagging (greTag)

### 2.1 Objectives:

By developing this project you will learn to

- program method elements in a modular manner,
- include a module in an existing method,
- use the Request Handling mechanism for a group of parameters.

### 2.2 Project description

We will develop a pulse program module for making a periodical pattern of stripes (tags) on the imaging object. Such a module may find application in dynamical studies of moving objects (e.g. in cardiac movies). The module will consist of a pair of identical non-selective RF pulses separated by a gradient pulse on the readout channel. Depending on the position along the readout axis, and on the precession caused by the gradient pulse, the flip angles of the RF pulses will either add or cancel. This will produce a periodic pattern of saturation tags perpendicular to the readout axis with the spacing given by the relation

$$\text{gradient\_amplitude} * \text{gradient\_time} * \text{tag\_spacing} = 1 \quad [1]$$

with the gradient amplitude expressed in Hz/mm. We will also add a spoiling gradient pulse to disperse the residual transverse magnetisation. The module will be inserted in the gre.ppg pulse program in front of the excitation pulse.





Fig. 1: Definition of Tag\_Spacing



**Fig. 2:** The tagging module.

The module will be controlled by a dedicated group of PVM parameters handled in the same way as the predefined parameter groups (Diffusion, contrast modules, etc.). This modular design of sequence elements has the advantage of being easy to transfer to different methods.

The parameters describing the gradient pulse and the tag spacing will be inter-related using the request handling mechanism. It will be possible to change the value of any of the three parameters and the method will react by adapting the remaining ones to keep the relation [1].

Your task will consist of modifying the *greTag* method, in which a simple form of the tagging module is already included. You will add a new parameter to the module and change the group behavior. The module will then be transferred to the *spinEcho* method.

### 2.3 The greTag method

The *greTag* method that you can find in the `/opt/PPC/PVM/src` directory already includes a simplified version of the Tagging module. The module allows specifying the gradient duration (*Tag\_GradientTime*) and the tag spacing (*Tag\_Spacing*) but the gradient amplitude is fixed at 5%. Your task will be to introduce a new parameter for the gradient amplitude and to produce specified changes in the module behavior.

The module is represented by the following files:

- The pulse program section (*Tag.mod*) to be included in the main ppg,



- TagPars.h containing the definitions of the module parameters,
- TagLayout.h containing the definition of the module parameters class Tagging,
- Tag.c containing the relations of the module parameters as well as the initialization, update and base-level-setting functions.

The \*.h files of the module are included in the corresponding \*.h of the method. The \*.c file is listed in the Makefile. You should remember this arrangement for transferring the module to a different method.

By working with this method you will notice that the parameters *Tag\_Spacing* and *Tag\_GradientTime* are related. This relation is programmed in the *Tag\_Update* function using the *UT\_GetRequest* function. This function allows checking which parameter has last been edited. The switch statement inside *Tag\_Spacing* provides different reactions to the editing of these parameters:

- **Tag\_Spacing changed:** Calculate *Tag\_GradientTime* correspondingly: if too short, restrict it and readjust *Tag\_Spacing*.
- **Tag\_GradientTime (or any other parameter) changed:** Calculate the corresponding *Tag\_Spacing*.

Note that the *UT\_GetRequest* function detects parameter changes only if the relations of these parameters include calls of *UT\_SetRequest*. You can also use this mechanism to detect edited structure fields and array elements, as described in ParaVision® Advanced Users Manual. This, however is out of scope of our project.

## 2.4 Specific method parameters

Parameters of the Tagging module:

|                    |                                                                                |
|--------------------|--------------------------------------------------------------------------------|
| Tag_RfPulse –      | RF pulse used (twice) in the module.                                           |
| Tag_Spacing –      | Tag pattern spacing (period) in mm.                                            |
| Tag_GradientTime – | Effective duration of the gradient pulse                                       |
| Tag_Gradient –     | Amplitude of the gradient pulse in % of max, to be introduced in this project. |
| Tag_ModuleTime –   | Duration of the tagging module                                                 |

## 2.5 Special features

When the tag spacing is set too narrow (below one pixel) the pattern disappears. The depth of the pattern (saturation strength) depends on the flip angle of the tag pulses, which should be typically in the range from 15 to 45 degree.

The value of *Tag\_GradientTime* must not be shorter than the gradient rise time.

The method must prevent zero value for any of the parameters involved in Eq. [1]: all of them are used in denominators in the code!

## 2.6 Your task

Review the code of the `greTag` method. In particular, look at the structure of `Tag_` parameters relations and at the usage of the request handling in the `Tag_Update` function.

Modify the `greTag` method (locally on your computer, not in PPC/PVM) in the following way:

- Add a new parameter `Tag_GradientTime` to control the amplitude of the tag gradient pulse
- Modify the `Tag_Update` function so that it reacts to the to the changes of all three parameters involved in Eq. [1] in the following way:
  - **Tag\_Gradient changed:** Calculate `Tag_GradientTime` to keep the current `Tag_Spacing`; if `Tag_GradientTime` becomes too short, restrict it and readjust `Tag_Spacing`.
  - **Tag\_Spacing changed:** Calculate `Tag_Gradient` to keep the current `Tag_Spacing`. Check the new value of `Tag_Gradient`. If too low or too high, restrict it and readjust `Tag_Spacing`.
  - **Tag\_GradientTime (or any other parameter) changed:** Calculate the corresponding `Tag_Spacing`.
- Transfer the modified tagging module to the `spinEcho` method (call it `spinEchoTag`).

## 2.7 Realisation hints

### 2.7.1 Starting point

The method should be developed based on the `greTag` method .

### 2.7.2 Steps

1. Copy `greTag` from the `/opt/methPool` to your local `methods/src` directory
2. Add a definition of `Tag_Gradient` in `TagPars.h`. Declare its relations as `Tag_GradientRels`.
3. Add `Tag_Gradient` to the `Tagging` parclass in `TagLayout.h`.
4. Write the relations (`Tag_GradientRels`) and range checker (`Tag_GradientRange`) of `Tag_Gradient`. You can take the ones of `TagGradientTime` as a basis. `Tag_Gradient` should be restricted to `[0.1 .. 100]`. Remember about `UT_SetRequest` in the relations.
5. Initialise `Tag_Gradient` by calling its range checker in `Tag_Init()`.
6. Modify `Tag_Update()`:
  - Replace the fixed amplitude (variable gradient) by the new parameter.
  - Add the parameter name to the list used by `UT_GetRequest`.

- Modify the switch statement: add a case for Tag\_Gradient (5, if you have placed it at the end of the list) and change the case for Tag\_Spacing, as described in 2.6.
7. Modify Tag\_BaseLevel() to set the gradient amplitude based on Tag\_Gradient.
  8. Test the modified greTag.
  9. **Optional:** You can restrict the tag spacing to at least twice the pixel resolution and at most a quarter of the FOV (so that the tagging pattern never disappears). This will require two additional arguments (fov, resolution) in the Tag\_Update function. There are different ways for this restriction. You can, for example, adjust the Tag\_Gradient to set the Tag\_Spacing within limits.
  10. Introduce the same module in spinEcho (call it spinEchoTag):
    - copy all Tag.\* files to spinEchoTag
    - add Tag\$(OBJEXT) to the OBJLIST definition in the Makefile
    - include Tag.mod in the pulse program in front of the excitation.
    - include the Tag\*.h files in the corresponding method's files.
    - call Tag\_Init() in initMeth()
    - call Tag\_Update() in the backbone() – when?
    - call Tag\_BaseLevel in SetBaseLevel().

## 2.8 Proposed solution

### File Tag.c (modifications in boldface)

```

/* relations of Tag_Gradient */
void Tag_GradientRels(void)
{
 DB_MSG(("--> Tag_GradientRels()"));

 UT_SetRequest("Tag_Gradient");
 Tag_GradientRange();
 ParxRelsParRelations("Tag_Handler", No);

 DB_MSG(("<-- Tag_GradientRels()"));
}

/* range check of Tag_Gradient */
void Tag_GradientRange(void)
{
 DB_MSG(("--> Tag_GradientRange()"));

 if(ParxRelsParHasValue("Tag_Gradient") == 0)
 Tag_Gradient = 5.0;

 Tag_Gradient = MAX_OF(Tag_Gradient, 0.1);
 Tag_Gradient = MIN_OF(Tag_Gradient, 100.0);

 DB_MSG(("<-- Tag_GradientRange()"));
}

```

```

/* -----
 * Tag_Init
 * initialisation of the Tag module
 *
 * ----- */
void Tag_Init(void)
{
 DB_MSG(("--> Tag_InitTagging()"));

 Tag_RfPulseRange();
 Tag_ModuleTimeRange();
 Tag_GradientTimeRange();
 Tag_GradientRange();
 Tag_SpacingRange();

 DB_MSG(("<-- Tag_InitTagging()"));
}

/* -----
 * Tag_Update
 * Updates the Tag module.
 * Arguments:
 * nucleus: name of the nucleus (PVM_Nucleus1)
 * deriveGains - if Yes, the gains for tagging pulses
 * will be derived.
 * ----- */
void Tag_Update(char *nucleus, YesNo deriveGains)
{
 YesNo referenceAviable;
 double referenceAttenuation=30;
 YesNo requestDetected;
 int dummy, parN;
 double gcc; /* gradient removed */
 const char *parList =
 "Tag_Spacing,"
 "Tag_GradientTime,"
 "Tag_RfPulse,"
 "Tag_ModuleTime,"
 "Tag_Gradient";

 DB_MSG(("--> Tag_Update()"));

 /* gradient calibration const, for later use */
 gcc = CFG_GradCalConst(nucleus);

 /* update the RF pulse */
 if(deriveGains == Yes)
 referenceAviable =
 CFG_GetGlobRefAtt(nucleus, &referenceAttenuation);
 else
 referenceAviable = No;
}

```

```

STB_UpdateRFPulse("Tag_RfPulse",
 &Tag_RfPulse,
 referenceAviable,
 referenceAttenuation);

/* react to requests */
/* we do not test structure fields, therefore "" and dummy: */
requestDetected = UT_GetRequest(parList, &parN, "", &dummy);

if(requestDetected == No)
 parN = 0; /* Lack of request is treated as default; see switch below */

DB_MSG(("request: %d",parN));

switch(parN)
{
 case 1: /* Tag_Spacing requested: Derive Tag_Gradient */
 Tag_Gradient = 1e5/(Tag_Spacing * Tag_GradientTime * gcc);
 if(Tag_Gradient < 0.1)
 {
 Tag_Gradient = 0.1;
 Tag_Spacing = 1e5 / (Tag_GradientTime * Tag_Gradient * gcc);
 }
 else if(Tag_Gradient > 100.0)
 {
 Tag_Gradient = 100.0;
 Tag_Spacing = 1e5 / (Tag_GradientTime * Tag_Gradient * gcc);
 }
 break;

 case 5: /* Tag_Gradient requested: Derive Tag_GradientTime */
 Tag_GradientTime = 1e5/(Tag_Spacing * Tag_Gradient * gcc);
 if(Tag_GradientTime < CFG_GradientRiseTime())
 {
 Tag_GradientTime = CFG_GradientRiseTime();
 Tag_Spacing = 1e5 / (Tag_GradientTime * Tag_Gradient * gcc);
 }
 break;

 case 2: /* Tag_GradientTime requested, */
 case 0: /* no request - */
 default: /* - derive Tag_Spacing: */
 Tag_Spacing = 1e5 / (Tag_GradientTime * Tag_Gradient * gcc);
}

/* module duration */
Tag_ModuleTime =
 0.01 +
 2 * Tag_RfPulse.Length +
 2 * CFG_AmplifierEnable() +
 Tag_GradientTime + /* tag gradient pulse */
 CFG_GradientRiseTime() +
 5.0 + /* spoiler */
 CFG_GradientRiseTime();

DB_MSG(("<-- Tag_Update()"));
}

/*-----

```

```

* Tag_BaseLevel
* Sets ACQP parameters used by the Tag module
-----/
void Tag_BaseLevel(void)
{
 /* delays */
 D[25] = CFG_AmplifierEnable()/1000.0;
 D[26] = Tag_GradientTime/1000.0;
 D[27] = CFG_GradientRiseTime()/1000.0;

 /* gradient amplitude */
 PVM_ppgGradAmp5 = Tag_Gradient/100.0;

 /* rf pulse length */
 P[5] = Tag_RfPulse.Length*1000.0;

 /* rf pulse */
 strcpy(TPQQ[5].name, Tag_RfPulse.Filename);
 if(PVM_DeriveGains == Yes)
 TPQQ[5].power = Tag_RfPulse.Attenuation;
 TPQQ[5].offset = 0.0;

 /* frequency list */
 ACQ_O2_list_size = -1;
 ParxRelsParRelations("ACQ_O2_list_size",Yes);
 ACQ_O2_list[0] = 0.0;
}

```

## 3 Accelerated Imaging (encGre)

### 3.1 Objectives

By developing this project you will learn how to

- implement a new parameter group called *Encoding*.

This new parameter group can handle all common encoding schemes including Parallel Imaging (PI), linear and centric k-space sampling, block-wise and interleaved segmentation, partial Fourier encoding and zero filling. Its purpose is to allow the user specify the required encoding scheme and to calculate the corresponding acquisition matrix size as well as the table of phase encoding steps.

With implementation of the *Encoding* group different mechanisms for accelerated imaging are provided.


### 3.2 Project description

We will start from the *gre* sequence and create a new method *encGre* using `copyMethod`. You will learn to perform the necessary steps to make the new parameter group visible in the *encGre* method with correct initialisation of all *Encoding* parameters. Afterwards you will learn how to update the *Encoding* parameters and how to set all relevant base-level parameters for controlling the phase encoding and the reconstruction of multi-channel acquisitions.

For experienced programmers there exists an additional exercise described in chapter 3.6. In this exercise you will use the implemented *Encoding* group to control the phase encoding schema of a segmented acquisition. For this purpose *PVM\_RareFactor* must be activated in the method and must be used as an update of *segment size* in the *Encoding* group.

### 3.3 Specific method parameters

#### 3.3.1 New parameter group: Encoding



- **Use Multiple Receivers (PVM\_EncUseMultiRec)** – A Yes/No parameter deciding whether the acquisition should work on a single or multiple channels. It can be switched to Yes only if the system contains multiple receive channels.
- **Active Receivers (PVM\_EncActiveReceivers)** – An array of On/Off parameters to activate/deactivate the data acquisition on selected channels.
- **PPI Acceleration (PVM\_EncPpiAccel1)** – An integer defining the acceleration by means of Parallel Imaging. This parameter introduces an undersampling of the k-space, i.e., it removes a periodic pattern of steps from the encoding scheme. The maximum value of PPI Acceleration is given by the number of active channels. Parallel Imaging preserves the resolution of the image but reduces the SNR by at least the square root of PVM\_EncPpiAccel1.
- **PPI Reference Lines (PVM\_EncPpiRefLines1)** – Number of nominally sampled lines in the k-space centre to be taken disregarding the PPI acceleration. The reference lines are needed for the auto-calibrating PI reconstruction methods, like the one currently implemented in ParaVision.
- **Zero-Fill Acceleration (PVM\_EncZfAccel1)** – A number telling how much the experiment is accelerated by a symmetric truncation of the k-space. Zero-filling should be used with caution since it reduces the effective resolution of the image (causes blurring). A similar parameter exists for the 3<sup>rd</sup> spatial dimension (Zero-Fill Acceleration 3d, PVM\_EncZfAccel2).



- **Phase Encoding Order (PVM\_EncOrder1)** – A symbol taking values ENC\_LINEAR or ENC\_CENTRIC and determining the order of acquisition of k-space lines (from edge to edge, or from centre outwards, respectively). The centric order allows achieving minimum effective echo times in multiple spin-echo methods such as RARE. It can also increase the signal-to-noise ratio in fast gradient echo methods (e.g. FISP) by effectively using the approach to the steady state. However, the centric scheme usually introduces a blurring effect.
- **Phase Encoding Start (PVM\_EncStart1)** – A number indicating the starting value of the phase encoding, ranging from -1 (lower edge of k-space) to +1 (upper edge). When the encoding start is higher than -1, the k-space is sampled cyclically. This parameter may be used to shorten the effective echo time in segmented experiments, such as RARE. This parameter is not available with the centric encoding and with partial-FT. A similar parameter exists for the third dimension (**Phase Encoding Start 3D**)
- **Acquisition Matrix Size (PVM\_EncMatrix)**. The effective size of the acquisition matrix taking into account the geometry parameters and all acceleration methods. Non-editable. When all acceleration parameters are set to 1.0 (no acceleration), the acquisition matrix size is a product of the image size (PVM\_Matrix, see In-plane geometry) and the anti-aliasing factor (PVM\_AntiAlias). Selecting one or more acceleration schemes (parallel, partial-FT, zero-filling) causes the acquisition matrix size to be reduced.
- **Encoding Output**. This subclass contains non-editable information parameters calculated by the Encoding mechanism.

### 3.4 Realization hints

#### 3.4.1 Starting point

The method should be developed based on the *gre* method.

#### 3.4.2 Steps

1. Use copyMethod to copy *gre* to *encGre*.
2. Like any other parameter group, *Encoding* must be included in the definition of MethodClass in the file **parsLayout.h** to be visible in the protocol.
3. In the **callbackDefs.h** file the group handler must be directed to the method's *backbone* function.
  - relations PVM\_EncodingHandler backbone;
4. The group should also be initialized in **initMeth.c** by call of
  - STB\_InitEncoding();
5. In the *backbone* function in **parsRelations.c** the *Encoding* group must be updated:

- `STB_UpdateEncoding(PTB_GetSpatDim(), /* dimension */  
                   PVM_Matrix,      /* image size */  
                   PVM_AntiAlias,   /* a-alias */  
                   &seg_size,      /* segment size */  
                   SEG_SEQUENTIAL,  /* segmenting mode */  
                   Yes,              /* PI in 2nd dim allowed */  
                   Yes, /* PI ref lines in 2nd dim allowed */  
                   No); /* partial ft in 2nd dim NOT allowed */`

Note that this update requires information about the dimension, image size and anti-aliasing and should therefore be made after the update of the in-plane geometry. The fourth argument of `STB_UpdateEncoding` is a pointer to a parameter defining the segment size in a multi-shot experiment, which typically is the RARE-factor (`PVM_RareFactor`). This parameter may change as a result of this function call. The next argument specifies the way of k-space segmenting. `SEG_SEQUENTIAL` means that consecutive scans should sample adjacent blocks of k-space lines. `SEG_INTERLEAVED` used for this argument would cause the consecutive k-space scans to be interleaved as in the RARE method. The final Yes/No arguments tell the group to activate or deactivate the possibility of PI, reference lines for PI, and partial (asymmetric) Fourier encoding.

Following the update, all parameters of the Encoding group have legal and consistent values.

6. Use `PVM_EncMatrix` for experiment time calculation in **parsRelations.c**, because `PVM_EncMatrix` describes the acquisition matrix and is defined as product of image size and anti-aliasing factor:
  - `TotalTime = PVM_EncMatrix[1] * PVM_RepetitionTime * PVM_NAverages;`
7. Use `PVM_EncMatrix` for setting `ACQ_size` in **BaseLevelRelations.c**
  - `ATB_SetAcqSize( Spatial, spatDim, PVM_EncMatrix, NULL, No );`
8. Set ACQP parameters for *User Encoding* in **BaseLevelRelations.c**
  - Use the switch statement: `switch(spatDim){....}`
9. Setting base-level parameters for controlling a multiple receiver acquisition:
 

This is done in the **BaseLevelRelations.c** file using the toolbox function `ATB_SetMultiRec()`. If this function returns *Yes*, a special 4-channel version of the pulse program needs to be declared. This takes place on older versions of the spectrometer hardware. Starting from Avance II this function always returns *No* and no special pulse program is needed.

  - `If(Yes == ATB_SetMultiRec())  
           ATB_SetPulprog("encGre.4ch");`
10. Multi receiver data require a special reconstruction mode which is defined in the **RecoRelations.c**:
  - To make things easier, copy `RecoRelations.c` from the master solution (`/opt/PPC/PVM/src/encGre/`) into your method directory.

### 3.5 Proposed Solution

#### parsLayout.h

```
parclass
{
 Method;
 PVM_EchoTime;
 PVM_MinEchoTime;
 PVM_RepetitionTime;
 PVM_NAverages;
 PVM_ScanTimeStr;
 PVM_DeriveGains;
 RF_Pulses;
 NuClei;
 Encoding;
 SequenceDetails;
 StandardInplaneGeometry;
 StandardSliceGeometry;
} MethodClass;
```

#### callbackDefs.h

```
/* *****
 * redirection of global parameter groups
 * *****/

/* Encoding */
relations PVM_EncodingHandler backbone;

/* inplane geometry */
relations PVM_InplaneGeometryHandler InplaneGeometryRel;
```

#### initMeth.c

```
STB_InitStandardInplaneGeoPars(2, dimRange, lowMat, upMat, No);

/*
 * init slice geometry parameters
 */
STB_InitSliceGeoPars(0, 0, 0);

/*
 * init encoding parameters
 */
STB_InitEncoding();

/*
 * init spectroscopy parameters (no csi)
 */
```

```
PTB_SetSpectroscopyDims(0, 0);
```

### parsRelations.c

```
STB_UpdateStandardInplaneGeoPars(minFov,2);

/* calculate gradients in logical directions */

if(dim == 3)
{
 /* Connect slice thickness to FOV in 3rd direction. */
 PVM_SliceThick = PVM_Fov[2];

 /* constrain maximum slices per package to 1 */
 STB_UpdateSliceGeoPars(0,0,1,minSliceThick);
}
else
{
 /* no constrain to slices in 2D mode */
 STB_UpdateSliceGeoPars(0,0,0,minSliceThick);
}
FreqEncodingGradients(PVM_Fov[0],readGradRatio,PVM_GradCalConst);
UpdatePhase2DGradients(PVM_Matrix[1],PVM_Fov[1],PVM_GradCalConst);

if(dim == 3)
{
 UpdatePhase3DGradients(PVM_Matrix[2],PVM_Fov[2],PVM_GradCalConst);
}
else
{
 Phase3DGrad =0.0;
}

/* once the image size is decided (after first update of inpl geo),
we can update the Encoding parameters, to get the acquisition size
(PVM_EncMatrix). */

int seg_size;
seg_size = 1;

STB_UpdateEncoding(PTB_GetSpatDim(), /* total dimensions */
 PVM_Matrix, /* image size */
 PVM_AntiAlias, /* a-alias */
 &seg_size, /* segment size */
 SEG_SEQUENTIAL, /* segmenting mode */
 Yes, /* ppi in 2nd dim allowed */
 Yes, /* ppi ref lines in 2nd dim allowed*/
 No); /* partial ft in 2nd dim allowed */

UpdateSliceSelectionGradients(PVM_SliceThick,
 sliceGradRatio,
 ExcPulse.Bandwidth*BwScale/100.0,
 PVM_GradCalConst);
```

```

(...)

PVM_RepetitionTime=MAX_OF(PVM_MinRepetitionTime,PVM_RepetitionTime);

dim = PTB_GetSpatDim();

TotalTime = PVM_RepetitionTime
 * PVM_EncMatrix[1]
 * PVM_NAverages;

if(dim == 3)
{
 TotalTime *= PVM_EncMatrix[2];
}

UT_ScanTimeStr(PVM_ScanTimeStr,TotalTime);
ParxRelShowInEditor("PVM_ScanTimeStr");
ParxRelMakeNonEditable("PVM_ScanTimeStr");

```

### BaseLevelRelations.c

```

SetMachineParameters();

if(PVM_ErrorDetected == Yes)
{
 UT_ReportError("SetBaseLevelParam: In function call!");
 return;
}

/* -----
 Sets parameters needed for multi-receiver acq. Overrides some
 previously set parameters such as NUCn Must be called at the end
 of SetBaseLevel.
 ----- */
if(Yes==ATB_SetMultiRec())
{
 ATB_SetPulprog("encGre.4ch");
}

PrintTimingInfo();

DB_MSG(("<--SetBaseLevelParam"));

(...)

/* ACQ_dim_desc */

ATB_SetAcqDimDesc(specDim, spatDim, NULL);
if(PVM_ErrorDetected == Yes)
{
 UT_ReportError("SetBasicParameters: In function call!");
 return;
}

```

```

 }

 /* ACQ_size */
 /* With the Encoding group, this call
 ATB_SetAcqSize(Spatial, spatDim, PVM_Matrix, PVM_AntiAlias, No);
 is replaced by: */
 ATB_SetAcqSize(Spatial, spatDim, PVM_EncMatrix, NULL, No);

 if(PVM_ErrorDetected == Yes)
 {
 UT_ReportError("SetBasicParameters: In function call!");
 return;
 }

 (...)

void SetGradientParameters(void)
{
 int spatDim, dim;

 DB_MSG(("-->SetGradientParameters"));

 ATB_SetAcqPhaseFactor(1);
 if(PVM_ErrorDetected == Yes)
 {
 UT_ReportError("SetGradientParameters: In function call!");
 return;
 }

 spatDim = PTB_GetSpatDim();

 dim = PARX_get_dim("ACQ_phase_encoding_mode", 1);
 PARX_change_dims("ACQ_phase_encoding_mode", spatDim);
 PARX_change_dims("ACQ_phase_enc_start", spatDim);
 switch(spatDim)
 {
 case 3:
 ACQ_phase_encoding_mode[2] = User_Defined_Encoding;
 ACQ_phase_enc_start[2] = -1; /* set, but no used */
 ACQ_spatial_size_2 = PVM_EncMatrix[2];
 ParxRelsCopyPar("ACQ_spatial_phase_2", "PVM_EncValues2");
 /* no break */
 case 2:
 ACQ_phase_encoding_mode[1] = User_Defined_Encoding;;
 ACQ_phase_enc_start[1] = -1.0; /* set, but no used */
 ACQ_spatial_size_1 = PVM_EncMatrix[1];
 ParxRelsCopyPar("ACQ_spatial_phase_1", "PVM_EncValues1");
 /* no break */
 default:
 ACQ_phase_encoding_mode[0] = Read;
 ACQ_phase_enc_start[0] = -1;
 }
}

```

```

ATB_SetAcqGradMatrix(PVM_NSpacks, PVM_SPackArrNSlices,
 PtrType3x3 PVM_SPackArrGradOrient[0],
 PVM_ObjOrderList);

if(PVM_ErrorDetected == Yes)
{
 UT_ReportError("SetGradientParameters: In function call!");
 return;
}

```

### RecoRelations.c

```

int dim,i;
int size,ftSize[3];

DB_MSG(("-->SetRecoParam"));

/* set baselevel reconstruction parameter */
/* default initialization of reco based on acqp pars allready set */

ATB_InitDefaultReco();

for(i=0; i<PTB_GetSpatDim(); i++)
 ftSize[i] = PVM_Matrix[i]*PVM_AntiAlias[i];

if(PVM_EncUseMultiRec == Yes)
{
 /* select method specific reconstruction method */
 RECO_mode = USER_MODE;
 ParxRelsParRelations("RECO_mode",Yes);
 ATB_InitUserModeReco(ACQ_dim, PVM_EncMatrix, ftSize, PVM_EncSteps1,
 PVM_EncSteps2,PVM_EncNReceivers, PVM_EncPpiAccell,
 PVM_EncPpiRefLines1, NI, ACQ_obj_order,
 ACQ_phase_factor, PVM_EchoPosition);
}

/* set reco rotate according to phase offsets */
dim = PTB_GetSpatDim();

/* set reco sizes and ft_mode (dim 2&3) */
/* (dim 1 is kept as it was set by ATB_InitDefaultReco) */
for(i=1; i<dim; i++)
{
 size = PVM_Matrix[i]*PVM_AntiAlias[i];
 RECO_ft_mode[i] = (size == PowerOfTwo(size)) ?
 COMPLEX_FFT:COMPLEX_FT;
 RECO_ft_size[i] = size;
 RECO_size[i] = PVM_Matrix[i];
}

ParxRelsParRelations("RECO_ft_mode",Yes);

```

```

ParxRelParRelations("RECO_ft_size", Yes);
ParxRelParRelations("RECO_size", Yes);

ATB_SetRecoRotate(PVM_EffPhase1Offset,
 PVM_Fov[1]*PVM_AntiAlias[1],
 NSLICES,
 PVM_NEchoImages,
 1) ; /* phase1 direction*/

(...)

ATB_SetRecoTranspositionFromLoops(PtrType3x3 ACQ_grad_matrix[0],
 NSLICES,
 1,
 NI,
 ACQ_obj_order);

RECO_bc_mode[0] = AUTO_OFFSET_BC;

DB_MSG(("<--SetRecoParam"));

}

int PowerOfTwo(int x)
/* returns next power of two */
{
 return (1<<(int)ceil(log((double)x)/log(2.0)));
}

```

## 3.6 Additional Exercise: Segmented acquisition

### 3.6.1 Objectives

- Use *PVM\_RareFactor* to define the segment size in the Encoding group.

The purpose of this exercise is to demonstrate the functionality of the new introduced *Encoding* group: only a few steps are necessary in order to implement the functionality of a segmented acquisition.

### 3.6.2 Project description

We will use *PVM\_RareFactor* as an argument of the update of the Encoding group. The fourth argument of *STB\_UpdateEncoding* is a pointer to a parameter defining the segment size in a multi-shot experiment, which typically is the *PVM\_RareFactor*.

To use *PVM\_RareFactor* in *encGre*, we have to make this parameter visible in *parsLayout.h*, to initialize it in *initMeth.c*, to use it as an update of the Encoding group in *parsRelations.c* and some base-level parameters must be set in *BaseLevelRelations.c*.



A new parameter *Seg\_DelayTime* will be defined in order to control segment repetition time.

### 3.6.3 Realization hints

1. Include *PVM\_RareFactor* in the definition of MethodClass in the file **parsLayout.h** to be visible in the protocol.
2. Initialize *PVM\_RareFactor* in **initMeth.c**:
  - `if(ParxRelsParHasValue("PVM_RareFactor") == No)`  
`PVM_RareFactor = 8;`
3. In **callbackDefs.h** the relation of *PVM\_RareFactor* must be redirected:
  - `relations PVM_RareFactor backbone;`
4. In **parsRelations.c** the Encoding group must be updated with a pointer argument of *PVM\_RareFactor* as segment size. Change the segmenting mode from `SEG_SEQUENTIAL` to `SEQ_INTERLEAVED`.
5. Define a new parameter *Seg\_DelayTime* in **parsDefinition.h** and make this parameter visible in the MethodClass in the file **parsLayout.h**.
6. *Seg\_DelayTime* should be restricted in **parsRelations.c**:
  - `Seg_DelayTime = MAX_OF(PVM_RepetitionTime, Seg_DelayTime);`
7. Set base-level parameters in the file **BaseLevelRelations.c**:
  - `ATB_SetAcqPhaseFactor( PVM_RareFactor);`
  - `ACQ_rare_factor = PVM_RareFactor;`
  - `L[0] = ACQ_size[1]/PVM_RareFactor;`
  - `L[4] = PVM_RareFactor;`
  - `D[8] = (Seg_DelayTime - PVM_RepetitionTime)/1000.0;`
8. Add *Seg\_DelayTime d8* and *loop l4* to the pulse program **encGre.ppg**.
9. Incorporate the *Seg\_DelayTime* into the calculation of the estimated total scan time (*PVM\_ScanTimeStr*) in **parsRelations.c**.

### 3.6.4 Proposed Solution

#### **parsLayout.h**

```
parclass
{
 Method;
 PVM_EchoTime;
 PVM_MinEchoTime;
```

```

PVM_RepetitionTime;
Seg_DelayTime;
PVM_NAverages;
PVM_ScanTimeStr;
PVM_DeriveGains;
RF_Pulses;
Nuclei;
Encoding;
PVM_RareFactor;
SequenceDetails;
StandardInplaneGeometry;
StandardSliceGeometry;
} MethodClass;

```

### initMeth.c

```

STB_InitEncoding();

if (ParxRelParHasValue("PVM_RareFactor") == No)
 PVM_RareFactor = 8;

/*
 * init spectroscopy parameters (no csi)
 */

```

### callbackDefs.h

```

/* *****
 *
 * relation redirection of single global parameters used in
 * this method:
 *
 * *****/

relations PVM_RepetitionTime RepetitionTimeRel;
relations PVM_NAverages AveragesRel;
relations PVM_EchoTime EchoTimeRel;
relations PVM_MinEchoTime backbone;

relations PVM_AcquisitionTime backbone;
relations PVM_DeriveGains DeriveGainsRel;
relations PVM_EffSWh EffSWhRel;
relations PVM_RareFactor backbone;

```

### parsRelations.c

```

/* once the image size is decided (after first update of inpl geo),
we can update the Encoding parameters, to get the acquisition size
(PVM_EncMatrix). */

```

```

STB_UpdateEncoding(PTB_GetSpatDim(), /* total dimensions */
 PVM_Matrix, /* image size */

```

```

PVM_AntiAlias, /* a-alias */
&PVM_RareFactor, /* segment size */
SEG_INTERLEAVED, /* segmenting mode */
Yes, /* ppi in 2nd dim allowed */
Yes, /* ppi ref lines in 2nd dim allowed */
No); /* partial ft in 2nd dim allowed */

```

(...)

```
PVM_RepetitionTime=MAX_OF(PVM_MinRepetitionTime,PVM_RepetitionTime);
```

```
Seg_DelayTime = MAX_OF(PVM_RepetitionTime, Seg_DelayTime);
```

```
dim = PTB_GetSpatDim();
```

```

TotalTime = PVM_RepetitionTime
 * PVM_EncMatrix[1]
 * PVM_NAverages;

```

### **parsDefinition.h**

```

int parameter
{
 display_name "Number of dummy scans";
 relations dsRelations;
} NDummyScans;

double parameter
{
 display_name "Segment Delay Time";
 units "ms";
 format "%.2f";
 relations backbone;
} Seg_DelayTime;

```

### **BaseLevelRelations.c**

```

void SetGradientParameters(void)
{
 int spatDim, dim;

 DB_MSG(("-->SetGradientParameters"));

 ATB_SetAcqPhaseFactor(PVM_RareFactor);
 if(PVM_ErrorDetected == Yes)
 {
 UT_ReportError("SetGradientParameters: In function call!");
 return;
 }
}

```

(...)

```

ACQ_scaling_read = 1.0;
ACQ_scaling_phase = 1.0;
ACQ_scaling_slice = 1.0;

ACQ_rare_factor = PVM_RareFactor;

ACQ_grad_str_X = 0.0;
ACQ_grad_str_Y = 0.0;
ACQ_grad_str_Z = 0.0;

```

(...)

```

sprintf(PULPROG, "encGre.ppg");

L[0] = ACQ_size[1]/PVM_RareFactor;
L[4] = PVM_RareFactor;

if(PTB_GetSpatDim() == 3)
{
 L[1] = ACQ_size[2];
}
else
{
 L[1] = 1;
}

slices = GTB_NumberOfSlices(PVM_NSpacks, PVM_SPackArrNSlices);
igwT = CFG_InterGradientWaitTime();
riseT = CFG_GradientRiseTime();

D[0] = ((PVM_RepetitionTime - PVM_MinRepetitionTime)/slices
 + igwT)/1000.0;
D[1] = (SliceGradStabTime + riseT)/1000.0;
D[2] = CFG_AmplifierEnable()/1000.0;
D[3] = (Phase2DGradDur - riseT) / 1000.0;
D[5] = (PVM_EchoTime - PVM_MinEchoTime + riseT + igwT)/1000.0;

D[4] = (riseT + PVM_AcqStartWaitTime)/1000.0;
D[6] = (ReadSpoilGradDur - Phase2DGradDur)/1000.0;
D[7] = riseT/1000.0;
D[8] = (Seg_DelayTime - PVM_RepetitionTime)/ 1000.0;

```

**encGre.ppg**

```

;=====
; definition of delays
;=====

define delay denab
"denab = d4 - de + depa"

"13 = 10 + ds"
"18 = 14 * ds"

```

```
=====
; declaration of 2d and 3d loop
=====

lgrad r2d<2d> = ACQ_size[1]
zgrad r2d
lgrad r3d<3d> = L[1]
zgrad r3d

if(DS>0)
{
 dsl, dgrad r2d
 lo to dsl times 18
}

(...)
lo to start times NSLICES
 lu ipp1 zslice
lo to start times NA
 lu rpp1 igrad r2d
lo to start times 13
 lu igrad r3d
 d8 ;Seg_DelayTime
lo to start times 14 ;PVM_RareFactor
lo to start times 11
lo to start times NAE
SETUP_GOTO(start)
exit
```





## Pipeline Filter

ParaVision Programming Course  
April 03-07, 2006

Arno Nauerth



## Pipeline Filter



- Introduction
- Overview of tutorial Examples
- Basics
- Examples: Step by Step
- Discussions




## Introduction: Usage



- Integration of an user program into an acquisition pipeline
- Sophisticated, user-defined raw data manipulations, such as
  - K-space corrections
  - K-space weighted accumulation
  - Keyhole imaging
  - Navigator techniques
  - ...




## Introduction: Acquisition Pipeline





## Introduction: User Filter Program



## Overview: Examples



- NoOperation – pass data without processing
- NoOpFast – optimized NoOperation
- SyncScan – Synchronized Experiment Start
- CartToPolar – Cartesian to Polar conversion of Scans
- KeyHole – Sequence acceleration
- MagControl – Navigator Corrected Magnetization



# Pipeline Filter: KeyHole






Image 1      Image 2      Image 3      Image 4      Image 5




# Pipeline Filter: MagControl



## Basics: Activating a Pipeline Filter



## Basics: Activating a Pipeline Filter



## Basics: General Form of a Pipeline Filter




```
Declarations (optional)
PIPE_INIT
PARX_INIT
Declarations (optional)
Program code (optional)
PIPE_PROC
Declarations (optional)
Program code (required)
PIPE_WRAPUP
Program code (optional)
PIPE_End
```



## Pipeline Filter: NoOperation




Open a ParaVision® Shell



Enter

```
cd $XWINNMRHOME/exp/stan/nmr/au/src
cp /opt/PPC/AU/src/NoOperation .
```

Edit AU NoOperation



## Pipeline Filter: NoOperation



```
PIPE_INIT
PARK_INIT
PIPE_PROC
int scan_cnt;
int dummy;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp; scan_cnt++)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy)
 memcpy(out_ptr, in_ptr, sizeof_scan);
 PIPE_SIGNAL_OUTPUT(1,1);
 }
 if (acq_scan_type == Scan_Experiment) break;
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
PIPE_END
```

Exit the Editor and let it be compiled




## Pipeline Filter: NoOperation



Load a Single Slice Gradient Echo Sequence with short TR


Edit ACQP Parameters



Enable Pipeline Filters



Activate NoOperation



Run Acquisition



## Pipeline Filter: NoOperation

BRUKER

Add a print statement to NoOperation

```
PIPE_INIT
PARX_INIT
PIPE_PROC
int scan_cnt;
int dummy;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp; scan_cnt++)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy)
 memcpy(out_ptr, in_ptr, sizeof_scan);
 PIPE_SIGNAL_OUTPUT(1,1);
 printf("%d scans transferred\n", scan_cnt);
 }
 if (acq_scan_type == Scan_Experiment) break;
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
PIPE_END
```

BRUKER  
BIOSPIN

## Pipeline Filter: NoOpFast

BRUKER


Enter within the  
ParaVision® Shell

```
{ cd $XWINNMRHOME/exp/stan/nmr/au/src }
cp /opt/PPC/AU/src/NoOpFast .
```

Activate NoOpFast

|                                 |                 |
|---------------------------------|-----------------|
| Add user's filter to pipeline ? | Yes             |
| User's filter program name      | NoOpFast        |
| Setup mode user's filters       | expand          |
| User's filter buffer memory     | For one PE step |
| User's filter sampling mode     | Standard        |

Edit AU NoOpFast



BRUKER  
BIOSPIN

## Pipeline Filter: NoOpFast

BRUKER

```
PIPE_INIT
PARX_INIT
PIPE_PROC
int scan_cnt;
int inCnt,outCnt,useCnt;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp;)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,inCnt,outCnt);
 useCnt = MIN_OF(inCnt,outCnt);
 memcpy(out_ptr, in_ptr, sizeof_scan*useCnt);
 PIPE_SIGNAL_OUTPUT(useCnt,useCnt);
 scan_cnt += useCnt;
 }
 if (acq_scan_type == Scan_Experiment) break;
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
PIPE_END
```

#scans input

#scans passed to

output(ACQ-out)

faster bc if more than 1 scan available,  
can transfer more than 1 scan.

BRUKER  
BIOSPIN

## Pipeline Filter: NoOpFast

BRUKER

```
FILE *fdTty=NULL;
PIPE_INIT
fdTty = fopen("/dev/tty","w");
if (NULL == fdTty) EXCEPT_printf("cannot fopen /dev/tty");
PARX_INIT
PIPE_PROC
int scan_cnt;
int inCnt,outCnt,useCnt;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp;)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,inCnt,outCnt);
 useCnt = MIN_OF(inCnt,outCnt);
 memcpy(out_ptr, in_ptr, sizeof_scan*useCnt);
 PIPE_SIGNAL_OUTPUT(useCnt,useCnt);
 scan_cnt += useCnt;
 fprintf(fdTty,"in=%d, out=%d, use=%d\n",inCnt,outCnt,
 useCnt);
 }
 if (acq_scan_type == Scan_Experiment) break;
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
if (NULL != fdTty) fclose(fdTty);
fdTty = NULL;
PIPE_END
```

BRUKER  
BIOSPIN

## Pipeline Filter: SyncScan

BRUKER

Enter within the  
ParaVision® Shell

```
cp /opt/PPC/AU/src/SyncScan .
pvWish
wm withdraw .
tk_messageBox -message Hello
exit
```

Activate NoOpFast for  
Setup an SyncScan for  
Scan experiments

|                                 |                                                                                                                                                                                                      |
|---------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Add user's filter to pipeline ? | Yes                                                                                                                                                                                                  |
| User's filter program name      | SyncScan                                                                                                                                                                                             |
| Setup mode user's filters       | <input checked="" type="checkbox"/> expand   p-1 < 2<br><input type="checkbox"/> [0] <input type="button" value="NoOpFast"/><br><input type="checkbox"/> [1] <input type="button" value="NoOpFast"/> |
| User's filter buffer memory     | For_one_PE_step                                                                                                                                                                                      |
| User's filter sampling mode     | Standard                                                                                                                                                                                             |

← traffic light Rxdj.  
won't do SyncScan.  
will do NoOpFast

BRUKER  
BIOSPIN

## Pipeline Filter: SyncScan

BRUKER

```
const char sysStr[] = {"echo \"wm withdraw .;tk_messageBox -message \

 \\\"Start Scan\\\"\" -type ok -icon info;exit\"|pvWish\"};
PIPE_INIT
PARX_INIT
system((const char *)sysStr);
PIPE_PROC
...
PIPE_WRAPUP
PIPE_END
```

prepares  
all pipelines

↑  
command prompt  
shell commands.

Use  and 

BRUKER  
BIOSPIN



## Pipeline Filter: CartToPolar

BRUKER

Load SETUP\_SHIM\_bas protocol, location B\_BASIC\_PVM

Enter within the  
ParaVision® Shell

cp /opt/PPC/AU/src/CartToPolar .

Activate CartToPolar  
for  
Setup Experiments

Add user's filter to pipeline? Yes

User's filter program name:  ← **CartToPolar**

Setup mode user's filters: # expand:  <1

User's filter buffer memory:

User's filter sampling mode:

BRUKER  
BIOSPIN

## Pipeline Filter: CartToPolar

BRUKER

```
double scale;
int dr;
PIPE_INIT
dr = GETIPAR("DR");
scale = (1<<(dr-1))/M_PI;
PARX_INIT
PIPE_PROC
int scan_cnt, dummy, *ip, *op;
int mag_int, phase_int;
double Re, Im;
while (1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp; scan_cnt++)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy)
 ip = (int *) in_ptr;
 op = (int *) out_ptr;
 for (i1=0; i1 < scan_word_size/2 ; i1++)
 {
 Re = *ip++ ;
 Im = *ip++ ;
 mag_int = sqrt(Re*Re + Im*Im);
 phase_int = scale * atan2(Im,Re);
 *op++ = mag_int;
 *op++ = phase_int;
 }
 PIPE_SIGNAL_OUTPUT(1,1);
 }
 if (acq_scan_type == Scan_Experiment) break;
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
PIPE_END
```

useful for  
preemph adj.

mag: ↓  
adj gain  
preemph


phase ↓  
adj B<sub>0</sub> comp.

BRUKER  
BIOSPIN

shift Scan by = -1 to remove dig filter

## Pipeline Filter: KeyHole

BRUKER



BRUKER  
BIOSPIN

## KeyHole: Parameters

BRUKER





**ArrayIndices[]** specifies array position for each Scan

**ACQ\_spatial\_phase\_1[]** specifies the K-space position for each Scan

arraySize = acq\_size[1] +  
(CenterProjections + NumberOfSegments) \* SegmentSize \* 2


BRUKER  
BIOSPIN




# MagControl: Results



ACQ Start at  $K = 0$ , Flip Angle 30 deg



MagControl<sup>led</sup>



uncontrolled





---

# Pipeline Filter

---

ParaVision Programming Course  
April 03-07, 2006

Author:  
Arno Nauerth

Bruker BioSpin MRI GmbH





## Table of Content

|            |                                          |           |
|------------|------------------------------------------|-----------|
| <b>1</b>   | <b>NoOperation</b>                       | <b>5</b>  |
| <b>1.1</b> | <b>Objectives</b>                        | <b>5</b>  |
| <b>1.2</b> | <b>Implementation</b>                    | <b>5</b>  |
| 1.2.1      | Parameters and their Settings            | 5         |
| 1.2.2      | Realization Hints                        | 5         |
| 1.2.3      | Creating the Pipeline Filter NoOperation | 6         |
| <b>1.3</b> | <b>Exercises</b>                         | <b>6</b>  |
| <b>2</b>   | <b>NoOpFast</b>                          | <b>6</b>  |
| <b>2.1</b> | <b>Objectives</b>                        | <b>6</b>  |
| <b>2.2</b> | <b>Implementation</b>                    | <b>7</b>  |
| 2.2.1      | Parameters and their Settings            | 7         |
| 2.2.2      | Realization Hints                        | 7         |
| 2.2.3      | The Pipeline Filter NoOpFast             | 7         |
| <b>2.3</b> | <b>Exercises</b>                         | <b>7</b>  |
| <b>3</b>   | <b>SyncScan</b>                          | <b>8</b>  |
| <b>3.1</b> | <b>Objectives</b>                        | <b>8</b>  |
| <b>3.2</b> | <b>Implementation</b>                    | <b>8</b>  |
| 3.2.1      | Parameters and their Settings            | 8         |
| 3.2.2      | Realization Hints                        | 8         |
| 3.2.3      | The Pipeline Filter SyncScan             | 9         |
| <b>3.3</b> | <b>Exercises</b>                         | <b>9</b>  |
| <b>4</b>   | <b>CartToPolar</b>                       | <b>9</b>  |
| <b>4.1</b> | <b>Objectives</b>                        | <b>9</b>  |
| <b>4.2</b> | <b>Implementation</b>                    | <b>10</b> |
| 4.2.1      | Parameters and their Settings            | 10        |
| 4.2.2      | Realization Hints                        | 10        |
| 4.2.3      | The Pipeline Filter CartToPolar          | 10        |
| <b>4.3</b> | <b>Exercises</b>                         | <b>11</b> |
| <b>5</b>   | <b>KeyHole Imaging</b>                   | <b>11</b> |
| <b>5.1</b> | <b>Objectives</b>                        | <b>11</b> |
| <b>5.2</b> | <b>Implementation</b>                    | <b>11</b> |
| 5.2.1      | Parameters and their Settings            | 11        |

|            |                                                                                      |           |
|------------|--------------------------------------------------------------------------------------|-----------|
| 5.2.2      | Realization Hints                                                                    | 11        |
| 5.2.3      | Creating the Measuring Method keyhole                                                | 12        |
| <b>5.3</b> | <b>Exercises</b>                                                                     | <b>19</b> |
| <b>6</b>   | <b>Magnetization Control</b>                                                         | <b>20</b> |
| <b>6.1</b> | <b>Objectives</b>                                                                    | <b>20</b> |
| <b>6.2</b> | <b>Implementation</b>                                                                | <b>20</b> |
| 6.2.1      | Parameters and their Settings                                                        | 20        |
| 6.2.2      | Realization Hints                                                                    | 20        |
| 6.2.3      | Creating the measuring method greMc                                                  | 21        |
| 6.2.4      | Exercises                                                                            | 25        |
| <b>7</b>   | <b>References</b>                                                                    | <b>25</b> |
| <b>7.1</b> | <b>Advanced Users Manual, Chapter D-9.3: AU filters for Pipelined Acquisition 25</b> |           |
| <b>7.2</b> | <b>Advanced Users Manual, Chapter D-9: Automation Programming</b>                    | <b>25</b> |
| <b>7.3</b> | <b>\$XWINNMRHOME/exp/stan/nmr/au/vorspann_f</b>                                      | <b>25</b> |
| <b>7.4</b> | <b>\$XWINNMRHOME/prog/include/aupipe.h</b>                                           | <b>25</b> |



# Pipeline Filter

## 1 NoOperation

*Pipeline Filter which passes through all Scans without any processing of them.*

### 1.1 Objectives

The purpose of this Pipeline Filter is to

- understand the basic principle of a Pipeline Filter
- have a simple template for the development of a new Pipeline Filter which does neither a data reduction nor a data expansion
- have a no-operation pipeline filter wherever it is needed

### 1.2 Implementation

#### 1.2.1 Parameters and their Settings

- ACQ\_user\_filter = Yes
- ACQ\_user\_filter\_name = ACQ\_user\_filter\_setup\_name = NoOperation
- ACQ\_user\_filter\_memory = For\_one\_scan
- ACQ\_user\_filter\_mode = Standard

#### 1.2.2 Realization Hints

- The predefined integer variable **scans\_per\_exp** is initialized to the total number of scans expected by the pipeline filters succeeding the user filter
- **PIPE\_WAIT\_INPUT\_OUTPUT(x1, x2, y1, y2)** – Wait for at least **x1** scans to be available in the input buffer, where at least **y1** scans are contiguous, and wait until there is space in the output buffer for at least **x2** scans, where at least **y2** scans are contiguous.
- After the return from **PIPE\_WAIT\_INPUT\_OUTPUT()**, the predefined char pointers **in\_ptr** and **out\_ptr** are pointing to the input buffer and the output buffer respectively
- **PIPE\_SIGNAL\_OUTPUT(x1,x2)** – Signal the system that **x1** input scans and **x2** output scans have been processed
- **void \*memcpy(void \*dest, const void \*src, size\_t n)** – copy **n** bytes from memory area **src** to memory area **dest**

### 1.2.3 Creating the Pipeline Filter NoOperation

Use the Menu function *Edit AUs* of the *Spectrometer Control Tool* in order to create and compile the program to be executable within the acquisition pipeline environment.

```
PIPE_INIT
PARK_INIT
PIPE_PROC
int scan_cnt;
int dummy;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp; scan_cnt++)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy)
 memcpy(out_ptr, in_ptr, sizeof_scan);
 PIPE_SIGNAL_OUTPUT(1,1);
 }
 if (acq_scan_type == Scan_Experiment) break;
} /* while (1) */
PIPE_WRAPUP
PIPE_END
```

The source of NoOperation will be found within the directory `$XWINNMRHOME/exp/stan/nmr/au/src`. During the development of a Pipeline Filter, it is often more efficient to have the source code permanently open with an arbitrary editor. When doing so, any change to the file needs to be manually compiled, e.g. by

```
$XWINNMRHOME/exp/stan/nmr/au/makeau NoOperation
```

## 1.3 Exercises

- comment out the `memcpy()` function
- add `printf()` statements for test purposes

## 2 NoOpFast

*Variation of the Pipeline Filter NoOperation for maximum throughput*

### 2.1 Objectives

The purpose of this Pipeline Filter is to

- optimize the data throughput for experiments with an very high effective sampling rate
- have a simple template for an optimized data throughput

- have an efficient no-operation pipeline filter wherever it is needed

## 2.2 Implementation

### 2.2.1 Parameters and their Settings

- ACQ\_user\_filter = Yes
- ACQ\_user\_filter\_name = ACQ\_user\_filter\_setup\_name = NoOpFast
- ACQ\_user\_filter\_memory = {For\_one\_PE\_step, For\_one\_Experiment}
- ACQ\_user\_filter\_mode = Standard

### 2.2.2 Realization Hints

The macro *MIN\_OF(a,b)* returns the minimum value of the two operands *a* and *b*.

### 2.2.3 The Pipeline Filter NoOpFast

The modified lines to NoOperation are printed in bold face:

```

PIPE_INIT
PARX_INIT
PIPE_PROC
int scan_cnt;
int inCnt,outCnt,useCnt;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp;)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,inCnt,outCnt)
 useCnt = MIN_OF(inCnt,outCnt);
 memcpy(out_ptr, in_ptr, sizeof_scan*useCnt);
 PIPE_SIGNAL_OUTPUT(useCnt,useCnt);
 scan_cnt += useCnt;
 }
 if (acq_scan_type == Scan_Experiment) break;
} /* while (1) */
PIPE_WRAPUP
PIPE_END

```

## 2.3 Exercises

- use FILE \*fdTty, fdTty=fopen("/dev/tty","w") and fprintf(fdTty,"...") in order to get a kind of real-time debug information within the ParaVision startup terminal window
- modify sequence in order to have a high effective sampling rate
- set NI > 1 and ACQ\_user\_filter\_size = For\_one\_PE\_step

## 3 SyncScan

*Pipeline Filter which raises a top-level window to control the start of the acquisition (modified NoOpFast)*

### 3.1 Objectives

The purpose of this Pipeline Filter is to

- understand the synchronization of the initialization- proceccing- and wrapup sections of Pipeline Filters respectively
- usage of `system()` commands within the different sections

### 3.2 Implementation

#### 3.2.1 Parameters and their Settings

- `ACQ_user_filter = Yes`
- `ACQ_user_filter_name = ACQ_user_filter_setup_name = SyncScan`
- `ACQ_user_filter_memory = x`
- `ACQ_user_filter_mode = Standard`

#### 3.2.2 Realization Hints

- `system(const char *string)` – executes a command in `string` by calling `/bin/sh -c string`
- `pvWish` – TclTk window shell for *ParaVision*<sup>®</sup>
  - `tk_messageBox` – display a toplevel message box within `pvWish`
  - enter the follwing commands within a *ParaVision*<sup>®</sup> shell

```
pvWish
wm withdraw .
tk_messageBox -message {Hello World}
exit
```

### 3.2.3 The Pipeline Filter SyncScan

```

const char sysStr[] = {"echo \\"wm withdraw .;tk_messageBox -message \
 \\\\\"Start Scan\\\\\\" -type ok -icon
info;exit\\"|pvWish"};

PIPE_INIT
PARX_INIT

system((const char *)sysStr);
PIPE_PROC
int scan_cnt;
int inCnt,outCnt,useCnt;
int *ip, *op;
while(1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp;)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,inCnt,outCnt)
 useCnt = MIN_OF(inCnt,outCnt);
 memcpy(out_ptr, in_ptr, sizeof_scan*useCnt);
 PIPE_SIGNAL_OUTPUT(useCnt,useCnt);
 scan_cnt += useCnt;
 }
 if (acq_scan_type == Scan_Experiment) break;
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
PIPE_END

```

### 3.3 Exercises

- Add a system call to signalize the end of the experiment
- What does it mean to put a system call outside the *while (1) {...}* body but before the *PIPE\_WRAPUP* statement?

## 4 CartToPolar

*Reorganize raw data from Cartesian Real/Imaginary pairs into Magnitude/Phase pairs.*

### 4.1 Objectives

The purpose of this Pipeline Filter is to

- get a better idea of the power of Pipeline Filters
- have a PREEMPHASIS tool available for separating gradient from B0 eddy current effects

## 4.2 Implementation

### 4.2.1 Parameters and their Settings

- ACQ\_user\_filter = Yes
- ACQ\_user\_filter\_name = ACQ\_user\_filter\_setup\_name = CartToPolar
- ACQ\_user\_filter\_memory = x
- ACQ\_user\_filter\_mode = Standard

### 4.2.2 Realization Hints

- **M\_PI** is a constant definition for  $\pi$
- **GETIPAR("parName")** – returns the value of the PARX integer parameter **parName**
- **double atan2(double Y, double X)** – calculates the arc tangent of the two variables **X** and **Y**

### 4.2.3 The Pipeline Filter CartToPolar

```

double scale;
int dr;

PIPE_INIT
dr = GETIPAR("DR");
scale = (1<<(dr-1))/M_PI;
PARX_INIT
PIPE_PROC
int scan_cnt, dummy, *ip, *op;
int mag_int, phase_int;
double Re, Im;

while (1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp; scan_cnt++)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy)
 ip = (int *) in_ptr;
 op = (int *) out_ptr;
 for (i1=0; i1 < scan_word_size/2 ; i1++)
 {
 Re = *ip++ ;
 Im = *ip++ ;
 mag_int = sqrt(Re*Re + Im*Im);
 phase_int = scale * atan2(Im,Re);
 *op++ = mag_int;
 *op++ = phase_int;
 }
 PIPE_SIGNAL_OUTPUT(1,1);
 }
 if (acq_scan_type == Scan_Experiment) break;
}

```

```

scans_in = 0;
} /* while (1) */
 PIPE_WRAPUP
PIPE_END

```

### 4.3 Exercises

- Remove the scaling by the variable `scale`

## 5 KeyHole Imaging

*Measuring method keyHole using a Pipeline Filter for speeding up an imaging sequence*

### 5.1 Objectives

The purpose of this Pipeline Filter is to

- understand the pipeline mechanism for data reduction or data expansion  
↳ needed to fill in the k-space lines that we kept from the previous scans (i.e. Centre projections)

### 5.2 Implementation

#### 5.2.1 Parameters and their Settings

- `ACQ_user_filter` = Yes
- `ACQ_user_filter_name` = `ACQ_user_filter_setup_name` = KeyHoleAu
- `ACQ_user_filter_memory` = x
- `ACQ_user_filter_mode` = Special
- `ACQ_user_filter_size[0]` = `ACQ_size[0]`
- `ACQ_user_filter_size[1]` < `scans_per_exp`

#### 5.2.2 Realization Hints

- `ACQ_user_filter_size[0] = 0` is a substitution for `ACQ_size[0]`
- `ACQ_phase_encoding_mode[1] = User_Defined_Encoding` can be used to have a user defined r2d phase encoding function
  - `ACQ_spatial_size_1` defines the size of the gradient function `r2d`
  - `ACQ_spatial_phase_1[]` contains the phase encoding values for `r2d`
- Apply the command `copyMethod` within a *ParaVision*<sup>®</sup> shell in order to create a private measuring method
- The environment variable `ParxMethodSearchPath` can be used to invoke private user-created measuring method
- `PIPE_WAIT_OUTPUT(x,y)` Wait for at least x scans to be available in the output buffer, where at least y scans are contiguous

## 5.2.3 Creating the Measuring Method keyhole

### 5.2.3.1 Using *copyMethod* to copy *gre* to *keyHole*

```
>copyMethod ($XWINMRHOME/prog/service/copyMethod)
Enter name of source method : gre
Enter name of target method : keyHole
Copy source method e.g. to /home/$USER/pv302/methods
>cd /home/$USER/pv302/methods/keyHole
```



### 5.2.3.2 Modification of the Pulse Program

There is no modification of the pulse program required when using `ACQ_phase_encoding_mode[1] = User_Defined_Encoding`.

### 5.2.3.3 Extension to `parsDefinition.h`

```

/*****
 * PARAMETER DEFINITION FOR KEY-HOLE
 *****/
int parameter
{
 display_name "nuber of segments";
 relations KeyHoleRelations;
}NumberOfSegments;
int parameter
{
 display_name "segment size";
 relations KeyHoleRelations;
}SegmentSize;
int parameter
{
 display_name "number of center projections";
 relations KeyHoleRelations;
}CenterProjections;
int parameter
{
 display_name "array indices";
}ArrayIndices[];

```

### 5.2.3.4 Extension to `parsLayout.h`

```

/*-----*
 * PV class...Extension for KeyHole
 -----/

parclass
{
 NumberOfSegments;
 CenterProjections;
 SegmentSize;
 ArrayIndices;
}KeyHole;

parclass
{
 Method;
 KeyHole;
...
}MethodClass;

```

### 5.2.3.5 Extension of *initMeth.c*

```
int dimRange[2] = {2,2}; //limit experiment to two dimensions
```

### 5.2.3.6 Extension of *parsRelations.c*

```
/*:=MPB:=*****
 *
 * Global Function: KeyHoleRelations()
 * Description: calculate gradient- and index array for Keyhole
 * imaging.
 * Preconditons: - linear phase encoding, starting at -1.0 */

void KeyHoleRelations (void)

/*:=MPE:=******/
{
 int i,j, idx;
 int arraySize;
 int matrixSize = ACQ_size[1];

 DB_MSG("-->KeyHoleRelations");
 if (NumberOfSegments < 1)
 NumberOfSegments = 2;
 if (NumberOfSegments >= matrixSize/4)
 NumberOfSegments = matrixSize/4;
 if (CenterProjections >= matrixSize/4)
 CenterProjections = matrixSize / 4;
 if (CenterProjections < 2)
 CenterProjections = matrixSize / 8;
 CenterProjections = (CenterProjections + 1) / 2 * 2;
 NumberOfSegments = (NumberOfSegments + 1) / 2 * 2;
 SegmentSize = (matrixSize - CenterProjections) / NumberOfSegments;
 CenterProjections = matrixSize - SegmentSize * NumberOfSegments;
 arraySize = matrixSize +
 (CenterProjections + SegmentSize) * NumberOfSegments;
 NR = NumberOfSegments + 1;
 ParxRelsParRelations("NR", Yes);
 PARX_change_dims("ArrayIndices",arraySize);
 ACQ_phase_encoding_mode[1] = User_Defined_Encoding;
 ParxRelsParRelations("ACQ_phase_encoding_mode", Yes);
 ACQ_spatial_size_1 = L[0] = arraySize;
 ParxRelsParRelations("ACQ_spatial_size_1", Yes);
 ACQ_user_filter = Yes;
 ParxRelsParRelations("ACQ_user_filter", Yes);
 ACQ_user_filter_mode = Special;
 ParxRelsParRelations("ACQ_user_filter_mode", Yes);
 ACQ_user_filter_memory = For_one_scan;
 sprintf(ACQ_user_filter_name,"KeyHoleAu");
 ACQ_user_filter_size[0] = 0;
 ACQ_user_filter_size[1] = arraySize;
 sprintf(ACQ_user_filter_setup_name,"NoOperation");
```

```

ACQ_user_filter_setup_size[0] = 1;
GO_reco_each_nr = Yes;

for (i=0; i<matrixSize; i++)
{
 ArrayIndices[i] = i - matrixSize/2;
}
idx = matrixSize;
for (i=0; i<NumberOfSegments; i++)
{
 for (j=0; j<SegmentSize/2; j++)
 {
 ArrayIndices[idx++] = j - CenterProjections/2 -
(i+1)*SegmentSize/2;
 }
 for (j=0; j<CenterProjections; j++)
 {
 ArrayIndices[idx++] = - CenterProjections/2 + j;
 }
 for (j=0; j<SegmentSize/2; j++)
 {
 ArrayIndices[idx++] = j + CenterProjections/2 +
(i)*SegmentSize/2;
 }
}
for (i=0; i<arraySize; i++)
 ACQ_spatial_phase_1[i] = 2.0 * ArrayIndices[i] / matrixSize;
DB_MSG(("<--KeyHoleRelations"));
...
void backbone(void)
{
 ...
 // make sure to get NI set to 1 just to make it fast
 STB_UpdateSliceGeoPars(1,1,1,minSliceThick);
 ...
 KeyHoleRelations();
}

```

### 5.2.3.7 The Pipeline Filter KeyHoleAu

```

#define DEBUG
#ifdef DEBUG
#define DB(P1) fprintf(fdTty,"%s",P1); (void)fflush(fdTty);
#define DB2(P1,P2) fprintf(fdTty,P1,P2); (void)fflush(fdTty);
#define DB3(P1,P2,P3) fprintf(fdTty,P1,P2,P3); (void)fflush(fdTty);
#else
#define DB(P1)
#define DB2(P1,P2)
#define DB3(P1,P2,P3)
#endif

/* global variables */

```

```

ACQ_SCAN_TYPE scan_type;
int *Image;
int *ArrayIndices;
int NumberOfArrayIndices;
int CenterProjections;
int SegmentSize;

#ifdef DEBUG
FILE *fdTty=NULL;
#endif

PIPE_INIT
PARX_INIT
{
#ifdef DEBUG
 fdTty = fopen("/dev/tty","w");
#endif

 /*-----*/
 /* allocate memory for image buffer and the array indices */
 /*-----*/
 DB2("allocating %d bytes for image buffer\n",scan_byte_size *
acq_size[1]);
 Image = malloc(scan_byte_size * acq_size[1]);
 if (NULL == Image)
 EXCEPT_printf("Cannot allocate memory for image buffer");
 NumberOfArrayIndices =
PARX_get_nth_dim(PARX_psid,"ArrayIndices",1);
 DB2("allocating %d integers for array
indices\n",NumberOfArrayIndices);
 ArrayIndices = malloc(NumberOfArrayIndices * sizeof(int));
 if (NULL == ArrayIndices)
 EXCEPT_printf("Cannot allocate memory for array indices");

 /*-----*/
 /* retrieve some parameters from parameter space */
 /*-----*/
 PARX_get_all_values(PARX_psid,
 "ArrayIndices", ArrayIndices,
 "CenterProjections", &CenterProjections,
 "SegmentSize", &SegmentSize,
 "ACQ_scan_type", &scan_type,
 NULL_PTR);
}
PIPE_PROC

int scanCnt;
int outCnt;
int idxImageIn;
int idxImageOut;
int dummy,pipeOutCnt;
idxImageOut = 0;

```

```

for (scanCnt=0; scanCnt <NumberOfArrayIndices ; scanCnt++)
{
 idxImageIn = *(ArrayIndices + scanCnt) + acq_size[1] / 2;
 outCnt = idxImageIn - idxImageOut;
 if (0 < outCnt)
 {
 /*-----*/
 DB3("\n transfer %d scans from 1st image at
%d\n",outCnt,idxImageOut);
 /*-----*/
 do
 {
 PIPE_WAIT_OUTPUT(outCnt,pipeOutCnt);
 memcpy(out_ptr,
 (char *) (Image + idxImageOut * acq_size[0]) ,
 pipeOutCnt * scan_byte_size);
 PIPE_SIGNAL_OUTPUT(0,pipeOutCnt);
 idxImageOut += pipeOutCnt;
 outCnt -= pipeOutCnt;
 } while (outCnt);
 }
 /*-----*/
 /* wait for next scan from acq_control and pass it to next filter
*/
 /*-----*/
 DB2(" %d ",idxImageOut);
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy);
 memcpy((char *) (Image + idxImageOut*acq_size[0]),in_ptr
,scan_byte_size);
 memcpy(out_ptr,(char *) (Image + idxImageOut*acq_size[0]),scan_byte_size);
 PIPE_SIGNAL_OUTPUT(1,1);
 idxImageOut++;
 if (scanCnt < (NumberOfArrayIndices - 1))
 {
 if (idxImageIn > *(ArrayIndices + scanCnt + 1) + acq_size[1] / 2)
 {
 outCnt = acq_size[1] - idxImageIn - 1;
 if (outCnt > 0)
 {
 /*-----*/
 DB3("\n transfer %d scans from 1st image at %d
\n",outCnt,idxImageOut);
 /*-----*/
 do
 {
 PIPE_WAIT_OUTPUT(outCnt,pipeOutCnt);
 memcpy(out_ptr,
 (char *) (Image + idxImageOut*acq_size[0]) ,
 pipeOutCnt * scan_byte_size);
 PIPE_SIGNAL_OUTPUT(0,pipeOutCnt);
 idxImageOut += pipeOutCnt;
 outCnt -= pipeOutCnt;
 } while (outCnt);
 }
 }
 }
}

```

```
 idxImageOut = 0;
 /*-----*/
 DB("\n **** Image transferred **** \n");
 /*-----*/
 }
}
}
PIPE_WRAPUP
If (NULL != fdTty)
{
 /*-----*/
 DB("\n **** Experiment Completed **** \n");
 /*-----*/
 (void)fclose(fdTty);
 fdTty = NULL;
}
PIPE_END
```

### 5.2.3.8 *Making keyHole ready to use*

```
make clean
make cproto
make depend
make install

cp keyhole.ppg $XWINNMRHOME/exp/stan/nmr/lists/pp
cp KeyHoleAu $XWINNMRHOME/exp/stan/nmr/au/src
$XWINNMRHOME/exp/stan/nmr/au/makeau KeyHoleAu
```

## 5.3 Exercises

- Setup and run keyhole experiments with different keyhole- setting

## 6 Magnetization Control

*Measuring method greMc using a Pipeline Filter for the online correction of scans having different magnetization amplitudes due to*

- *variable TR's (e.g. due to ECG-gating)*
- *non-steady-state (e.g. DS = 0)*
- *variable RF power*

### 6.1 Objectives

The purpose of this Pipeline Filter is to

- understand the pipeline mechanism for data reduction or data expansion
- implement Navigator techniques by the help of a Pipeline Filter

### 6.2 Implementation

#### 6.2.1 Parameters and their Settings

- ACQ\_user\_filter = Yes
- ACQ\_user\_filter\_name = ACQ\_user\_filter\_setup\_name = MagControl
- ACQ\_user\_filter\_memory = x
- ACQ\_user\_filter\_mode = Special
- ACQ\_user\_filter\_size[0] > ACQ\_size[0]
- ACQ\_user\_filter\_size[1] = 0

#### 6.2.2 Realization Hints

- The number of complex data points acquired during a delay is given by the equation  $nx = SW\_h * delay$
- NI is copied to the predefined parameter ni
- ACQ\_size[] is copied into predefined array parameter acq\_size[]
- The array parameter filter\_size[] is adapted to ACQ\_user\_filter\_size[]
- ACQ\_scan\_type is copied to the predefined parameter acq\_scan\_type
- ACQ\_scan\_type may have the values Scan\_Experiment or Setup\_Experiment during the execution time of an acquisition
- Use the construct

```

If (ACQ_scan_type == Setup_Experiment)
 { ...
 }

```

in order to have a special pulse program section for a setup experiment.



## 6.2.3 Creating the measuring method greMc

### 6.2.3.1 Using copyMethod to copy gre to greMC

```
>copyMethod ($XWINMRHOME/prog/service/copyMethod)
Enter name of source method : gre
Enter name of target method : greMc
Copy source method e.g. to /home/$USER/pv302/methods
>cd /home/$USER/pv302/methods/greMc
```

### 6.2.3.2 Modification of the Pulse Program

Original ppg section:

```
...
 (p0:sp0 ph1):f1
d3 grad{(t2) | r2d(t3) | (t1)+r3d(t4)}
d5 groff
denab REC_ENABLE grad{(t5) | (0) | (0)}
 ADC_INIT_B(ph1, ph0)
aqq ADC_START
d3 grad{(t8) | r2d(t6) | r3d(t7) }
d6 grad{(t8) | (0) | (0) }
d7 groff
d0 ADC_END
...
```

Modified ppg section

```
define delay aqScan
"aqScan = aqq - d5 - d3 - d4"
...
 (p0:sp0 ph1):f1
d3 grad{(0) | (0) | (t1)+r3d(t4)}
d3 REC_ENABLE groff
 ADC_INIT_B(ph1, ph0)
d5 ADC_START
d3 grad{(t2) | r2d(t3) | (0)}
denab grad{(t5) | (0) | (0)}
aqScan
d3 grad{(t8) | r2d(t6) | r3d(t7)}
d6 grad{(t8) | (0) | (0) }
d7 groff
d0 ADC_END
if (ACQ_scan_type == Setup_Experiment)
{
 3s ; additional T1 recovery for Setup and Auto_RG
}
...
```

The scan time has been extended by d3+d5+d4, since denab is just a compensation for d4 for the hidden delay de-depa.

### 6.2.3.3 Extensions to *parsDefinition.h* and *parsLayout.h*

None.

### 6.2.3.4 Extension of *initMeth.c*

```
int dimRange[2] = {2,2}; //limit experiment to two dimensions
```

### 6.2.3.5 Extensions to *parsRelations.c*

```
void backbone(void)
{
...
// make sure to get NI set to 1
STB_UpdateSliceGeoPars(1,1,1,minSliceThick);
...
}
```

### 6.2.3.6 Extensions to *BaseLevelRelations.c*

```
void SetBaseLevelParam()
{
...
SetPipeFilterParameters();
DB_MSG(("<--SetBaseLevelParam"));
}
void SetPipeFilterParameters(void)
{
ACQ_user_filter = Yes;
ParxRelsParRelations("ACQ_user_filter", Yes);
sprintf(ACQ_user_filter_name,"MagControl");
sprintf(ACQ_user_filter_setup_name,"");
ACQ_user_filter_memory = For_one_scan;
ACQ_user_filter_mode = Special;
ParxRelsParRelations("ACQ_user_filter_mode", Yes);
ACQ_user_filter_size[0] = ACQ_size[0] +
(int)(SW h * (D[3] + D[5] + D[4]) + 0.5)*2;
ACQ_user_filter_size[1] = 0;
DS = NDummyScans = 0; // Disable dummy scans at all
}
```

to round up  
additional points

### 6.2.3.7 The Pipeline Filter *MagControl*

```
// Magnetization Control for non-steady-state acquisition

double d5, sw_h;
double *meanRef;
int magSize;

PIPE_INIT
PARX_INIT
/*-----*/
```

```

/* retrieve some parameters from parameter space */
/*-----*/
PARX_get_all_values(PARX_psid,
 "D[5]", &d5,
 "SW_h", &sw_h,
 NULL_PTR);

magSize = (int)(sw_h * d5)*2; // # of points recording the
Magnetization
if (magSize < 2) magSize = 2;
if (acq_scan_type == Setup_Experiment)
 if (magSize > acq_size[0]/4) magSize = acq_size[0]/4;
meanRef = malloc(ni * sizeof(double));
if (NULL == meanRef) EXCEPT_printf("Cannot allocate memory for mean
buffer");

PIPE_PROC
int niStep, niIdx, curNi, scan_cnt, dummy, *ip, *op;
double mean, Re, Im;

while (1)
{
 for (scan_cnt=0; scan_cnt < scans_per_exp; scan_cnt++)
 {
 PIPE_WAIT_INPUT_OUTPUT(1,1,dummy,dummy);
 ip = (int *) in_ptr;
 op = (int *) out_ptr;
 if (acq_scan_type == Setup_Experiment)
 {
 // Pass Magnetization during
 Setup
 for (il=0; il < magSize; il++)
 *op++ = *(ip+il);
 ip = ip + filter_size[0] - acq_size[0];
 for (il=0; il < scan_word_size - magSize; il++)
 *op++ = *ip++;
 }
 else
 {
 // Calculate mean during Scan
 niStep = scan_cnt / ni;
 niIdx = scan_cnt % ni;
 mean = 0.0;
 for (il=0; il < magSize / 2; il++)
 {
 Re = *(ip+il*2);
 Im = *(ip+il*2+1);
 mean += sqrt(Re*Re + Im*Im);
 }
 mean = mean / magSize / 2;
 if (niStep == 0)
 meanRef[niIdx] = mean; // meanRef[niIdx] = mean of
first NI Scans
 ip = ip + filter_size[0] - acq_size[0];
 for (il=0; il < scan_word_size; il++) // Pass normalized Scan
 *op++ = (int)(meanRef[niIdx] / mean * *ip++);
 }
 PIPE_SIGNAL_OUTPUT(1,1);
 }
 if (acq_scan_type == Scan_Experiment) break;
}

```

```
 scans_in = 0;
} /* while (1) */
PIPE_WRAPUP
PIPE_END
```

### 6.2.3.8 Making greMc ready to use

```
make clean
make cproto
make depend
make install

cp greMc.ppg $XWINNMRHOME/exp/stan/nmr/lists/pp
cp MagControl $XWINNMRHOME/exp/stan/nmr/au/src
$XWINNMRHOME/exp/stan/nmr/au/makeau MagControl
```

### 6.2.4 Exercises

- Run a GSP experiment and make a SnapShot of the acquisition display
- Set manually `ACQ_size[0] = ACQ_user_filter_size[0]`, deactivate the pipeline filter by setting `ACQ_user_filter=No` and run a GSP experiment in order to have the whole Scan displayed
- Compare the SnapShot with the whole Scan
- **Experiment 1:** Setup a **short TR** and a **30 degree excitation pulse** and run the sequence **without dummy scans** by setting `PVM_RepetitionTime = min` and `ExcPulse.FlipAngle = 30.0 deg`
- **Experiment 2:** **Move the K-space center to the beginning of the acquisition** by setting `ACQ_phase_enc_start[1] = 0`
- **Experiment 3:** Apply a **90 degree excitation pulse** instead of an alpha pulse by setting `ExcPulse.FlipAngle = 90.0 deg`
- Disable the normalization in **MagControl** and run the experiments again
- Compare the scaled with the un-scaled experiments

## 7 References

7.1 **Advanced Users Manual, Chapter D-9.3: AU filters for Pipelined Acquisition**


7.2 **Advanced Users Manual, Chapter D-9: Automation Programming**

7.3 `$XWINNMRHOME/exp/stan/nmr/au/vorspann_f`

Definition of the "Convenience Parameters"

7.4 `$XWINNMRHOME/prog/include/aupipe.h`

Macro definitions and initialization of the "Convenience Parameters"



*adventurment*

**Method "ltGre"**

